REF: CG0351 10 x 1 ml
CG0352 10 x 1 ml

Coagnos[®] Routine Control GB

INTENDED PURPOSE:

The Coagnos[®] Routine Control N and Coagnos[®] Routine Control P are intended for use as a quality control material. The Coagnos[®] Routine Control kit is indicated for professional use only.

Coagnos[®] Routine Control N and Coagnos[®] Routine Control P are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assayed for Fibrinogen, TT and ATXa, and are prepared from normal human plasma.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody, HCV antibody. However they should be handled with the same precautions as a human patient sample. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
Routine Control N	10 x 1 mL	Prepared from pooled normal plasma. Each vial contains 1 mL of buffered, lyophilised human plasma.	Reconstitute each vial of the appropriate control with 1 mL of distilled or deionised water. Swirl gently. Allow to stand for 10 minutes for complete dissolution and mix well before use.
Routine Control P	10 x 1 mL	Prepared from adsorbed human plasma. Each vial contains 1 mL of buffered, lyophilised human plasma.	Reconstitute each vial of the appropriate control with 1 mL of distilled or deionised water. Swirl gently. Allow to stand for 10 minutes for complete dissolution and mix well before use.

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 8 hours when kept at +2 –+8°C or 4 weeks at -20°C when flash frozen. Keep covered.

SAMPLE COLLECTION AND PREPARATION:

Not applicable.

PROCEDURE:

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS:

Routine Control N should give values within the laboratory normal range for PT, aPTT and fibrinogen assays. Routine Control P has been standardised to give prolonged and markedly prolonged PT and aPTT times respectively. Lot and instrument specific expected values are provided with each pack of controls.

LIMITATIONS:

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{1,2,3}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives. Each laboratory should establish its own performance data.

Precision

Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

Coagnos[®] PT

	Routine Control N			Routine Control P		
	Mean	SD	CV [%]	Mean	SD	CV [%]
Clots [s]	11.8	0.31	2.6	22.3	0.86	3.87
%PT	114.1	2.67	2.34	38.6	1.79	4.65
INR	0.92	0.02	1.86	1.94	0.05	2.72

Coagnos[®] APTT

	Routine Control N			Routine Control P		
	Mean	SD	CV [%]	Mean	SD	CV [%]
Clots [s]	29.8	0.54	1.81	48.7	0.94	1.92
Ratio	1.11	0.03	2.47	1.80	0.04	2.13

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !
	Batch identification number
	Item number
	Manufacturer

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