REF: CG0361 10 x 1 ml
CG0362 10 x 1 ml

Coagnos® Special Control GB

INTENDED PURPOSE:

The Coagnos® Special Control N and Coagnos® Special Control P are intended for use as a quality control material. Coagnos® Special Control N may be used as a normal control, and Coagnos® Special Control P may be used as an abnormal control when assaying for Factors II, V, VII, VIII, IX, X¹, XI, XII², fibrinogen³, von Willebrand Factor, antigenic and functional Protein C and Protein S (total and free), as well as the chromogenic assays including Antithrombin III, Protein C, Factor VIII and Plasminogen. The Coagnos® Special Control kit is indicated for professional use only. Coagnos® Special Controls are used in factor assays and other testing in the same manner as a fresh plasma pool. Factor II, VII, VIII, IX and X values and the chromogenic Factor VIII, Antithrombin Xa and Protein C values are traceable to World Health Organisation standards to ensure the utmost credibility in values stated⁴. The control plasma should be used to gauge internal factors in each laboratory's system.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody, HCV antibody. However they should be handled with the same precautions as a human patient sample. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
Special Control N	10 x 1 mL	Prepared from a frozen pool of citrated plasma and is buffered and lyophilised to ensure stability of all plasma constituents ⁵ .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.
Special Control P	10 x 1 mL	Prepared from a frozen pool of citrated plasma and is buffered and lyophilised to ensure stability of all plasma constituents ⁵ .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Component	Stability	Signs of Deterioration
Special Control N	Values for Factor VIII, von Willebrand factor and ristocetin co-factor are stable for 2 hours at +2 → +8°C. All other factors are stable for 4 hours at +2 → +8°C or 4 weeks at -20°C when flash frozen.	Unreconstituted Special Controls must appear as a light yellow, dry plug. If any unusual conditions are noted, the customer should notify Analyticon Biotechnologies AG before using.
Special Control P		

SAMPLE COLLECTION AND PREPARATION:

Not applicable.

PROCEDURE:

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS:

The assayed percent activities of the various coagulation factors should be used as a guideline. If the product does not perform as expected, the patient results should be considered invalid. Ensure the lot number printed on this assay sheet is the same as that on the vial of Coagnos® Special Control to be used.

LIMITATIONS:

The results obtained with Special Control N and Special Control P depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{6,7,8}. Each laboratory should establish an expected range for the particular instrument-reagent system used.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES:

Refer to the reference values insert for lot specific reference ranges. Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Parameter	Special Control N		Special Control P	
	Mean	CV (%)	Mean	CV (%)
Fibrinogen (g/L)	2.9	6.9	1.24	3.35
Factor IX (%)	132.7	4.3	50.3	5.4
Protein S (%)	110.4	5.0	42.4	3.1

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !
	Batch identification number
	Item number
	Manufacturer

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REF: CG0361 10 x 1 ml
CG0362 10 x 1 ml

**BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/
LITERATURVERZEICHNIS:**

1. Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, AJCP, 64:817-819
2. Hardisty RM, et al. (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma, Thrombosis et Diathesis Haemorrhagica, 7:215-229
3. Morse EE, et al. (1971) Automated Fibrinogen Determination, AJCP, 55:671-676
4. Elodi S, Katalin V, Hollan S (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, Haemostasis, 7:1-9
5. Thelin M (1968) Preparation and Standardization of a Stable AHF Plasm, Thrombosis et Diathesis Haemorrhagica, 19:423
6. Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568
7. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564
8. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235