Coagnos® Thrombin Time GB

INTENDED PURPOSE:

The Coagnos® Thrombin Time kit is intended for carrying out clot based haemostasis assays in human plasma. The Coagnos® Thrombin Time kit is indicated for professional use only.

The Thrombin Time reagent contains bovine thrombin. The reagent can be used manually or on semi-automated and automated instruments. The test is commonly applied to detect various sources of interference with normal blood coagulation. Prolongation of the thrombin clotting time can be taken as a qualitative indication of abnormal fibrinogen levels (high or low), or the presence of interfering substances such as Fibrin Degradation Products (FDPs) or heparin. Quantitative evaluation of the possible causes of prolonged thrombin clotting time should be performed as follow-up studies, such as aPTT or chromogenic assay for heparin, Clauss fibrinogen, FDP determinations, heparin neutralisation by protamine sulphate or polybrene¹, normal plasma mixing studies² or reptilase assay³ to distinguish between hypofibrinogenaemia and FDP effects.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

| Component | Content | Description | Preparation |
|---------------|-----------|--|--|
| Thrombin Time | 10 x 2 mL | Each vial contains a lyophilized preparation of bovine thrombin with buffers and stabilisers. The reconstituted reagent contains ≤10 NIH units/mL of thrombin. The reagent should be a white lyophilised plug. | Reconstitute the reagent with the recommended volume of purified water: 1 mL - REF CG0311 Allow to stand for 5 minutes then mix gently by inversion and transfer to a plastic tube. Reconstituted reagent should be a clear, colourless solution. |

STORAGE, SHELF-LIFE AND STABILITY:

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

| | |
|---------------|--|
| Thrombin Time | Reconstituted reagent should be stored at +2 –+8°C and is stable for 14 days or at -20°C for 1 month. Reconstituted reagent can also be stored on-board a Sysmex CA1500 photo-optical instrument for up to 7 days. |
|---------------|--|

SAMPLE COLLECTION AND PREPARATION:

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at +2 –+8°C or +18 –+24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at +37°C prior to testing. Do not keep at +37°C for more than 5 minutes⁴.

PROCEDURE:

1. Collect and prepare the blood sample according to the directions outlined in "Sample Collection And Preparation".
2. Reconstitute the control plasma according to the package insert provided with each control. Prepare the reagent for use in the procedure according to the reconstitution instructions in the "Composition" section.
3. Perform all tests in duplicate. Calculate the mean clotting time of the duplicate determinations to the nearest 0.1 seconds. Individual values should be within ±5% of the mean value.

Manual Method

1. Pipette 0.2 mL of the patient plasma or control plasma into a reaction tube.
2. Incubate at +2. 37°C for 3 minutes.
3. Pipette 0.1 mL of Thrombin Time into the reaction tube containing patient or control plasma whilst simultaneously starting a timer.
4. Record the time taken for clot formation to the nearest 0.1 second.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Analyticon Biotechnologies AG for instrument specific application guides.

INTERPRETATION OF RESULTS:

The results of the Thrombin Time test should be reported to the nearest 0.1 seconds. The normal range (usually the mean ±2 standard deviations) should be established by each laboratory. Results outside the normal range should be considered abnormal and follow-up testing performed.

LIMITATIONS:

Expected values for the Thrombin Time test will vary from one laboratory to another depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and

method of sample storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional testing. As well as the cause of elongated thrombin times indicated, a report has suggested that many systemic amyloidosis patients with bleeding complications may have a circulating inhibitor which prolongs the Thrombin Time⁵. Also, therapeutic levels of heparin may entirely abolish clotting in the Thrombin Time test, although neutralisation with protamine sulphate or polybrene should correct the Thrombin Time¹.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Analyticon Biotechnologies AG supplies the following controls available for use with this product:

REF CG0351 Coagnos® Routine Control N
REF CG0352 Coagnos® Routine Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using a Coagulyzer® Auto Pro instrument. Each laboratory should establish its own performance data.

Precision

Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

| Sample | Mean | Repeatability | | Within-device | |
|-----------|------|---------------|--------|---------------|--------|
| | | SD | CV [%] | SD | CV [%] |
| Clots [s] | 11.0 | 0.27 | 2.49 | 0.56 | 5.12 |
| | 18.0 | 0.48 | 2.67 | 1.02 | 5.66 |
| | 24.9 | 1.00 | 4.03 | 1.66 | 6.67 |
| %PT | 0.95 | 0.02 | 2.49 | 0.05 | 5.12 |
| | 1.56 | 0.04 | 2.67 | 0.09 | 5.66 |
| | 2.15 | 0.09 | 4.03 | 0.14 | 6.67 |

Interferences

Coagnos® Thrombin Time is insensitive to Heparin levels of up to 0.12 U/ml. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 99 mg/dL. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 40 mg/dL for Coagnos® PT. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 300 mg/dL. Lipid concentrations in excess of this prevent clot detection.

Method comparison

Comparison of clot time in seconds were determined using Coagnos® Thrombin Time on Coagulyzer® Auto Pro and Siemens® Thrombin Time with Thromboclotin® on a Sysmex® CS-5100 system on 110 samples. The following correlations were obtained:

Coagnos® Thrombin Time [s] = 1.050x – 3.421 r2 = 0.59 n = 110

SYMBOLS

| | |
|---|---|
|  | In vitro diagnostics product |
|  | The product complies with European legislation |
|  | Follow the instructions for use ! |
|  | Use by |
|  | Permitted storage temperature range |
|  | Read warnings and precautions in instructions for use ! |
|  | Batch identification number |
|  | Item number |
|  | Manufacturer |

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BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Laposata et al. The Clinical Haemostasis Handbook, Yearbook Medical Publishers Inc., p219, 1989
2. Thompson AR and Harker LA. Manual of Haemostasis and Thrombosis. 3rd Ed., F.A. Davis Co., p62, 1983
3. DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990.
4. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
5. Gastineau DA et al. (1991) Inhibitor of the Thrombin Time in Systemic Amyloidosis: A Common Coagulation Abnormality, Blood, 77:2637-40