### Order information:

Catalog No.	Contents
5799	R1 1 x 500 ml

### Intended use:

Solution for deproteinisation of human blood, serum, plasma and liquor in glucose assays.

It is also suitable for sample pre-treatment in alcohol and lactate assays.

### Summary:

Glucose is a fast decomposed compound in samples. In order to stabilize Glucose the decomposing enzymes should be removed from the sample by deproteinisation and zentrifugation.

### Test principle:

Proteins like enzymes in the sample are denaturated by the deproteinisation solution and are sedimented by centrifugation. The supernatant can be used as sample in the assay.

### Reagent concentration:

<b>R1:</b>	
Perchloric acid	369 mmol/l
Na-Perchlorate	157 mmol/l

### Preparation and Stability:

R1: Ready for use

Reagent is stable up to the expiration date at +2 to +25°C.

### Specimen:

Liquor  
Capillary blood, EDTA- / Heparin-Plasma, Serum  
Deproteinise samples directly after collecting. Separate supernatant immediately after centrifugation.

Stability:	1 days	at +20°C to +25°C
	7 days	at +4°C to +8°C

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Testing procedure:

#### Materials provided

- Working solutions as described above

#### Additional materials required

- Centrifuge
- Centrifuge tubes

<b>Manual procedure:</b>		
	Macro	Semi micro
R1	1000 µl	500 µl
Sample	100 µl	50 µl
Mix well and centrifuge at 10000-13000rpm for 10min.		

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2. Test Fibel, Boehringer Mannheim GmbH 1977, Mannheim