Ery-Solution

GOWERS' METHOD

Order information:

Catalog-No.	Contents
6944	R1 1 x 2500 ml
6945	R1 1 x 500 ml

Intended use:

Test for photometrical determination of erythrocyte number in blood. Test can be performed in round and square bottom cuvettes.

Test principle:

The number of erythrocytes in a defined blood volume is measured photometrically via turbidity.

Reagents – contents and concentrations:

Ery-solution. Based on GOWERS' method:

Sodium sulphate	0.19 mol/l
Acetic acid	2.90 mol/l

Additional material required:

Capillary tubes 10µl for test procedure in round bottom cuvettes
Capillary tubes 5µl for test procedure in square bottom cuvettes.

Preparation and stability:

R1: ready to use

Stability:

Up to the expiration date at +15°C to +25°C

Close bottles tightly after use.

Specimen:

Use fresh capillary blood or blood with di-Potassium-EDTA as anticoagulant.

Stability of EDTA blood:

At room temperature: 8 hours

At +2° C to +8° C 24 hours

Note:

Keep away from children. After the stated expiry, the reagent should not be used. Do not use cuvettes with cloudy or brownish coloured reagents. Cuvettes should only be touched at ruffled sides or cap. For measurement, place smooth sides in measurement direction. Hemolysis can interfere.

Avoid strong squeezing of the finger pad during retrieval of capillary blood because sample could be diluted by lymph. It is recommended to determine own laboratory reference values.

There can be deviations in case of pathologically altered blood (e.g. deformation of erythrocytes, extreme changes of size and shape of erythrocytes and pathological decreased cell counts under 3.5 mio./µL). In these cases a microscopic count is recommended.

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Use of round bottom cuvettes, 1cm diameter:

Reagent	2.45 ml
Sample	10 µl

1. open the cuvette
2. pipette 10 µl blood into the cuvette (e.g. use capillary tube)
3. rinse pipette tip
4. close the cuvette and mix well. Avoid formation of foam.
5. measure within 3 to 15 min.

Use of square bottom cuvettes, 1cm:

Reagent	2.15 ml
Sample	5 µl

1. open the cuvette
2. place the blood filled capillary tube inside the cuvette and close the cap.
3. directly remove blood out of capillary tube by rotation of cuvette, avoid the formation of foam.
4. measure within 3 to 15 min. Capillary tube has to be in one corner of cuvette during measurement.

Measurement:

Wavelength:	546 nm
Temperature:	ambient temperature (+15°C bis +25°C)
Cuvette (light path):	1 cm

Calculation:

Using a COMPUR or DR.LANGE photometer the measured value is directly shown as Mio. Ery/µL. (For further information please refer to user manual)

Using other types of photometers the result is shown as Mio. Ery/µL based on a calibration curve or a photometer-specific factor.

Normal values:

Men	4.5 - 5.9	Mio./µL
Women	4.1 - 5.1	Mio./µL
Children (age-depend.)	3.2 - 5.8	Mio./µL

Quality control:

All commercially available blood controls can be used.

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