Ery-cuvettes square bottom

GOWERS' METHOD

Order information:

Catalog No.	Contents
6956	40 Single-use square cuvettes

Intended use:

Test for photometrical determination of erythrocyte number in blood.

Test principle:

The number of erythrocytes in a defined blood volume is measured photometrically via turbidity.

Reagents – contents and concentrations:

40 single-use square bottom cuvettes filled with ready to use Ery-solution.

Based on GOWERS' method:

Sodium sulphate	0.19 mol/l
Acetic acid	2.90 mol/l
40 capillary tubes 5µl, heparinised	

Preparation and stability:

Cuvettes are ready for use

Stability:

Sealed packages	Up to the expiration date	at +15°C to +25°C
Opened packages	24 hours	at +15°C to +25°C
	8 days	at +2°C to +8°C

Close aluminium bags tightly with a clip for storage.

Specimen:

Use fresh capillary blood or blood with di-Potassium-EDTA as anticoagulant.

Stability of EDTA blood:

At room temperature:	8 hours
At +2° C to +8° C	24 hours

Note:

Keep away from children. After the stated expiry, the reagent should not be used. Do not use cuvettes with cloudy or brownish coloured reagents. Cuvettes should only be touched at ruffled sides or cap. For measurement, place smooth sides in measurement direction. Hemolysis can interfere.

Avoid strong squeezing of the finger pad during retrieval of capillary blood because sample could be diluted by lymph. It is recommended to determine own laboratory reference values.

There can be deviations in case of pathological modified blood (e.g. deformation of erythrocytes, extreme changes of size and shape of erythrocytes and pathological decreased bloods under 3.5 mio./µL). In these cases, a microscopic count is recommended.

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

1. open the cuvette
2. place the blood filled capillary tube inside the cuvette and close the cap.
3. directly remove blood out of capillary tube by rotation of cuvette, avoid the formation of foam.
4. measure within 3 to 15 min. Capillary tube has to be in one corner of cuvette during measurement.

Measurement:

Wavelength:	546 nm
Temperature:	ambient temperature (+15°C to +25°C)
Cuvette (light path):	1 cm

Calculation:

Using a COMPUR photometer the measured value is directly shown as Mio. Ery/µL. (For further information please refer to user manual)

Using other types of photometers the result is shown as Mio. Ery/µL based on a calibration curve or a photometer-specific factor.

Normal values:

Men	4.5 - 5.9	Mio./µL
Women	4.1 - 5.1	Mio./µL
Children (age-depend.)	3.2 - 5.8	Mio./µL

Quality control:

All commercially available blood controls can be used.

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