### BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B5931	200 / 600	R1	3 x	20 ml
		R2	3 x	15 ml
B5903	300 / 600*	R1	2 x	10 ml
		R2	2 x	6 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

### Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of ferritin in human serum and plasma.

### Summary:

Ferritin is the iron storage protein. It has a molecular weight of >440000 daltons, depending upon the iron content, and consists of a protein shell (apoferritin) of 24 subunits and an iron core containing an average of approx. 2500 Fe<sup>3+</sup> ions (in the basic isoforms). Common to all isoforms is their construction from two separate subunits, the acidic H (heavy)-type subunit and the weakly basic L (light)-type subunit. The basic isoforms are responsible for the long-term iron storage function and are mainly detectable in the liver, spleen and bone marrow. Acid isoforms are found mainly in the myocardium, placenta, tumor tissue and – to a lesser extent – in the depot organs.

The determination of ferritin is necessary above all in iron metabolism diagnosis, monitoring iron therapy, ascertaining the iron reserves in groups at risk and in the differential diagnosis of anaemia. It encompasses prelatent and latent iron deficiency as well as iron overloading. It is also used to distinguish between hypoferric anaemia and hypochromic anaemia (chronic infection and tumor anaemias, sideroblastic anaemia or thalassemia).

Ferritin determinations are particularly suitable for monitoring renal anaemia when iron utilization and distribution disorders are present during therapy with erythropoietin. The ferritin detectable in blood is in equilibrium with the body's depot iron and hence acts as an indicator for the level of the iron stores.

The Analyticon ferritin assay is based on the immunological agglutination principle with enhancement of the reaction by latex.

### Test principle:

Particle-enhanced immunoturbidimetric assay. Anti-ferritin antibodies bound to latex micro particles react with the antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically.

### Reagent concentration:

**R1:**  
Glycin-buffer pH 9.0 180 mmol/l  
NaCl 100 mmol/l  
Preservative

**R2:**  
Latex particles coated with polyclonal anti-human ferritin antibodies 0.12 %  
Preservative

### Preparation and stability:

R1: Ready for use.  
R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C.

Onboard stability at +2°C to +8°C:  
R1: 60 days  
R2: 60 days

### Note:

Mix the latex reagent well before placing on the analyzer and once every 5 days thereafter.

### Specimen:

Collect serum using standard sampling tubes.  
Li-heparin, EDTA- or citrate plasma.  
Stability: 7 days at +20°C to +25°C  
7 days at +4°C to +8°C

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.  
Icterus: No significant interference up to an index I of 100 (approximate conjugated and unconjugated bilirubin concentration: 100 mg/dl).  
Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/dl).  
Lipemia (Intralipid): No significant interference up to an index L of 145 (approximate triglycerides concentration: 290 mg/dl).  
There is poor correlation between turbidity and triglycerides concentration.  
Rheumatoid factors up to 500 IU/ml do not interfere.  
A high-dose hook-effect may occur at ferritin concentrations above 1100 ng/ml.  
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

*Materials provided*  
• Working solutions as described above  
*Additional materials required*  
• Calibrators and controls as indicated below  
• 0.9% NaCl-solution

### Measuring / reportable range:

Measuring range 2 - 1000 ng/ml (µg/L)  
Determine higher ferritin concentration via the rerun function using 0.9% NaCl.

### Expected values:

Men	30-400 ng/ml
Women	15-150 ng/ml
Children 3 months-16years	20-200 ng/ml
Children 2 <sup>nd</sup> -3 <sup>rd</sup> month	80-500 ng/ml
Children 1 <sup>st</sup> month	150-450 ng/ml

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the ferritin results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit):

Detection limit: 2 ng/ml (2 µg/l)  
The lower detection limit represents the lowest measurable ferritin- concentration that can be distinguished from zero.

### Imprecision:

Reproducibility within run was determined using human control sera (n = 20). The following results were obtained:

Sample	Within run		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	260.7	5.05	1.9
Sample 2	210.8	4.94	2.3
Sample 3	83.9	2.93	3.5

Sample	Run to run		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	244.5	5.64	2.3
Sample 2	115.3	3.12	2.7
Sample 3	393.3	12.86	3.3

### Method comparison:

A comparison of the Analyticon Turbitex® Ferritin (y) with a commercial obtainable assay (x) gave the following result:  
y = 0.9162x + 2.1541; r = 0.9957

### Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Traceability: This method has been standardized against WHO reference material (1<sup>st</sup> International Standard for Ferritin WHO80/602).

Calibration Type: Logit-Log  
S1: NaCl 0.9%  
S2-S5: Bio Cal® FER Calibration set 4 x 2 ml #14555

### Calibration frequency

Full calibration is recommended:  
• after reagent lot change  
• as required following quality control procedures

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