Order information:

Catalog No.	Contents					
4560	R1	2 x	10 ml	R2	2 x	6 ml
H5901 Hit I (ILab*)	R1	3 x	20 ml	R2	3 x	10 ml
H5903 Hit 917 (AU*)	R1	2 x	10 ml	R2	2 x	5 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 361-400 (user defined method)
 Hitachi 917: ACN/BCN 901-905 (user defined method)
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support

Intended use:

Immuno-turbidimetric assay for the in vitro quantitative determination of ferritin in human serum and plasma.

Summary:

Ferritin is the iron storage protein. It has a molecular weight of >440000 daltons, depending upon the iron content, and consists of a protein shell (apoferritin) of 24 subunits and an iron core containing an average of approx. 2500 Fe³⁺ ions (in the basic isoforms). Common to all isoforms is their construction from two separate subunits, the acidic H (heavy)-type subunit and the weakly basic L (light)-type subunit. The basic iso-ferritins are responsible for the long-term iron storage function and are mainly detectable in the liver, spleen and bone marrow. Acid iso-ferritins are found mainly in the myocardium, placenta, tumor tissue and – to a lesser extent – in the depot organs.

The determination of ferritin is necessary above all in iron metabolism diagnosis, monitoring iron therapy, ascertaining the iron reserves in groups at risk and in the differential diagnosis of anaemia. It encompasses prelatent and latent iron deficiency as well as iron overloading. It is also used to distinguish between hypoferric anaemia and hypochromic amenity (chronic infection and tumor anaemias, sideroblastic anaemia or thalassemia).

Ferritin determinations are particularly suitable for monitoring renal anaemia when iron utilization ad distribution disorders are present during therapy with erythropoietin. The ferritin detectable in blood is in equilibrium with the body's depot iron and hence acts as an indicator for the level of the iron stores. A variety of methods are available for determining ferritin, e.g. radioimmunoassay (RIA), enzyme-linked immunosorbent assay (ELISA), fluorescence immunoassay (FIA), luminescence immunoassay (LIA) and nephelometric immunoassay.

The Analyticon ferritin assay is based on the immunological agglutination principle with enhancement of the reaction by latex.

Test principle:

Particle-enhanced immunoturbidimetric assay.

- Sample and addition of R1 (buffer)
 - Addition of R2 (Anti-ferritin antibody latex) and start of reaction:
- Anti-ferritin antibodies bound to latex micro particles react with the antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically.

Reagent concentration:

R1:	
Glycin-buffer pH 9,0	180 mmol/l
NaCl	100 mmol/l
Preservative	
R2:	
Latex particles coated with polyclonal anti-human ferritin antibodies	0.12 %
Preservative	

Preparation and stability:

R1: Ready for use.
 R2: Ready for use.
 Unopened kit components: Up to the expiration date at +2°C to +8°C.
 Onboard stability at +2°C to +8°C:
 R1: 60 days
 R2: 60 days

Note:

Mix the starter reagent well before placing on the analyzer and once every 5 days thereafter.

Specimen:

Collect serum using standard sampling tubes.
 Li-heparin, EDTA- or citrate plasma.
 Stability: 7 days at +20°C to +25°C
 7 days at +4°C to +8°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
 Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
 Icterus: No significant interference up to an index I of 60 (approximate conjugated and unconjugated bilirubin concentration: 60mg/dl).
 Hemolysis: No significant interference up to an index H of 800 (approximate haemoglobin concentration: 800 mg/dl).
 Lipemia (Intralipid): No significant interference up to an index L of 750 (approximate triglycerides concentration: 1500 mg/dl).
 There is poor correlation between turbidity and triglycerides concentration.
 Rheumatoid factors up to 500 IU/ml do not interfere.
 A high-dose hook-effect may occur at ferritin concentrations above 1100 ng/ml.
 The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl-solution

For performing the assay on Hitachi-analyzers refer to the appropriate operator's manual for user defined methods.

Hitachi 917: Prior to testing, put the barcodes according to your chosen application code on the bottles.

e.g. ACN 901 put R1 barcode 901 on bottle R1 and R2 barcode 901 on bottle R2.
 Place R1 (buffer) on reagent plate 1; R2 (Ferritin-antibody latex) on reagent plate 2.

Manual procedure:

Wavelength:	575 nm
Temperature:	+37°C
Zero adjustment:	against blank
R1 Blank/Sample/Calibrator	Blank/Sample/Calibr. 800 µl 40 µl
Incubate 5 min. Then add:	
R2	400 µl
Incubate 45s and read initial absorbance (blank). Read absorbance again after 2min (Endpoint).	

Note: Reagent and sample volumes may have to be adjusted to the specific system used.

Measuring / reportable range:

Measuring range 8-1000 ng/ml (µg/L)

When the ferritin concentration of the sample is above the measuring range, manually dilute with 0.9% NaCl (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. 3) or determine with rerun-function.

Expected values:

Men	30-400 ng/ml
Women	15-150 ng/ml
Children 3 months-16years	20-200 ng/ml
Children 2 nd -3 rd month	80-500 ng/ml
Children 1 st month	150-450 ng/ml

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the ferritin results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 5 ng/ml (5 µg/l)

The lower detection limit represents the lowest measurable ferritin- concentration that can be distinguished from zero.

Imprecision:

Reproducibility within run was determined using human control sera (n = 21). The following results were obtained:

Sample	Within run		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	14.90	0.6	4.0
Sample 2	100.00	0.65	0.6
Sample 3	431.05	2.20	0.5

Sample	Between day		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	16.47	0.87	5.3
Sample 2	105.18	1.60	1.5
Sample 3	428.71	3.52	0.8

Method comparison:

A comparison of the Analyticon Turbitex® Ferritin (y) with a commercial obtainable assay (x) with 64 samples gave the following result:

$$y = 0.89 x - 9.4; r = 0.997$$

Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Traceability: This method has been standardized against WHO reference material (1st International Standard for Ferritin WHO80/602).

S1:	NaCl 0.9%		
S2-S5:	Bio Cal® FER	1 x 2 ml	#14550
	Bio Cal® FER Calibration Set	4 x 2 ml	#14555

Calibration frequency

Full calibration is recommended:

- after reagent lot change
- as required following quality control procedures

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