### Order information:

Catalog No.	Contents
H5601 Hit 1 / 917 (ILab* / AU*)	R1 6 x 14 ml   R1a 6 x 6 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911/917: ACN 667  
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

In vitro test for the quantitative determination of glycated protein (fructosamine) in human serum and plasma on automated clinical chemistry analyzers. Measurement of fructosamine may be a suitable alternative to measurement of HbA1c in patients with Hb variants.

### Summary:

Fructosamine is a time-averaged indicator of blood glucose levels and is used to assess the glycemic status of diabetics. The concentration of glycated proteins such as glycohemoglobin, glycoalbumin or glycated total protein is generally recognized to be valuable in evaluating the glycemic status of diabetic patients. Various other methods used for such determinations such as affinity chromatography and the thiobarbituric acid method are labour-intensive and time consuming and results obtained from different laboratories are difficult to compare. This assay is based on the Nitrotetrazolium-blue method and yields a precise and easily automatable determination of the non-enzymatic glycation of serum proteins. The turnover of serum proteins (albumin has a half-life of 19 days) is less than that of hemoglobin (lifespan of erythrocytes is approx. 120 days), and therefore fructosamine determinations provide a means of monitoring patient blood glucose status over a shorter period (1-3 weeks) than glycohemoglobin (6-8 weeks). As a result, changes in fructosamine values alert the physician to deteriorating glycemic control earlier than changes in HbA1c values. In addition, fructosamine levels decrease more quickly than HbA1c when diabetic patients are brought under better control.

### Test principle:

Colorimetric assay  
This colorimetric assay is based on the ability of ketoamines to reduce nitro-tetrazolium-blue (NB) to formazan in an alkaline solution. The rate of formation of formazan is directly proportional to the concentration of fructosamine. Uric acid interference is eliminated by Uricase and detergent eliminates matrix effects. The rate of reaction is measured photometrically at 546nm.

### Reagent Concentration:

<b>R1:</b> Potassium carbonate buffer pH 10.3	250 mmol/l
<b>R1a:</b> Nitrotetrazolium-blue buffer	0.57 mmol/l
Sodium cholate	4.9 mmol/l
Potassium chloride	49 mmol/l
Potassium phosphate	49 mmol/l
Uricase (Arthrobacter species)	> 2.8 kU/l
Detergent	2.1 %

### Preparation and stability:

R1: Add the contents of one bottle 1a careful to one bottle 1. Mix by swirling gently. A slight discoloration of R1 does not interfere with the performance of the assay.

Unopened kit components: up to the expiration date at +2 - +8°C

Store protected from light after opening.

Onboard stability: 28 days.

### Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.  
Disposal of all waste material should be in accordance with local guidelines.

### Specimen:

Collect serum using standard sampling tubes.  
Centrifuge samples containing precipitate before performing the assay.  
Heparinized or EDTA plasma.  
Stability: 3 days at +20°C to +25°C  
2 weeks at + 4°C to + 8°C  
2 month at -20°C

Avoid repeated freezing and thawing. Mix samples well after thawing.

### Limitations - interference:

Criterion: Recovery within ± 10% of initial values.  
Icterus: Significant positive interference at an I index of > 5 (approximate conjugated and unconjugated bilirubin concentration: > 5mg/dl).  
Hemolysis: Significant positive interference at an H index of > 500 (approximate hemoglobin concentration > 500mg/dl).  
Lipemia (Intralipid): No significant interference up to an L index of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.  
Glycemia: No significant interference up to glucose concentrations of 900 mg/dl (50 mmol/l) and therefore fasting samples are not required. Glucose levels of more than 900 mg/dl falsely elevate fructosamine values to a significant extent in samples having normal fructosamine values.  
Ascorbic acid concentrations of up to 4 mg/dl (220 µmol/l), exceeded only after infusion of high doses, do not interfere with the assay.  
No significant interference from uric acid concentrations up to 20 mg/dl (1200 µmol/l).  
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

- Working solutions as described above
- Additional materials required
- Calibrators and controls
- 0.9% NaCl
- General laboratory equipment

### Calculation:

Hitachi systems automatically calculate the fructosamine concentration of each sample.

### Measuring / reportable range:

10-1000 µmol/l

Determine samples with higher activities via the rerun function. On instruments without rerun function, manually dilute with 0.9% NaCl or distilled/deionized water (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2).

### Expected values:

Fructosamine concentrations were determined in 555 apparently healthy subjects between the ages of 20 and 60. A reference range of 205 to 285 µmol/l for adults without diabetes was determined in this study.  
In a poorly controlled diabetic patient population, mean fructosamine values were reported to be 396 µmol/l (range 228-563 µmol/l). A fructosamine concentration above the established expected values is an indicator for hyperglycemia during the preceding 1-3 weeks or longer.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, fructosamine results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Note:

In hydraemic states (e. g. during pregnancy) it is recommended to relate fructosamine to total protein using the following formula:

$$\text{Fructosamine}_{\text{related to total protein}} = \frac{\text{fructosamine}_{\text{measured in } \mu\text{mol/l}} \times 7.2}{\text{total protein}_{\text{measured in g/l}}} [\mu\text{mol/l}]$$

Correction for serum albumin is not recommended. Dysproteinemic states may produce erroneous fructosamine values.

### Analytical sensitivity (lower detection limit)

10 µmol/l

The lower detection limit represents the lowest measurable fructosamine activity that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using samples in an internal protocol. The following results were obtained.

	With in run		
Sample	Mean (µmol/l)	SD (µmol/l)	CV %
Human serum	288	2.58	0.9
Control 1	272	1.88	0.7
Control 2	512	4.12	0.8
	Between day		
Sample	Mean (µmol/l)	SD (µmol/l)	CV %
Human serum	296	8.69	2.9
Control 1	273	3.89	1.4
Control 2	521	9.01	1.7

### **Method comparison:**

A comparison of the Analyticon FRUC (y) with a commercial obtainable assay (x) resulted in the following correlation ( $\mu\text{mol/l}$ ):  
 $y = -8.171 + 1.019 x$ ;  $r = 0.996$

### **Quality Control:**

For quality control, use suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### **Calibration:**

Standardization: the fructosamine assay was standardized against glycated poly-L-lysine and C-glucose.

S1: 0.9% NaCl

S2: Bio Cal® FRUC                      3 x 1 ml                      #14560

### **Calibration frequency**

Two-point calibration is recommended:

- every 7 days if reagent bottles are onboard the analyzer for more than 7 days
- after reagent bottle change if previous reagent bottles were onboard the analyzer for more than 7 days
- after reagent lot change
- as required following quality control procedures

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