Order information:

Catalog No.		Contents					
4551		R1	1 x	100ml	R2	1 x	20 ml
		R4	1 x	5 ml			
H5551	Hit I (ILab*)	R1	6 x	47 ml	R2	6 x	11 ml
H5553	Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	14 ml
AU5553	AU	R1	6 x	60 ml	R2	6 x	14 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 661
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

In vitro test for the quantitative determination of iron in human serum and plasma.

Summary:

Ingested iron is mainly absorbed in the form of Fe²⁺ in the duodenum and upper jejunum. The trivalent form and the heme-bound Fe²⁺-component of iron in food has to be reduced by vitamin C. About 1 mg of iron is assimilated daily. Upon reaching the mucosal cells, Fe²⁺ ions become bound to transport substances. Before passing into the plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bound to trans-ferrin in this form. The transport of Fe ions in blood plasma takes place via trans-ferrin-iron complexes. A maximum of 2 Fe³⁺ ions per protein molecule can be trans-ported. Serum iron is almost completely bound to transferrin. Iron (non-heme) measurements are used in the diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis (a disease associated with widespread deposit in the tissue of the two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin), and chronic renal disease. Iron determinations are performed for the diagnosis and monitoring of microcytic anemia (e.g. due to iron metabolism disorders and hemoglobinopathy), macrocytic anemia (e.g. due to vitamin B12 deficiency, folic acid deficiency and drug-induced metabolic disorders of unknown origin) as well as normocytic anemias such as renal anemia (erythropoietin deficiency), hemolytic anemia, hemoglobinopathy, bone marrow disease and toxic bone marrow damage.


Numerous photometric methods have been described for the determination of iron. All have the following in common:

- Liberation of Fe³⁺ ions from the transferrin complex using acids or detergents.
- Reduction of Fe³⁺ ions to Fe²⁺ ions.
- Reaction of the Fe²⁺ ions to give a colored complex.

The method described here is based on the FerroZine® method without deproteinization.

Test principle:

Colorimetric assay

Under acidic conditions, iron is liberated from transferrin. Lipemic samples are clarified by the detergent. Ascorbate reduces the released Fe³⁺ ions to Fe²⁺ ions which react with FerroZine to form a colored complex. The color intensity is directly proportional to the iron concentration and can be measured photometrically.

Reagent Concentration:

R1:	
Citric acid pH 1.8	200 mmol/l
Thiourea	115 mmol/l
Na-ascorbate	150 mmol/l
detergent	< 1 %
R2:	
Citric acid pH 6.4	200 mmol/l
FerroZine®	6 mmol/l
Stabilizers/ detergents	< 1 %
R4: (#4551)	
Iron	166 µg/dl

Preparation and stability:

R1: Ready for use
 R2: Ready for use
 R4: Ready for use
 Unopened kit components: Up to the expiration date at +2°C to +8°C.
 Onboard stability:
 R1: 28 days at +2°C to +8°C.
 R2: 28 days at +2°C to +8°C.
 Store protected from light.

Specimen:

Collect serum using standard sampling tubes. Heparinized plasma.

Stability:	7 days	at +20°C to +25°C
	3 weeks	at +4°C to +8°C
	several years	at -20°C

Separate serum or plasma from the clot or cells within 1 hour. EDTA and oxalate plasma cause decreased values.
 Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial values.

Icterus: No significant interference up to an index I of 90 (approximate conjugated and unconjugated bilirubin concentration: 90 mg/dl). Hemolysis: No significant interference up to an index H of 70 (approximate hemoglobin concentration: 70 mg/dl). Higher hemoglobin concentrations lead to false-positive values due to contamination of the sample with hemoglobin-bound iron.

Lipemia (Intralipid): No significant interference up to an index L of 500 (approximate triglycerides concentration: 1000 mg/dl). There is poor correlation between turbidity and triglycerides concentration. In patients treated with Deferoxamine, the drug-bound serum iron does not react in the test, resulting in falsely lowered values.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual Testing Procedure:

Wavelength:	Hg 578 nm (560nm)	
Reaction temperature:	+37°C	
Cuvette:	1 cm light path	
Zero adjustment:	One reagent blank for each series	

	Blank	Sample/ Calibrator
R1	1000 µl	1000 µl
Aqua dest.	100 µl	---
Sample/ Calibrator	---	100 µl

Mix and incubate 5 minutes. Read the absorbance A₁. Then add:

	Blank	Sample/ Calibrator
R2	200 µl	200 µl

Mix and incubate 5 minutes. Read the absorbance A₂

Calculation:

$$\frac{A2 \text{ Sample} - A1 \text{ Sample}}{A2 \text{ Calibrator} - A1 \text{ Calibrator}} \times \text{Calibrator conc.} = \text{Iron in } \mu\text{g/dl}$$

SI Units: (µg/dl) x 0.1791 = µmol/l

Measuring range:

5 -1000 µg/dl (0.895-179 µmol/l)

Determine samples with higher concentrations µg/dl via the rerun function. On instruments without rerun function, manually dilute these samples with 0.9% NaCl (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. 2).

Reference value:

Serum/plasma
 Men: 59-158 µg/dl (10.6-28.3 µmol/l)
 Women: 37-145 µg/dl (6.6-26.0 µmol/l)

The concentration of iron in serum/plasma is dependent on the diet and is subject to circadian variations.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the iron results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 1.5 µg/dl (0.269 µmol/l)

The lower detection limit represents the lowest measurable iron concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls within run (n = 20). The following results were obtained:

Within run			
Sample	Mean µmol/l	SD µmol/l	CV %
Sample 1	120.97	2.36	1.95
Sample 2	183.32	1.75	0.95
Sample 3	191.58	3.02	1.58

Reproducibility was determined using controls in an internal protocol between day (n = 10). The following results were obtained:

Between day			
Sample	Mean µmol/l	SD µmol/l	CV %
Sample 1	112.82	1.72	1.52
Sample 2	178.01	4.63	2.60
Sample 3	190.22	3.50	1.84

Method comparison:

A comparison of the Analyticon Fluitest Fe-FZ (y) with a commercial obtainable assay (x) gave with 78 samples the following result (µmol/l):

$$y = 1.001x - 0.493; \quad r = 0.997$$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Multicalibration

S1: 0.9% NaCl		
S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420

Calibration frequency:

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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