BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B5531	200 / 600	R1	6 x	47 ml
		R2	6 x	13 ml
B5553	300 / 600*	R1	6 x	20 ml
		R2	6 x	6 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:


In vitro test for the quantitative determination of iron in human serum and plasma.

Summary:

Ingested iron is mainly absorbed in the form of Fe²⁺ in the duodenum and upper jejunum. The trivalent form and the heme-bound Fe²⁺-component of iron in food has to be reduced by vitamin C. About 1 mg of iron is assimilated daily. Upon reaching the mucosal cells, Fe²⁺ ions become bound to transport substances. Before passing into the plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bound to transferrin in this form. The transport of Fe ions in blood plasma takes place via transferrin-iron complexes. A maximum of 2 Fe³⁺ ions per protein molecule can be transported. Serum iron is almost completely bound to transferrin. Iron (non-heme) measurements are used in the diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis (a disease associated with widespread deposit in the tissue of the two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin), and chronic renal disease. Iron determinations are performed for the diagnosis and monitoring of microcytic anemia (e.g. due to iron metabolism disorders and hemoglobinopathy), macrocytic anemia (e.g. due to vitamin B12 deficiency, folic acid deficiency and drug-induced metabolic disorders of unknown origin) as well as normocytic anemias such as renal anemia (erythropoietin deficiency), hemolytic anemia, hemoglobinopathy, bone marrow disease and toxic bone marrow damage. Numerous photometric methods have been described for the determination of iron. All have the following in common:
 - Liberation of Fe³⁺ ions from the transferrin complex using acids or detergents.
 - Reduction of Fe³⁺ ions to Fe²⁺ ions.
 - Reaction of the Fe²⁺ ions to give a coloured complex.
 The method described here is based on the FerroZine® method without deproteinization.

Test principle:

Colorimetric assay

Under acidic conditions, iron is liberated from transferrin. Lipemic samples are clarified by the detergent. Ascorbate reduces the released Fe³⁺ ions to Fe²⁺ ions which react with FerroZine to form a coloured complex. The colour intensity is directly proportional to the iron concentration and can be measured photometrically.

Reagent Concentration:

R1:	
Citric acid pH 1.8	200 mmol/l
Thiourea	115 mmol/l
Na-ascorbate	150 mmol/l
detergent	< 1 %
R2:	
Citric acid pH 6.4	200 mmol/l
FerroZine®	6 mmol/l
Stabilizers/ detergents	< 1 %

Preparation and stability:

R1: Ready for use
 R2: Ready for use
 Unopened kit components: Up to the expiration date at +2°C to +8°C.
 Onboard stability: R1: 60 days
 R2: 60 days
 Store protected from light.

Specimen:

Collect serum using standard sampling tubes. Heparinized plasma.
 Stability: 7 days at +20°C to +25°C
 3 weeks at +4°C to +8°C
 several years at -20°C

Separate serum or plasma from the clot or cells within 1 hour. EDTA and oxalate plasma cause decreased values.
 Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
 Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial values.
 Icterus: No significant interference up to an index I of 100 (approximate conjugated and unconjugated bilirubin concentration: 100 mg/dl). Hemolysis: No significant interference up to an index H of 75 (approximate hemoglobin concentration: 75 mg/dl). Higher hemoglobin concentrations lead to false-positive values due to contamination of the sample with hemoglobin-bound iron.
 Lipemia (Intralipid): No significant interference up to an index L of 300 (approximate triglycerides concentration: 600 mg/dl). There is poor correlation between turbidity and triglycerides concentration. In patients treated with Deferoxamine, the drug-bound serum iron does not react in the test, resulting in falsely lowered values.
 The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided
 • Working solutions as described above
Additional materials required
 • Calibrators and controls as indicated below
 • 0.9% NaCl

Measuring range:

35 - 1000 µg/dl (6.3 - 179 µmol/l)
 Determine samples with higher concentrations via the rerun function using 0.9% NaCl as diluents.

Reference value:

Serum/plasma
 Men: 59-158 µg/dl (10.6-28.3 µmol/l)
 Women: 37-145 µg/dl (6.6-26.0 µmol/l)
 The concentration of iron in serum/plasma is dependent on the diet and is subject to circadian variations.
 Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the iron results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 1.15 µg/dl (0.21 µmol/l)
 The lower detection limit represents the lowest measurable iron concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls within run (n = 20). The following results were obtained:

Sample	Within run		
	MW µmol/l	SD µmol/l	CV %
Sample 1	111.31	3.439	3.1
Sample 2	180.01	1.018	0.6
Sample 3	190.71	2.247	1.2

Reproducibility was determined using controls in an internal protocol between day (n = 10). The following results were obtained:

Sample	Run to run		
	MW µmol/l	SD µmol/l	CV %
Sample 1	116.06	5.171	4.5
Sample 2	182.27	2.457	1.3
Sample 3	202.64	2.043	1.0

Method comparison:

A comparison of the Analyticon Fluitest Fe-FZ (y) with a commercial obtainable assay (x) gave with 65 samples the following result (µmol/l):
 $y = 1.1503x + 2.5992$; $r = 0.9941$

Quality Control:

Human Control Serum		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Fluitest® FE FZ

IRON

Calibration:

Multicalibration

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E

10 x 3 ml

#1430

Calibration frequency:

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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