Fluitest® GLDH

GLUTAMATE DEHYDROGENASE

Order information:

Catalog No.		Contents			
H7401	Hit I (ILab*)	R1	5 x 20 ml	R2	1 x 20 ml
AU7401	AU	R1	5 x 20 ml	R2	1 x 20 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 361-400 (user defined method)
 Hitachi 917: ACN 588
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Test for the quantitative in vitro determination of glutamate dehydrogenase (GLDH) in human serum und plasma with clinical chemistry analyzers.

Summary:

Glutamate dehydrogenase (GLDH) is a mitochondrial enzyme which is present in many tissues. Significant elevations of GLDH activity are measured in necrosis of hepatocytes, as in acute toxic liver necrosis and in hypoxic liver diseases. The measurement of GLDH is used to evaluate the extent of parenchymal liver damage and in conjunction with the transaminases ALAT/GPT and ASAT/GOT in the differential diagnosis of liver disorders. The calculation of the (ALAT+ASAT)/GLDH ratio enables to differentiate between inflammatory liver diseases and liver necrosis due to intoxication or ischemia.

Test principle:

Optimised UV-Test according to DGKC
 $2\text{-Oxoglutarate} + \text{NADH} + \text{NH}_4^+ \xrightarrow{\text{GLDH}} \text{L-Glutamate} + \text{NAD}^+ + \text{H}_2\text{O}$

Reagents – contents and concentrations:

R1:	
Triethanolamine pH 8.0	75 mmol/l
α -Ketoglutarate	10 mmol/l
Ammoniumacetate	150 mmol/l
EDTA	3.75 mmol/l
ADP	1.5 mmol/l
LDH	> 2.3 kU/l
R2:	
NADH	1.3 mmol/l

Preparation and stability:

Reagents are stable up to the indicated expiry stored at +2°C to +8°C. Protect from light. Avoid contamination!
 R1: ready to use
 R2: ready to use

Specimen:

Serum, Heparin or EDTA plasma.
 Stability: 7 days at +20°C to +25°C
 7 days at +2°C to +8°C
 4 weeks at -20°C

Centrifuge samples containing precipitates before performing the assay.

Note:

For in vitro diagnostic use. The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.
 Ascorbic acid: no interference up to 30 mg/dl
 Icterus: no interference up to 60 mg/dl
 Hemolysis: no interference up to 500 mg/dl.
 Lipids interfere. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Measuring/reportable range:

2-120 U/l
 Determine samples having higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl solution (e.g. 1 + 5). Multiply the result by the appropriate dilution factor (e.g. 6).

Expected values:

Women < 5,0 U/L
 Men < 7,0 U/L
 Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, GLDH results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

2 U/l
 The lower detection limit represents the lowest measurable GLDH activity that can be distinguished from zero.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Reagents as described above
- Additional materials required
- Calibrators and controls as indicated below.
- 0.9% NaCl solution
- General laboratory equipment

Manual procedure:	
Wavelength:	340nm (Hg 334nm)
Temperature:	+37°C
Cuvette:	1 cm light path
Zero adjustment:	Air or Aqua dest.
Serum/Plasma	
R1	1000 μ l
Sample	120 μ l
Mix and incubate for approx. 3min. Then add:	
R2	200 μ l
Mix, read initial absorbance after 30sec. Start stopwatch simultaneously. Read absorbance again after 1, 2 und 3min. Calculate Δ A/min.	
Calculation:	
<i>With factor:</i>	
340 nm	Δ A/min. x 1111 = GLDH activity (U/l) ^a
Hg 334 nm	Δ A/min. x 1133 = GLDH activity (U/l)
<i>With calibrator:</i>	
$\frac{\Delta$ A/min Sample}{\DeltaA/min calibrator}	X Calibrator conc. = GLDH activity (U/l)

Imprecision:

Reproducibility according to an internal protocol gave the following results:

Sample	Between Day		
	MW U/l	SD U/l	% VK
Sample 1	5.77	0.51	6.98
Sample 2	18.3	0.39	2.11
Sample 3	32.0	0.78	2.43

Sample	Within run		
	MW U/l	SD U/l	% VK
Sample 1	6.18	0.43	8.78
Sample 2	16.1	0.49	3.02
Sample 3	33.2	0.80	2.40

Method comparison:

A comparison of the Analyticon Fluitest® GLDH (y) with a commercial obtainable assay (x) gave the following result (U/l):
 $y = 1.034x + 0.006$; $r = 0.999$

Quality control:

Human control serum		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: NaCl (0.9%)		
S2: Bio Cal® E	20 x 3 ml	#1430

Calibration frequency:

- A two-point calibration is recommended:
- after bottle change
 - quality control requirements

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