### BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B5831	200 / 600	R1	6 x	47 ml
		R2	6 x	13 ml
B5833	300 / 600*	R1	6 x	54 ml
		R2	6 x	14 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

### Intended use:


Enzymatic in vitro test for the quantitative regulation (determination) of glucose in human serum and plasma.

### Summary:

Carbohydrates supply the body with glucose. Glucose is the most important monosaccharide in the blood; the postprandial concentration is 5mmol/L of glucose. Glucose substrate is an indispensable energy supplier which supports cellular function. Glucose degradation occurs in glycolysis. Glucose measurements are used in the diagnosis and monitoring of carbohydrate metabolism disorders including diabetes mellitus, neonatal hypoglycaemia, idiopathic hypoglycaemia, and pancreatic islet cell carcinoma. The hexokinase method, based on the work of Schmidt, Peterson and Young, is a recognized reference method.

### Test principle:

Determination of glucose concentration according to the following reaction using hexokinase and glucose-6-phosphate-dehydrogenase:

### Reagent concentration:

<b>R1:</b>	
Tris buffer, pH 7.8	100 mmol/l
Magnesiumacetate	4 mmol/l
NADP	1.2 mmol/l
ATP	1.2 mmol/l
<b>R2:</b>	
HEPES buffer, pH 7.0	30 mmol/l
Magnesiumacetate	4 mmol/l
Hexokinase	> 900 U/l
G-6-P-DH	> 1800 U/l

### Preparation and stability:

R1: Ready for use  
R2: Ready for use

Unopened kit components: Up to the expiration date at +2°C to +8°C.

Onboard stability: R1: 60 days  
R2: 60 days

### Specimen:

Blood, serum, plasma, urine and CSF  
Cellular constituents from sera or plasma have to be separated within 1/2 h after collecting blood specimen. By adding an inhibitor for glycolysis (NaF, KF) samples can be stored

up to 24h at +15°C to 25°C or  
up to 7 days at +2°C to +4°C

Liquor, urine

up to 24h at +15°C to 25°C or  
up to 7 days at +2°C to +4°C

Centrifuge samples containing precipitate before performing the assay.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.  
Icterus: No significant interference up to an index I of 100 (approximate conjugated bilirubin: 100 mg/dl)  
Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/dl).  
Lipemia (Intralipid): No significant interference up to an index L of 2100 (approximate triglycerides concentration: 4200 mg/dl) There is poor correlation between turbidity and triglycerides concentration.  
Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.  
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

**Note:** Glucose values in the external reference material which were determined in a method comparison against GOD/oxygen-electrode, show on average an approx. 3% higher abnormality.

### Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Testing procedure:

Materials provided  
• Working solutions as described above  
*Additional materials required*  
• Calibrators and controls as indicated below  
• 0.9% NaCl

### Measuring /reportable range:

2.5 - 500 mg/dl or 0.14 - 27.8 mmol/l  
Determine samples with higher concentrations via the rerun function using 0.9% NaCl as diluent.

### Expected values:

Thomas, Krieg and Colombo:  
**Serum/plasma**  
55 - 115 mg/dl bzw. 3.05 - 6.38 mmol/l  
**24 hour urine:**  
20 - 90 mg/24h or 0.11 - 0.50 mmol/24h  
1.3 - 6.0 mg/dl or 0.07 - 0.33 mmol/l  
**Morning urine:**  
6 - 20 mg/dl bzw. 0.3 - 1.1 mmol/l

Expected values according to Tietz:

**Liquor**  
Children: 60 - 80 mg/dl or 3.33 - 4.44 mmol/l  
Adult 40 - 70 mg/dl or 2.22 - 3.89 mmol/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the glucose results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 0.22 mg/dl  
The lower detection limit represents the lowest measurable glucose concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	57.75	0.593	1.0
Sample 2	146.99	2.296	1.6
Sample 3	107.55	2.073	1.9

Sample	Run to run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	261.08	9.931	3.8
Sample 2	93.90	3.492	3.7
Sample 3	207.35	7.282	3.5

### Method comparison:

A comparison of the Analyticon Fluitest GLU HK (y) with a commercial obtainable assay (x) gave the following result (mg/dl):  
 $y = 0.91x + 4.155$ ;  $r = 0.9940$

### Quality Control:

Human Control Serum:  
Contronorm® Plus 5 x 5 ml #1205  
20 x 5 ml #1220  
ControPath® Plus 5 x 5 ml #1305  
20 x 5 ml #1320  
Human Urine Control:  
Urine Control 8 x 5 ml #1507

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

# Fluitest® GLU HK

GLUCOSE HEXOKINASE 5+1

## **Calibration:**

Standardization: The glucose method was standardized against ID/MS as a reference method.

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Text passages with grey background were changed in the latest revision of this package insert.

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