### Order information:

Catalog No.	Contents					
4981	R1	4 x	100 ml	R2	1 x	5 ml
	R4	1 x	5 ml			

### Intended use:

Enzymatic in vitro test for the quantitative regulation (determination) of glucose in human serum, -plasma, -urine and liquor.


### Summary:

Glucose is the central energy source of the cells in the organism. The most common supply follows hydrolytic cleavage of polymeric carbohydrates, in general starch. Glucose is a monosaccharide with a postprandial concentration of 5 mmol/l in the blood and serves as an indispensable energy-supply for cellular functions. The glucose catabolism takes place via the glycolysis as the first step, followed by the citric acid cycle and oxidative phosphorylation.

Glucose regulations become executed for the diagnosis and course control of carbohydrate metabolism illnesses like the diabetes mellitus, neonatal hypoglycemia, idiopathic hypoglycemia and with insulinoma. The present, on works of Schmidt as well as from Peterson and Young being based Hexokinase-method serves as a recognized reference method.

### Test principle:

Determination of glucose concentration according to the following reaction using hexokinase and glucose-6-phosphate-dehydrogenase:

### Reagent concentration:

<b>R1:</b>	
Bis-Tris buffer, pH 8.0	50 mmol/l
MgCl <sub>2</sub>	3.8 mmol/l
NADP	1.2 mmol/l
ATP	1.2 mmol/l
<b>R2:</b>	
Hexokinase	> 200 KU/l
G-6-P-DH	> 200 KU/l
<b>R4:</b>	
Glucose	100 mg/dl (5,55 mmol/l)

### Preparation and stability:

Mix 100 parts of buffer/R1 with 1 part of enzyme reagent/R2.

R4: Ready for use.

The working solution R2 mixed with R1 is stable for:

7 days	at +2°C to +8°C
1 day	at +20°C to +25°C

### Specimen:

Blood, serum, plasma, urine and CSF

Cellular constituents from sera or plasma have to be separated within 1/2 h after collecting blood specimen. By adding an inhibitor for glycolysis (NaF, KF) samples can be stored

up to 24h	at +15°C to 25°C or
up to 7 days	at +2°C to +4°C

Liquor, urine

up to 24h	at +15°C to 25°C or
up to 7 days	at +2°C to +4°C

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 15

(approximate conjugated bilirubin: 15 mg/dl)

Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 250 (approximate triglycerides concentration: 250 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

**Note:** Glucose values in the external reference material which were determined in a method comparison against GOD/oxygen-electrode, show on average an approx. 3% higher abnormality.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

*Additional materials required*

• Calibrators and controls as indicated below

• 0.9% NaCl

### Manual procedure:

For serum or plasma without deproteinization:

Wavelength: Hg 334 nm, 340 nm, Hg 365 nm

Temperature: +25 / +37°C

Cuvette: 1 cm light path

Zero adjustment: against reagent blank

	sample	reagent blank
Working solution	2000 µl	2000 µl
saline (0.9%NaCl)	---	20 µl
serum / standard	20 µl	---

Mix and incubate 10-15 min. at +20 to +25°C or 5 min. at +37°C. Determine absorbance of sample (A<sub>sample</sub>) and blank (A<sub>blank</sub>) against reagent blank.

### Calculation:

$$A_{\text{sample}} - A_{\text{blank}} = \Delta A_{\text{sample}}$$

By factor:

	<b>mg/dl</b>	<b>mmol/l</b>
Hg 334 nm	$\Delta A_{\text{sample}} \times 291$	$\Delta A_{\text{sample}} \times 16.1$
340 nm	$\Delta A_{\text{sample}} \times 286$	$\Delta A_{\text{sample}} \times 15.8$
Hg 365 nm	$\Delta A_{\text{sample}} \times 515$	$\Delta A_{\text{sample}} \times 28.6$

by standard:

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times \text{Conc. of Standard} = \text{Glucose Conc.}$$

### Measuring /reportable range:

Up to 500 mg /dl or 27.8 mmol/l

At higher concentrations, dilute the sample 1:3 with 0.9 % NaCl. Multiply the result by factor 3.

### Expected values:

Thomas, Krieg and Colombo:

#### Serum/plasma

55 – 115 mg/dl or 3.05 – 6.38 mmol/l

#### 24 hour urine:

20 – 90 mg/24h or 0.11 – 0.50 mmol/24h

1.3 – 6.0 mg/dl or 0.07 – 0.33 mmol/l

#### Morning urine:

6 – 20 mg/dl or 0.3 – 1.1 mmol/l

Expected values according to Tietz:

#### Liquor

Children: 60 – 80 mg/dl or 3.33 – 4.44 mmol/l

Adult: 40 – 70 mg/dl or 2.22 – 3.89 mmol/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the glucose results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 2 mg/dl

The lower detection limit represents the lowest measurable glucose concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	97.9	1.41	1.44
Sample 2	240	2.16	0.90
Sample 3	454	3.67	0.81

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	98.9	2.62	2.64
Sample 2	232	7.30	3.15
Sample 3	473	18.9	4.00

### Method comparison:

A comparison of the Analyticon Fluitest GLU HK (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 0.999x + 1.423; \quad r = 1.008$$

### Quality Control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Urine Control:

Urine Control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Standardization: The glucose method was standardized against ID/MS as a reference method.

Multicalibration

S1: 0.9% NaCl

S2: Bio Cal® E                      10 x 3 ml                      #1430

Bio Cal®                              20 x 3 ml                      #1420

R4: Calibrator provided in kit

### Disposal:

Please note the legal regulations.

### Literature:

1. Bablock W et al. A General Regression Procedure for Method Transformation: J Clin Chem Clin Biochem 1988;26:783-790.
2. Bassing H, Bablok W.A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
3. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986; 24:863-869.
4. Greiling H, Gressner AM ed. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed., Stuttgart / New York; Schattauer Verlag: 1995.
5. Krieg M et al. Vergleichende quantitative Analytik klinisch-chemischer Kenngrößen im 24 Stunden- Urin und Morgenurin. J Clin Chem Clin Biochem 1985;24:863-869.
6. Peterson JI, Young DS. Anal Biochemistry 1985;23:301.
7. Schmidt FH, Klein Wschr 1961;39:1244
8. Thomas L ed. Labor und Diagnose, 4. ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992.
9. Tietz NW ed Clinical Guide to Laboratory Tests, 2<sup>nd</sup> ed Philadelphia, PA: WB Saunders Company, 1990:246-250.
10. Tietz NW ed. Clinical Guide to Laboratory Tests, 3. Auflage. Philadelphia. PA: WB Saunders Company: 1995:266-273.

Text passages with grey background were changed in the latest revision of this package insert.

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