Order information:

Catalog No.	Contents	
4611	R1	4 x 100 ml
	R2	4 x for 100 ml
458	R4	1 x 5 ml
	R2	4 x for 250 ml
	R1	4 x 250 ml
	R4	1 x 10 ml

Intended use:

Enzyme in vitro test for the quantitative determination of glucose in human serum or plasma.

Summary:


Glucose is the central energy source of the cells in the organism. The most common supply follows hydrolytic cleavage of polymeric carbohydrates, in general starch. Glucose is a monosaccharide with a postprandial concentration of 5 mmol/l in the blood and serves as an indispensable energy-supply for cellular functions. The glucose catabolism takes place via the glycolysis as the first step, followed by the citric acid cycle and oxidative phosphorylation.

Glucose regulations become executed in the diagnosis and course control of carbohydrate metabolism illnesses like the diabetes mellitus, neonatal hypoglycemia, idiopathic hypoglycemia and with insulinoma.

The test bases on the coupling of the enzymatic oxidation of glucose by glucose oxidase resulting in hydrogen peroxide, which is subsequently used for the generation of a coloured product by peroxidase. In the Trinder method the carcinogenic ortho-dianisidine used in earlier formulations has been replaced by phenole and 4-aminoantipyrine.

Test principle:

Enzymatic colorimetric test on basis of Trinder-Reaction:

Reagent Concentration:

R1:	
Phosphate buffer pH 7.5	150 mmol/l
Phenole	7.5 mmol/l
R2:	
GOD	12000 U/l
POD	660 U/l
4-aminoantipyrine	0.40 mmol/l
R4:	
Glucose	100 mg/dl (5.55 mmol/l)

Preparation and stability:

Dilute contents of enzyme reagent /R2 with the corresponding volume of buffer/R1.

The working solution is stable for:

3 months at	+2°C to +8°C or
4 weeks at	+20°C to +25°C

Keep protected from light!

R4: Ready for use.

The standard is refrigerated stable up to the expiry date.

Specimen:

Collect serum, heparinized plasma or EDTA plasma using standard sampling tubes. The separation of the cells of the blood test should take place not later than half an hour after the decline. Non-hemolyzed serum or plasma can be stored refrigerated for up to 12 hours before the determination.

Tests which contain precipitates must be centrifuged before performing the test. By adding of NaF or KF the specimen can be stored

24 hours	at +15°C to +25°C or
7 days	at +4°C.

Blood samples should only be drawn prior to the administration of Metaminazole.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 40 (approximate conjugated bilirubin: 40mg/dl)

Hemolysis: No significant interference up to an index H of 600 (approximate haemoglobin concentration: 600 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 80 (approximate triglycerides concentration: 160 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminazole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:			
Wavelength:	Hg 546 nm (492-550 nm)		
Temperature:	+25°C / +30°C / +37°C		
Cuvette:	1 cm light path		
Zero adjustment:	against reagent blank		
	Blank	Standard	Sample
Working reagent	1000 µl	1000 µl	1000 µl
Standard R4	----	10 µl	----
Sample	----	----	10 µl
Mix, measure (A) after incubating 15 minutes at +37°C for or 30 minutes at +25°C.			
Within 60 minutes read absorbance of sample and standard against reagent blank.			
Determine the absorbance change as $\Delta A \text{ sample} = (A \text{ sample} - A \text{ blank})$ $\Delta A \text{ standard} = (A \text{ standard} - A \text{ blank})$ and use this for the calculation.			
Calculation:			
$\frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times \text{standard conc.} = \text{Glucose conc. (mg/dl)}$			

Measuring /reportable range:

Up to 400 mg/dl

At higher concentrations, dilute the sample with 0.9 % NaCl (e.g. 1+1).

Multiply the result by the appropriate dilution factor (e.g. 2).

Expected values:

Thomas, Krieg and Colombo

Serum/plasma

55 – 115 mg/dl or 3.05 – 6.38 mmol/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the glucose results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 2 mg/dl

The lower detection limit represents the lowest measurable glucose concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	95.0	1.18	1.24
Sample 2	231	2.06	0.89
Sample 3	486	5.11	1.05

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	96.5	2.24	2.32
Sample 2	223	5.19	2.33
Sample 3	485	7.24	1.49

Method comparison:

A comparison of the Analyticon GLU (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 0.984 x - 0.561;$$

$$r = 0.999$$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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