Order information:

Catalog No.	Contents			
4341	R1	4 x	100 ml	R4 1 x 5 ml
4342	R1	4 x	250 ml	R4 1 x 5 ml
4347	R1	4 x	500 ml	R4 1 x 5 ml
H5701	Hit I (ILab*)	R1	9 x	100 ml
H5703	Hit 917 (AU*)	R1	12 x	60 ml
AU5703	AU	R1	12 x	60 ml
5799	Deproteinization-solution	R1	1 x	500 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 249
Hitachi 917: ACN 525

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Enzyme in vitro test for the quantitative regulation (determination) of glucose in human serum, -plasma.

Summary:



Glucose is the central energy source of the cells in the organism. The most common supply follows hydrolytic cleavage of polymeric carbohydrates, in general starch. Glucose is a monosaccharide with a postprandial concentration of 5 mmol/l in the blood and serves as an indispensable energy-supply for cellular functions. The glucose catabolism takes place via the glycolysis as the first step, followed by the citric acid cycle and oxidative phosphorylation.

Glucose regulations become executed the diagnosis and course control of carbohydrate metabolism illnesses like the diabetes mellitus, neonatal hypoglycemia, idiopathic hypoglycemia and with insulinoma.

The test bases on the coupling of the enzymatic oxidation of glucose by glucose oxidase resulting in hydrogen peroxide, which is subsequently used for the generation of a coloured product by peroxidase. In the Trinder method the carcinogenic ortho-dianisidine used in earlier formulations has been replaced by phenole and 4-amino-antipyrine.

Test principle:

Enzymatic colorimetric test on basis of Trinder – Reaction:

Reagent concentration:

R1:	
Phosphate buffer, pH 7.5	0.5 mol/l
Phenol	7.5 mmol/l
GOD	12000 U/l
POD	660 U/l
4-amino-antipyrine	0.40 mmol/l
R4: (#4341, #4342, #4347)	
Glucose	100 mg/dl (5.55 mmol/l)

Preparation and stability:

R1: Ready for use

R4: Ready for use. Single reagent, liquid stable up to the expiry date even after opening:

until expiry date	at +2°C to +8°C
3 weeks	at +20°C to +25°C

Keep protected from light!

Coloration of the reagent (reagent blank at 546 nm, 1 cm > 0.2) indicates a contamination or damage by storage at higher temperatures.

Onboard stability:	R1	28 days
	R4	28 days

Specimen:

Capillary blood*, serum, Heparin-or EDTA-plasma.

The separation of the cells of the blood test should take place to not later than one half an hour after the decline. Non-hemolyzed serum or plasma can be stored refrigerated for up to 12 hours before the determination.

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metaminazole.

*Capillary blood has to be liberated from protein (deproteinization)

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 55 (approximate conjugated bilirubin: 55 mg/dl)

Hemolysis: No significant interference up to an index H of 450 (approximate haemoglobin concentration: 450 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 2000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminazole (Dipyron). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

• for procedure with deproteinization: #5799 deproteinization solution

Manual procedure for kinetic:

Wavelength:	Hg 546 nm (492 – 550 nm)
Temperature:	+37°C
Cuvette:	1 cm light path
Zero adjustment:	against reagent blank

	Blank	w/o deproteinization	With deproteinization
R1	1000 µl	1000 µl	1000 µl
Calibrator/sample	-----	10 µl	
Deproteinized	-----	-----	100 µl
Sample/Calibrator			

Mix, and incubate 90 seconds at +37°C, then read absorbance (A₁) and start stopwatch at the same time. Repeat readings after exactly 5 minutes (A₂).

Determine the absorbance change as:

$$\Delta A \text{ sample} = [A_2 (\text{sample}) - A_1 (\text{sample})]$$

$$\Delta A \text{ standard} = [A_2 (\text{standard}) - A_1 (\text{standard})]$$

and use this for the calculation.

Manual procedure for endpoint:

Wavelength:	Hg 546 nm (492 – 550 nm)
Temperature:	+37°C
Cuvette:	1 cm light path
Zero adjustment:	against reagent blank

	Blank	w/o deproteinization	With deproteinization
R1	1000 µl	1000 µl	1000 µl
Calibrator/sample	-----	10 µl	
Deproteinized	-----	-----	100 µl
Sample/Calibrator			

Mix, measure (A) after incubating 15 minutes at +37°C or 30 minutes at +25°C.

Within 60 minutes read absorbance of sample and standard against reagent blank. Determine the absorbance change as:

$$\Delta A \text{ sample} = (A \text{ sample} - A \text{ blank})$$

$$\Delta A \text{ standard} = (A \text{ standard} - A \text{ blank})$$

and use this for the calculation.

Calculation:

Without deproteinization:

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Callibrator}}} \times \text{Calibrator conc.} = \text{Glucose conc.}$$

With deproteinization:

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Callibrator}}} \times \text{Calibrator conc.} \times 10 = \text{Glucose conc.}$$

Note: In assays with deproteinization the calibrator has to be pre-treated like samples.

Measuring /reportable range:

The endpoint method is linear up to 400 mg /dl (22.2 mmol/l)
 The kinetic method up to 700 mg /dl (38.9 mmol/l).
 At higher concentrations, dilute the sample 1 + 1 with 0.9 % NaCl. Multiply the result by factor 2.
 Alternatively, using an automatic analyzer, determine again with rerun-function.

Expected values:

Serum, Plasma: 75 – 115 mg/dl (4.1 – 6.4 mmol/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the glucose results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 2 mg/dl

The lower detection limit represents the lowest measurable glucose concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls between day. The following results were obtained:

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	92.5	1.74	1.88
Sample 2	221	3.94	1.78
Sample 3	483	7.04	1.46

Within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	95.7	0.43	0.45
Sample 2	235	1.28	0.55
Sample 3	501	3.34	0.67

Method comparison:

A comparison of the Analyticon Fluitest Glucose GOD-PAP (y) with a commercial obtainable assay (x) gave the following result (mg/dl):
 $y = 0.990 x - 1.001$; $r = 0.999$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl		
S2: Bio Cal® E	10 x 3 ml	#1430
BioCal®	20 x 3 ml	#1420
R4: Glucose standard provided in kit		

Calibration frequency

Two-point calibration is recommended:

- every 14 days if the reagent bottles are onboard the analyzer for more than 14 days.
- after reagent bottle change if the previous reagent bottles were onboard the analyzer for more than 14 days.
- after lot change
- as required following quality control procedures

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