Fluitest® GOT AST

ASPARTATE AMINOTRANSFERASE

BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B7001	200 / 600	R1	6 x	47 ml
		R2	6 x	11 ml
B7033	300 / 600*	R1	6 x	57 ml
		R2	6 x	14 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Fluitest GOT AST without pyridoxal phosphate activation

Intended use:

In vitro test for the quantitative determination of aspartate amino-transferase (AST) in human serum and plasma.


Summary:

Aspartate aminotransferase (glutamate-oxaloacetate-transaminase) belongs to the transaminases, which catalyze the interconversion of amino acids and α -ketoacids by transfer of amino groups. Aspartate aminotransferase is commonly found in human tissue. Although heart muscle is found to have the most activity of the enzyme, significant activity has also been seen in the brain, liver, gastric mucosa, adipose tissue, skeletal muscle, and kidneys. AST is present in both the cytoplasm and mitochondria of cells. In cases involving mild tissue injury, the predominant form of AST is that from the cytoplasm, with a smaller amount coming from the mitochondria. Severe tissue damage results in more of the mitochondrial enzyme being released. Elevated levels of the transaminases can signal myocardial infarction, hepatic disease, muscular dystrophy, and organ damage. In 1955, Karmen et al described the first kinetic determination of AST activity in serum. The International Federation of Clinical Chemistry (IFCC) recommended in 1977 and 1980 standardized procedures for AST determination, including optimization of substrate concentrations, employment of TRIS* buffers, preincubation of combined buffer and serum to allow side reactions with NADH to occur, substrate start, and optional pyridoxal phosphate activation. In 2002 the IFCC confirmed their recommendation and extended it to 37°C. This method is derived from the IFCC reference method.


*TRIS = Tris(hydroxymethyl)-aminomethane

Test principle:

UV test according to a standardized method:

AST is the enzyme which catalyzes this equilibrium reaction. The oxalacetate increase is measured in a subsequent indicator reaction which is catalyzed by malate dehydrogenase.

In the second reaction, NADH is oxidized to NAD. The rate of decrease in NADH (Measured photometrically) is directly proportional to the rate of formation of oxalacetate, and thus the AST activity.

Reagent Concentration:

R1:	
Tris buffer pH 7.8	100 mmol/l
L-Aspartate	200 mmol/l
LDH	800 U/l
MDH	600 U/l
R2:	
NADH ₂	0.18 mmol/l
2-Oxoglutarate	12 mmol/l

Preparation and stability:

Unopened kit components: Up to the expiration date at +2°C to +8°C

R1: Ready for use.
R2: Ready for use.

Onboard stability: R1 60 days
R2 60 days

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma.

Stability: 24 hours at +20°C to +25°C
7 days at +2°C to +8°C

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at +2°C to +8°C.

Centrifuge samples containing precipitate before performing the assay.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.

Icterus: No significant interference up to an index I of 70 (approximate conjugated and unconjugated bilirubin: 70 mg/dl)

Hemolysis: No significant interference up to an index H of 425 (approximate haemoglobin concentration: 425 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 425 (approximate triglycerides concentration: 850 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Blood samples should only be drawn prior to the administration of Sulfasalazine and Sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Notes:

For in vitro diagnostic use.

Warning! R1 contains hazardous components.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring /reportable range:

11 - 400 U/l or 0.18 - 6.67 μ kat/l.

Determine samples having higher activities via the rerun function using 0.9% NaCl solution as diluents.

Expected values:

According to the IFCC method.

	U/l	μ kat/l
Men	<50	< 0.85
Women	< 35	< 0.60

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the AST results should always be assessed in conjunction with the patients medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 4.5 U/l or 0.075 μ kat/l

The lower detection limit represents the lowest measurable AST concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol between day (n = 20). The following results were obtained:

Run to run			
Sample	Mean (U/l)	SD (U/l)	CV %
sample 1	48.16	1.817	3.8
sample 2	151.22	4.923	3.3
sample 3	139.82	6.827	4.9

Reproducibility was determined using controls within protocol between day (n = 20). The following results were obtained:

Within run			
Sample	Mean (U/l)	SD (U/l)	CV %
sample 1	47.27	1.047	2.2
sample 2	146.91	1.568	1.1
sample 3	101.93	1.204	1.2

Method comparison:

A comparison of the Analyticon Fluitest® AST (y) with a commercial obtainable assay (x) gave the following result with 165 samples:

$y = 0.943x - 0.355$; $r = 0.998$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
ControPath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Fluitest® GOT AST

ASPARTATE AMINOTRANSFERASE

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency:

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedures

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Text passages with grey background were changed in the latest revision of this package insert.

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* TRIS = Tris(hydroxymethyl)-aminomethan

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