GOT AST

ASPARTATE AMINOTRANSFERASE

Order information:

Catalog No.	Contents
11703	R1 1 x 65 ml R2 20 x for 3 ml
117	R1 1 x 110 ml R2 10 x for 10 ml
1174	R1 4 x 100 ml R2 4 x for 100 ml

Intended use:

In vitro test for the quantitative determination of aspartate amino-transferase (AST) in human serum and plasma.

Summary:

Aspartate aminotransferase (glutamate oxaloacetate transaminase) belongs to the transaminases, which catalyze the interconversion of amino acids and α -ketoacids by transfer of amino groups. Aspartate aminotransferase is commonly found in human tissue. Although heart muscle is found to have the most activity of the enzyme, significant activity has also been seen in the brain, liver, gastric mucosa, adipose tissue, skeletal muscle, and kidneys.

AST is present in both the cytoplasm and mitochondria of cells. In cases involving mild tissue injury, the predominant form of AST is that from the cytoplasm, with a smaller amount coming from the mitochondria. Severe tissue damage results in more of the mitochondrial enzyme being released. Elevated levels of the transaminases can signal myocardial infarction, hepatic disease, muscular dystrophy, and organ damage.

In 1955, Karmen et al described the first kinetic determination of AST activity in serum. The International Federation of Clinical Chemistry (IFCC) recommended in 1977 and 1980 standardized procedures for AST determination, including optimization of substrate concentrations, employment of TRIS* buffers, preincubation of combined buffer and serum to allow side reactions with NADH to occur, substrate start, and optional pyridoxal phosphate activation.

This method is derived from the IFCC reference method.

*TRIS = Tris(hydroxymethyl)-aminomethane

Test principle:

UV-Method recommended by the IFCC:

2-Oxoglutarate + L-Aspartate $\xrightarrow{\text{GOT}}$ Glutamate + Oxaloacetate

Oxaloacetate + NADH + H⁺ $\xrightarrow{\text{MDH}}$ Malate + NAD⁺

Reagent concentration:

R1:	
Tris buffer pH 7.8	80 mmol/l
L-Aspartate	200 mmol/l
R2:	
NADH	0.18 mmol/l
LDH	800 U/l
MDH	600 U/l
Oxoglutarate	12 mmol/l

Preparation and stability:

Dissolve enzyme reagent /R2 with the corresponding volume of buffer /R1.

Before first use wait at least 15 min.

This solution is stable: 1 week at +2 to +8°C

Unopened vials are stable: up to the expiration date at +2 to +8°C

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma.

Stability: 24 hours at +20°C to +25°C
7 days at +2°C to +8°C

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at +2°C to +8°C.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

 Warning! R1 contains hazardous components.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within \pm 10% of initial value.

Icterus: No significant interference up to an index I of 70 (approximate conjugated and unconjugated bilirubin: 70 mg/dl)

Hemolysis: No significant interference up to an index H of 900 (approximate haemoglobin concentration: 900 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 450 (approximate triglycerides concentration: 900 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Lipemia may cause absorbance flagging as a result of an absorbance increase.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Manual procedure:

Wavelength: Hg 334 nm, Hg 340 nm or Hg 365 nm
Temperature: +25 / +30 / +37°C
Cuvette: 1 cm light path
Zero adjustment: against air or distilled water

	Macro	Semi	Micro
Working reagent	2500 μ l	1000 μ l	500 μ l
Sample	250 μ l	100 μ l	50 μ l

Mix, incubate for 1 min. at assay temperature and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes.

If the $\Delta A/\text{min}$. is between 0.11 and 0.25 at Hg 334 nm/ 340 nm or 0.06 and 0.130 at Hg 365 nm use only the values for the first 2 minutes for the calculation.

Calculation:

	Macro	Semi	Micro
Hg 340 nm $\Delta A/\text{min} \times$	1750	1750	1750
Hg 334 nm $\Delta A/\text{min} \times$	1790	1790	1790
Hg 365 nm $\Delta A/\text{min} \times$	3240	3240	3240

Measuring /reportable range:

Up to 440 U/l

$\Delta E/\text{min}$ 0.250 at 340 nm or $\Delta E/\text{min}$ 0.130 at 365 nm.

At higher activities dilute the sample with 0.9% NaCl or distilled/deionized water (e.g. 1+9). Multiply the result by the appropriate dilution factor (e.g. factor 10).

Expected values:

According to the IFCC method.

	25°C	30°C	37°C
Women	up to 16 U/l	up to 22 U/l	up to 35 U/l
Men	up to 19 U/l	up to 26 U/l	up to 50 U/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the AST results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Imprecision:

Reproducibility was determined using controls (n = 20). The following results were obtained:

Between day			
Sample	Mean U/l	SD U/l	CV %
sample 1	33.2	0.67	2.02
sample 2	82.9	1.18	1.43
sample 3	148	2.06	1.40

Within run			
Sample	MW U/l	SD U/l	CV %
sample 1	33.4	0.76	2.29
sample 2	116	1.56	1.34

Analytical sensitivity (lower detection limit)

Detection limit: 4 U/l or 0.07 μ kat/l

The lower detection limit represents the lowest measurable AST concentration that can be distinguished from zero.

Method comparison:

A comparison of the Analyticon AST (y) with a commercial obtainable assay (x) gave the following result with 48 samples:

$y = 0.992x - 0.022$; $r = 0.995$

GOT AST

ASPARTATE AMINOTRANSFERASE

Quality control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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Text passages with grey background were changed in the latest revision of this package insert.

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* TRIS = Tris(hydroxymethyl)-aminomethan

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