Fluitest® GOT AST

ASPARTATE AMINOTRANSFERASE

Order information:

Catalog No.	Contents						
1176	R1	8 x	15 ml	R2	1 x	25 ml	
1177	R1	4 x	50 ml	R2	4 x	10 ml	
1178	R1	4 x	100 ml	R2	2 x	40 ml	
H7001	Hit I (ILab*)	R1	6 x	47 ml	R2	6 x	11 ml
H7003	Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	14 ml
AU7003	AU	R1	6 x	60 ml	R2	6 x	14 ml
1101	P-5-P	R0	6 x	3 ml			

(*) Kit contains only reagent barcode for Hitachi system

System information:

Without pyridoxal phosphate activation:

Hitachi 911: ACN 111
Hitachi 917: ACN 457, ACN 111 (STAT)

With pyridoxal phosphate activation:

Hitachi 911: ACN 143
Hitachi 917: ACN 143

Please note: For measurement with pyridoxal phosphate on Hitachi 911/912 analyzers special labels are required. These labels are available on request. For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

In vitro test for the quantitative determination of aspartate amino-transferase (AST) in human serum and plasma, with or without pyridoxal phosphate activation.

Summary:

Aspartate aminotransferase (glutamate-oxaloacetate-transaminase) belongs to the transaminases, which catalyze the interconversion of amino acids and α -ketoacids by transfer of amino groups. Aspartate aminotransferase is commonly found in human tissue. Although heart muscle is found to have the most activity of the enzyme, significant activity has also been seen in the brain, liver, gastric mucosa, adipose tissue, skeletal muscle, and kidneys. AST is present in both the cytoplasm and mitochondria of cells. In cases involving mild tissue injury, the predominant form of AST is that from the cytoplasm, with a smaller amount coming from the mitochondria. Severe tissue damage results in more of the mitochondrial enzyme being released. Elevated levels of the transaminases can signal myocardial infarction, hepatic disease, muscular dystrophy, and organ damage. In 1955, Karmen et al described the first kinetic determination of AST activity in serum. The International Federation of Clinical Chemistry (IFCC) recommended in 1977 and 1980 standardized procedures for AST determination, including optimization of substrate concentrations, employment of TRIS* buffers, preincubation of combined buffer and serum to allow side reactions with NADH to occur, substrate start, and optional pyridoxal phosphate activation. In 2002 the IFCC confirmed their recommendation and extended it to 37°C. This method is derived from the IFCC reference method.


*TRIS = Tris(hydroxymethyl)-aminomethane

Test principle:


Fluitest GOT AST without pyridoxal phosphate activation

UV test according to a standardized method:

- Sample and addition of R1 (buffer)
- Addition of R2 and start of reaction:

AST is the enzyme which catalyzes this equilibrium reaction. The oxaloacetate increase is measured in a subsequent indicator reaction which is catalyzed by malate dehydrogenase.



In the second reaction, NADH is oxidized to NAD. The rate of decrease in NADH (Measured photometrically) is directly proportional to the rate of formation of oxaloacetate, and thus the AST activity.


Fluitest GOT AST with pyridoxal phosphate activation

UV test according to a standardized method.

- Sample and addition of R1 (buffer/pyridoxal phosphate)
- Addition of R2 and start of reaction:

Pyridoxal phosphate is a coenzyme which transfers an amino group from one amino acid (aspartate) to another amino acid (glutamate), via the corresponding α -ketoacid. Addition of pyridoxal phosphate to the reaction mixture ensures maximum catalytic activity, while addition of α -ketoglutarate causes pyridoxal phosphate to be liberated, and thus reactivates the transaminase reaction. The increase in oxaloacetate is determined in an indicator reaction catalyzed by malate dehydrogenase.

NADH is oxidized to NAD⁺. The rate of the photometrically determined NADH decrease is directly proportional to the rate of formation of oxaloacetate and thus the AST activity.

Reagent Concentration:

R1:

Tris buffer pH 7.8 100 mmol/l
L-Aspartate 200 mmol/l
LDH 800 U/l
MDH 600 U/l

R2:

NADH₂ 0.18 mmol/l
2-Oxoglutarate 12 mmol/l

Preparation and stability:

Fluitest GOT AST without pyridoxal phosphate activation

Unopened kit components: Up to the expiration date at +2°C to +8°C

Serum start:

Mix 5 volumes of R1 with 1 volume of R2.
This solution is stable up to 10 days at +2°C to +8°C or
up to 1 day at +20°C to +25°C.

Substrate start:

R1: Ready for use.
R2: Ready for use.
Onboard stability: R1 28 days
R2 90 days

Fluitest GOT AST with pyridoxal phosphate activation

Unopened kit components: Up to the expiration date at +2°C to +8°C

Substrate start:

R1: Mix well 100 parts R1 and 1 part P-5-P-solution.
(e.g. 10ml R1 + 100µl P-5-P).
R2: Ready for use.

Onboard stability P-5-P activated working solution:

R1 6 days
R2 90 days

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma.

Stability: 24 hours at +20°C to +25°C
7 days at +2°C to +8°C

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at +2°C to +8°C.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Warning! R1 contains hazardous components.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure serum start:

Wavelength: Hg 334 nm, Hg 340 nm or Hg 365 nm
Temperature: +25 / +30 / +37°C
Cuvette: 1 cm light path
Zero adjustment: against air

Working solution 1000 µl
Serum/Plasma 100 µl

Mix, incubate for 1 min. at assay temperature and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes.

Calculation:

Hg 365 nm 3235 x ΔA /min
Hg 340 nm 1746 x ΔA /min
Hg 334 nm 1780 x ΔA /min

Fluitest® GOT AST

ASPARTATE AMINOTRANSFERASE

Manual procedure for substrate start:	
Wavelength:	Hg 334 nm, 340 nm or Hg 365 nm
Temperature:	+25 / +30 / +37°C
Cuvette:	1 cm light path
Zero adjustment:	against air
R1	1000 µl
Serum/Plasma	100 µl
Mix, incubate for 1 min. at assay temperature, then add:	
R2	200 µl
Mix, read initial absorbance and start stopwatch simultaneously. Repeat reading after exactly 1, 2 and 3 minutes.	
Calculation:	
Hg 365 nm	3823 x ΔA/min
Hg 340 nm	2063 x ΔA/min
Hg 334 nm	2103 x ΔA/min

Limitations - interference:

Criterion: Recovery within ± 10% of initial value.
 Hemolysis interferes due to AST activity from erythrocytes.
 Icterus: No significant interference up to an index I of 20 (approximate conjugated and unconjugated bilirubin: 20 mg/dl)
 Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).
 Lipemia (Intralipid): No significant interference up to an index L of 325 (approximate triglycerides concentration: 650 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
 Lipemia may cause absorbance flagging as a result of an absorbance increase.
 Blood samples should only be drawn prior to the administration of Sulfasalazine and Sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.
 The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Measuring range:

If the change of absorbance per minute is higher than ΔE/min 0.160 at 340 nm or ΔE/min 0.080 at 365 nm the sample has to be diluted. On instruments without rerun function, dilute the samples manually with 0.9% NaCl or distilled/deionized water (e.g. 1 + 3). Multiply the result by the appropriate dilution factor (e.g. factor 4).
 Dilution Limit Hitachi: 800 U/l or 13.34 µkat/l

Reference values:

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the AST results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Fluitest GOT AST without pyridoxal phosphate activation

According to the IFCC method.

	U/l	µkat/l
Men	<50	< 0.85
Women	< 35	< 0.60

Fluitest GOT AST with pyridoxal phosphate activation

According to DGKC:

+37°C	U/L	µkat/L
GOT AST/P5P men	10 – 50	0,17 – 0,83
GOT AST/P5P women	10 – 35	0,17 – 0,58

According to Tietz:

+37°C	U/L	µkat/L
GOT AST/P5P men	15 – 40	0,26 – 0,68
GOT AST/P5P women	13 – 35	0,22 – 0,60

Analytical sensitivity (lower detection limit):

Detection limit: 4 U/l or 0.07 µkat/l

The lower detection limit represents the lowest measurable AST concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined on Hitachi 911 using controls in an internal protocol between day (n = 20). The following results were obtained:

Between day			
Sample	Mean U/l	SD U/l	CV %
sample 1	43.5	2.05	4.71
sample 2	83.4	1.95	2.34
sample 3	127	2.92	2.30

Reproducibility was determined using controls within protocol between day (n = 20). The following results were obtained:

Within run			
Sample	Mean U/l	SD U/l	CV %
sample 1	39.8	1.13	2.83
sample 2	86.0	1.04	1.21
sample 3	135	1.51	1.12

Method comparison:

A comparison of the Analyticon Fluitest® AST (y) with a commercial obtainable assay (x) gave the following result with 58 samples:

$$y = 0.946x + 1.385; \quad r = 0.998$$

Quality Control:

Human Control Serum		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl		
S2: Bio Cal® E	10 x 3 ml	#1430

Calibration frequency:

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedures

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Text passages with grey background were changed in the latest revision of this package insert.

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