Fluitest® GPT ALT

ALANINE AMINOTRANSFERASE

BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B6501	200 / 600	R1	6 x	47 ml
		R2	6 x	11 ml
B6533	300 / 600*	R1	6 x	57 ml
		R2	6 x	14 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Fluitest GPT ALT without pyridoxal phosphate activation

Intended use:

In vitro test for the quantitative determination of alanine aminotransferase (ALT) in human serum and plasma.

Summary:

Alanine aminotransferase (glutamate-pyruvate-transaminase) belongs to the group of transaminases which catalyze the conversion of amino acids to the corresponding α -ketoacids via the transfer of amino groups; they also catalyze the reverse process. Although higher activities exist in the liver, minor activity can also be detected in the kidneys, heart, skeletal muscle, pancreas, spleen, and lungs. Elevated levels of transaminases are indicative of myocardial infarction, hepatopathies, muscular dystrophy, and damage to internal organs. Increased ALT activity in the serum, however, is a rather specific indicator of damage to the liver parenchyma, while AST is not necessarily a liver-specific parameter.

In 1956, Wroblewski and LaDue described the first kinetic determination of ALT in serum. In 1977 and 1980, the International Federation of Clinical Chemistry (IFCC) recommended standardized methods for the determination of ALT with optimized substrate concentrations, use of TRIS buffer*, simultaneous preincubation of serum with buffer (to avoid competing reactions with NADH), substrate start, and pyridoxal phosphate activation.

The method described here is derived from the IFCC reference method.

*TRIS = Tris(hydroxymethyl)-aminomethane


Test principle:

UV test according to the IFCC method

- Sample and addition of R1 (buffer)
- Addition of R2 and start of reaction:

ALT is the enzyme which catalyzes this equilibrium reaction. The pyruvate increase is measured in a subsequent indicator reaction which is catalyzed by lactate dehydrogenase.

In the second reaction, NADH is oxidized to NAD. The rate of decrease in NADH (Measured photometrically) is directly proportional to the rate of formation of pyruvate, and thus the ALT activity.

Reagent concentration:

R1:	
Tris buffer pH 7.8	100 mmol/l
L-Alanine	500 mmol/l
LDH	1200 U/l
R2:	
NADH ₂	0.18 mmol/l
2-Oxoglutarate	15 mmol/l

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability: R1: 60 days

R2: 60 days

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma.

Stability: 24 hours at +20°C to +25°C

3 days at +2°C to +8°C

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at +2°C to +8°C.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.

Hemolysis interferes due to ALT activity from erythrocytes.

Icterus: No significant interference up to an index I of 35 (approximate conjugated and unconjugated bilirubin: 35 mg/dl)

Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/dl).

Lipemia (Intralipid): No significant interference up to a triglyceride concentration of 950 mg/dl. There is poor correlation between turbidity and triglycerides concentration.

Blood samples should only be drawn prior to the administration of Sulfasalazine and Sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Measuring range:

3.0-300 U/l or 0.05-5 μ kat/l

Determine samples with higher activities via the rerun function.

Reference values:

Bases of the IFCC-method (37°C) and optimised standard-method of the DGKC (1972):

	U/l	μ kat/l
Men	<50	< 0.85
Women	< 35	< 0.60

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the ALT results should always be assessed in conjunction with the patient medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 3.69 U/l or 0.062 μ kat/l

The lower detection limit represents the lowest measurable ALT concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol within day (n = 20). The following results were obtained:

sample	Run to run		
	Mean (U/l)	SD (U/l)	CV %
sample 1	123.59	3.297	2.7
sample 2	20.80	0.852	4.1
sample 3	109.03	1.392	1.3

sample	Within run		
	Mean (U/l)	SD (U/l)	CV %
sample 1	121.06	0.982	0.8
sample 2	19.81	0.886	4.5
sample 3	105.88	0.646	0.6

Method comparison:

A comparison of the Analyticon Fluitest GPT ALT (y) with a commercial obtainable assay (x) gave with 166 samples the following result:

$$y = 0.845x + 1.798; \quad r = 0.995$$

Quality Control:

Human Control Serum		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits

Fluitest® GPT ALT

ALANINE AMINOTRANSFERASE

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency:

Two point calibration is recommended;

- after lot change
- as required following quality control procedures

Calibration verification: Not necessary.

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Text passages with grey background were changed in the latest revision of this package insert.

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* TRIS = Tris(hydroxymethyl)-aminomethan

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