GPT ALT

ALANINE AMINOTRANSFERASE

Order information:

Catalog No.	Contents		
11803	R1 1 x	65 ml	R2 20 x for 3 ml
118	R1 1 x	110 ml	R2 10 x for 10 ml
1184	R1 4 x	100 ml	R2 4 x for 100 ml

Intended use:

In vitro test for the quantitative determination of alanine amino-transferase (ALT) in human serum and plasma.

Summary:

Alanine aminotransferase (glutamate pyruvate transaminase) belongs to the group of transaminases which catalyze the conversion of amino acids to the corresponding α -keto acids via the transfer of amino groups; they also catalyze the reverse process. Although higher activity exists in the liver, minor activity can also be detected in the kidneys, heart, skeletal muscle, pancreas, spleen, and lungs. Elevated levels of transaminases are indicative of myocardial infarction, hepatopathies, muscular dystrophy, and damage to internal organs. Increased ALT activity in the serum, however, is a rather specific indicator of damage to the liver parenchyma, while AST is not necessarily a liver-specific parameter.


In 1956, Wroblewski and LaDue described the first kinetic determination of ALT in serum. In 1977 and 1980, the International Federation of Clinical Chemistry (IFCC) recommended standardized methods for the determination of ALT with optimized substrate concentrations, use of TRIS buffer*, simultaneous preincubation of serum with buffer (to avoid competing reactions with NADH), substrate start, and pyridoxal phosphate activation.

The method described here is derived from the IFCC Reference method.

*TRIS = Tris(hydroxymethyl)-aminomethane

Test principle:

Method recommended by the IFCC:

Reagent concentration:

R1:	
Tris buffer pH 7.8	100 mmol/l
L-Alanine	500 mmol/l
R2:	
NADH	0.18 mmol/l
LDH	1200 U/l
Oxoglutarate	15 mmol/l

Preparation and stability:

Dissolve enzyme reagent /R2 with the corresponding volume of buffer /R1. Before first use wait at least 15 min.

Reconstituted reagent:	2 weeks	at +2 to +8°C
Unopened vials:	up to the expiration date	at +2 to +8°C

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma.

Stability:	24 hours	at +20°C to +25°C
	3 days	at +2°C to +8°C

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at 2-8°C.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

 Warning! R2 contains hazardous material. The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within \pm 10% of initial value.

Icterus: No significant interference up to an index I of 37

(approximate conjugated and unconjugated bilirubin: 37 mg/dl)

Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 250 (approximate triglycerides concentration: 500 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

Lipemia may cause absorbance flagging as a result of an absorbance increase.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Manual procedure:

Wavelength:	Hg 334 nm, Hg 340 nm or Hg 365 nm
Temperature:	+25 / +30 / +37°C
Cuvette:	1 cm light path
Zero adjustment:	against air or distilled water

	Macro	Semi	Micro
Working reagent	2500 μ l	1000 μ l	500 μ l
Sample	250 μ l	100 μ l	50 μ l

Mix, incubate for 1 min. at assay temperature and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes and calculate the $\Delta A/\text{min}$.

If the $\Delta A/\text{min}$. is between 0.11 and 0.25 at Hg 334 nm/ 340 nm or 0.06 and 0.130 at Hg 365 nm use only the values for the first 2 minutes for the calculation.

Calculation:

	Macro	Semi	Micro
Hg 340 nm $\Delta A/\text{min} \times$	1750	1750	1750
Hg 334 nm $\Delta A/\text{min} \times$	1790	1790	1790
Hg 365 nm $\Delta A/\text{min} \times$	3240	3240	3240

Measuring /reportable range:

$\Delta E/\text{min}$ 0.250 at 340 nm or $\Delta E/\text{min}$ 0.130 at 365 nm

At higher activities dilute the sample with 0.9% NaCl solution (e.g. 1+9). Multiply the result by the appropriate dilution factor (e.g. 10)

Expected values:

Bases of the IFCC-method and optimised standard-method of the DGKC (1972)

		25°C	30°C	37°C
U/l	Men	up to 22	up to 29	up to 50
	Women	up to 17	up to 22	up to 35
$\mu\text{kat/l}$	Men	up to 0.37	up to 0.48	up to 0.85
	Women	up to 0.28	up to 0.37	up to 0.60

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the ALT results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 4 U/l or 0.07 $\mu\text{kat/l}$

The lower detection limit represents the lowest measurable ALT concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

within run			
Sample	Mean U/l	SD U/l	CV %
sample 1	30.6	0.59	1.92
sample 2	79.5	1.28	1.61
sample 3	119	1.08	0.91

Between day			
Sample	Mean U/l	SD U/l	CV %
sample 1	33.1	0.78	2.37
sample 2	93.3	1.25	1.34
sample 3	137	2.20	1.61

GPT ALT

ALANINE AMINOTRANSFERASE

Method comparison:

A comparison of the Analyticon ALT (y) with a commercial obtainable assay (x)

gave with 50 samples the following result:

$$y = 1.063 x + 0.263; \quad r = 0.996$$

Quality control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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Text passages with grey background were changed in the latest revision of this package insert.

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