H3-Check

Hemolyzer® 3 NG Control Set

Order information:

Catalog-No.	Article name	Content
HE3504	H3-Check Complete	6 x 2.5 mL (L, N, H)
HE3505	H3-Check Normal	6 x 2.5 mL (N)

Intended Use:

H3-Check Control Sets are designed for quality control measurements on the Analyticon hematology analyzer Hemolyzer® 3 NG (open and closed mode) (HE3100), and for validation of their measurement results.

Reagents:

The H3-Check Control Set contains human erythrocytes, simulated leukocytes, and mammalian platelets suspended in a plasma-like fluid with preservatives. For lot-specific target values, please refer to the table below.

This reagent contains non-hazardous concentrations of sodium azide (< 0.1 % (w/w)).

Warning and Precautions:

For in vitro diagnostic use only.

For trained laboratory professional only.

Sodium azide can form explosive compounds with heavy metals (e.g. copper and lead) which are commonly found in plumbing drains. Repeated contact of these reactants may result in the buildup of a shock sensitive layer.

All body fluid samples should be considered potentially infectious materials. The product was tested negative against anti-HIV 1 and 2, anti-HCV and HBsAg, but since no test method can rule out the danger of infection with absolute certainty, treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.

Exercise the normal precautions required for handling all laboratory reagents. Disposal of the product should be in accordance with local regulations.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Storage and Stability:

Unopened tube	Until the expiration date	2 – 8 °C
Opened tube	14 days	2 – 8 °C

Store the tubes in upright position. Do not freeze!

Indications of Deterioration:

After mixing, product should be similar in appearance to fresh whole blood. In unmixed vials, the supernatant may appear cloudy and reddish, which is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. Do not use the product if deterioration is suspected.

Procedure:

1. Allow the tubes to warm to room temperature (15 to 30 °C) for at least 15 minutes.
2. Invert the tubes for at least 20–30 seconds to mix the control material thoroughly. The red precipitate at the bottom of the tube has to be completely suspended and the red color of the fluid needs to be evenly distributed. Tubes stored for a long time might require longer mixing.
3. Gently invert the tube 8–10 times immediately before sampling and analyze the sample as instructed in the Quality Control section of the Operator's Manual for your instrument.
4. If the tube has been open for sampling, clean residual material from the cap and tube rim with a lint-free tissue after the sampling. Close the cap tightly.
5. Return the tubes to refrigerator (2–8 °C) within 30 minutes of use.

Results and Expected Values:

Verify that the lot number on the tube matches the lot number on the table of target values.

Target values are determined on well-maintained, properly calibrated instruments using the instrument manufacturer's recommended reagents. The "Expected Range" takes inherent imprecision of the method and expected biological variability of the control material into account.

Minor variations may derive from different reagent batches, maintenance status, operating technique, and calibration and might contribute to inter-laboratory variation.

For greater control sensitivity, it is recommended that each laboratory establishes its own target value and expected range, and periodically reevaluate the target value. The expected range of the laboratory may include values outside of the expected range given by the manufacturer.

Limitations of the Procedure:

The performance of this product is assured only if it is properly stored and handled as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

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