H3-Compact

Hemolyzer® 3 NG Reagent Pack

Order information:

Catalog-No.	Content	Volume	Total
HE3701	H3-DIL	1100 mL	100 Tests
	H3-LYSE	350 mL	
	H3/5-CLEAN	400 mL	
HE3705	H3-DIL	3500 mL	500 Tests
	H3-LYSE	750 mL	
	H3/5-CLEAN	1150 mL	

Intended use:

H3-Compact is a set of reagents intended to be used on the Analyticon hematology analyzer Hemolyzer® 3 NG (HE3100).

The reagents of the H3-Compact should only be replaced with Analyticon reagents. Erroneous results may be obtained if the reagent pack is applied with third party reagents. Please refer to the instrument's User Manual for further information.

Summary and Explanation:

HE3701: This reagent pack is designed to supply suitable amount of reagent to run 100 tests on the Analyticon hematology analyzer Hemolyzer® 3 NG under the condition that the user runs 3–5 measurements per day (open sampling mode). Further conditions that need to be fulfilled are listed under "Reagent Preparation and Replacement".

HE3705: This reagent pack is designed to supply suitable amount of reagent to run 500 tests on the Analyticon hematology analyzer Hemolyzer® 3 NG under the condition that the user runs 15 measurements per day (open sampling mode). Further conditions that need to be fulfilled are listed under "Reagent Preparation and Replacement".

Reagent Concentration:	H3-DIL	H3-LYSE	H3/5-CLEAN
Sodium chloride	< 1.5 %		< 1.5 %
Buffers	< 1.0 %	< 1.0 %	< 1.0 %
Stabilizers	< 0.5 %	< 0.5 %	< 0.5 %
Preservatives	< 0.5 %	< 0.5 %	< 0.5 %
Surfactants		< 3.5 %	

in ion-free water

Warning and Precautions:

For in vitro diagnostic use only.

For trained laboratory professional only.


The reagents are environmental friendly, azide-free and do not contain harmful ingredients. Exercise the normal precautions required for handling all laboratory reagents. Disposal of the product should be in accordance with local regulations.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Reagent Preparation and Replacement:

The reagent is ready for use.

1. Leave the reagent to reach room temperature before using it.
2. Open the caps of the bottles inside the reagent pack. Cut the sealing foil on the new reagents to allow the tubes go through. Immerse the caps with tubes into the reagents and tighten the caps. Make sure that the color on each tube, reagent container label and connector in the back of the instrument match (green=H3-DIL; yellow=H3-LYSE; white=H3/5-CLEAN).

Always keep the reagent in its original container and do not mix the remains of two reagent containers.

Place QR Code

3. Do not touch the tubes that go into the reagents. Avoid any dust or microbial contamination of the tubing and reagents. If you find the tubes dirty, you can clean them using a dampened lint-free wiping paper. Clean the tubes by swiping them top-down.
4. **REF HE3701:** Remove small tear flap and align the tubes so they are going through the upper edge of the pack. Close the pack lid but do not pinch the tubes.
REF HE3705: Do not close the pack lid to not pinch the tubes!
5. Prime thoroughly the new reagent and measure backgrounds according to the instrument's User Manual.
6. When installing a new lot of reagent, verify calibration by running control material. Recalibrate the instrument if necessary as specified in your User Manual.

7. Scan the QR code which can be found on the inside on the package insert of the pack (refer to the User Manual).

The reagent pack can be used for the designated number of tests, as long as the following conditions are fulfilled:

- The user measures more than the minimum number of tests per day (refer to "Summary and Explanation").
- The analyzer performs two Fast Blank measurements each day.
- The analyzer is powered off every day.
- The user performs necessary maintenance actions prompted by the analyzer's operating software.

The analyzer software will stop automatically after the specified number of tests were performed. Depending on usage, the amount of leftover reagent inside the package may vary.

Quality Control samples, Calibration samples, Open mode and Closed mode Samples will equally reduce the number of tests that can be performed with the package.

Running a lower number of tests per day than the conditions listed above will use relatively more reagents, and the pack will run out before the designated number of tests.

Fast Blank measurements consume less reagent than sample measurements, however, running Fast Blank measurements beyond the conditions listed above will also decrease the available number of tests.

In case the conditions described above are not fulfilled, there is no warranty that the designated number of tests that can be performed with the volume of the reagent pack. No liability is accepted for leftover reagent in the pack after all available tests have been performed.

Storage and Stability:

Unopened container	Up to the expiration date	2 – 35°C
Open container	6 months	2 – 35°C

Keep the reagent containers closed when not in use.

Do not freeze!

Specimen Collection and Preparation:

H3-Compact is intended for use with blood specimens collected by venipuncture in EDTA anticoagulant.

Specimens for hematological analysis may be stored for up to 6 hours at 15–30°C or up to 10 hours after collection when refrigerated at 2–8°C.

If the specimen is refrigerated at 2–8°C, please wait 30 minutes for the specimens to warm up before proceeding to the measurement.

All body fluid samples should be considered potentially infectious materials. Treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.

Text passages with grey background were changed in the latest revision of this package insert.

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