# H5-LYSE

Hemolyzer® 5 NG Lyse

## Order information:

Catalog-No.	Content
HE5111	1 L

## Intended use:

H5-LYSE is a stabilized and micro-filtered lysing agent for stromatolysis of erythrocytes (RBC), for quantitative determination of leukocytes (WBC), leukocyte five-part differentiation (LYM, MID, NEU, EOS, BAS) and hemoglobin (HGB) concentration measurement of human blood on Analyticon hematology analyzers Hemolyzer® 5 NG (closed mode and open mode) (HE5100 and HE5200).

H5-LYSE should only be used with Analyticon reagents. Erroneous results may be obtained if the reagent is applied with third party reagents. Please refer to the instrument's User Manual for further information.

## Reagent Concentration:

Surfactants	< 0.5 %
Buffers	< 1.8 %
Preservatives	< 1.0 %
Stabilizers	< 1.1 %

in ion-free water

## Warning and Precautions:

For in vitro diagnostic use only.

For trained laboratory professional only.

The reagents are environmental friendly, azide-free and do not contain harmful ingredients. Exercise the normal precautions required for handling all laboratory reagents. Disposal of the product should be in accordance with local regulations.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

## Reagent Preparation and Replacement:

The reagent is ready for use.

1. Leave the reagent to reach room temperature before using it.
2. Remove the inlet cap from the reagent container to be replaced. Cut the sealing foil on the new reagent to allow the tubes go through. Connect the reagent inlet to the new reagent container. Make sure that the color on each tube, reagent container label and connector in the back of the instrument match.
3. Avoid any dust or microbial contamination of the tubing and reagents. Always keep the reagent in its original container and do not mix the remains of two reagent containers.
4. Prime thoroughly the new reagent and measure backgrounds according to the instrument's User Manual.
5. When installing a new lot of reagent, verify calibration by running control material. Recalibrate the instrument if necessary as specified in your User Manual.

## Storage and Stability:

Unopened container	Up to the expiration date	2 – 35°C
Open container	6 months	2 – 35°C

Keep the reagent container closed when not in use.  
Do not freeze!

## Specimen Collection and Preparation:

H5-LYSE is intended for use with blood specimens collected by venipuncture in EDTA anticoagulant.

Specimens for hematological analysis may be stored for up to 6 hours at 15–30°C or up to 10 hours after collection when refrigerated at 2–8°C.

If the specimen is refrigerated at 2–8°C, please wait 30 minutes for the specimens to warm up before proceeding to the measurement.

All body fluid samples should be considered potentially infectious materials. Treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.

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