Fluitest® HDL-CHOL

HDL-CHOLESTOL PRECIPITATION REAGENT

Order information:

Catalog-No.	Contents
410	R1 6 x 40 ml

Intended use:

Precipitation reagent for sample pre-treatment for the in vitro quantitative determination of HDL-cholesterol in human serum and plasma.

Summary:

HDL (High Density Lipoproteins) are responsible for the reverse transport of cholesterol from the peripheral cells to the liver. Here, cholesterol is transformed to bile acids which are excreted into the intestine via the biliary tract. Monitoring of HDL-cholesterol in serum is of clinical importance since an inverse correlation exists between serum HDL-cholesterol concentrations and the risk of arteriosclerotic disease. Elevated HDL-cholesterol concentrations are protective against coronary heart disease, while reduced HDL-cholesterol concentrations, particularly in conjunction with elevated triglycerides, increase the cardiovascular risk.

A variety of methods are available to determine HDL-cholesterol, including ultracentrifugation, electrophoresis, HPLC, and precipitation-based methods. Of these precipitation-based methods are used routinely. HDL-cholesterol is first separated by precipitating apoprotein B-containing lipoproteins from serum by using a combination of a polyanion and a divalent cation, such as dextran sulfate/magnesium chloride or phosphotungstate/magnesium chloride. Such precipitation-based methods are, however, time-consuming and not amenable to automated analysis. Thus, there is a great clinical need for a convenient and reliable method for measuring HDL-cholesterol in serum without any pre-treatment. Several approaches for direct measurement of HDL-cholesterol in serum have been proposed, including the use of magnetically responsive particles as polyanion-metal combinations and the use of polyethylene glycol (PEG) with anti-apoprotein B and anti-apoprotein CIII antibodies.

Test principle:

The chylomicrons, VLDL (very low density lipoproteins) and LDL (low density lipoproteins) are precipitated by addition of phosphotungstic acid and magnesium chloride. After centrifugation the supernatant fluid contains the HDL (high density lipoproteins)-fraction, their cholesterol content is determined enzymatically by the CHOD-PAP test.

Reagent Concentration:

R1	
Phosphotungstic acid	0.55 mmol/l
Magnesium chloride	25 mmol/l

Preparation and Stability:

- Precipitant for macro assays. Use contents undiluted.
- Precipitant for semi-micro assays:
Dilute 4 parts of precipitating reagent with 1 part of redistilled water (e.g. 80 ml + 20 ml)

The HDL reagent is stable up to the expiry date when stored at +15 to +25°C.

Specimen:

Collect serum using standard sampling tubes

Serum, heparin- or EDTA-plasma

Do not use citrate-, oxalate- or fluoride-plasma.

Stability: 7 days at +2°C - +8°C

3 months at -20°C

Fasting and non-fasting samples can be used. EDTA plasma causes decreased results.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Interferences with the Chol measurement.

See manual instruction of Analyticon Chol.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Analytical sensitivity (lower detection limit)

Detection limit: 3 mg/dl (0.08 mmol/l)

The lower detection limit represents the lowest measurable HDL-cholesterol concentration that can be distinguished from zero.

Testing procedure:

Materials provided

- Working solutions as described above

Additional materials required

- Analyticon CHOL Test-Kit #4046 or #4010
- Calibrators and controls as indicated below
- centrifuge
- Appropriate centrifuge tubes

Manual testing procedure for PRECIPITATION:

Pipette into centrifuge tubes.

	Macro	Semi Micro
Sample	500 µl	200 µl
HDL reagent undiluted	1000 µl	---
HDL reagent diluted	---	500 µl

Mix well, allow to stand for 10 min at +15°C to +25°C and centrifuge for 2 min at 1000g or 10 min at 4000g. After centrifugation separate the clear supernatant from the precipitate within 1 hour and determine the cholesterol concentration.

Manual testing procedure for CHOLESTEROL DETERMINATION:

Wavelength:	Hg 546nm, 500nm
Temperature:	+25 / +30 / +37°C
Cuvette:	1cm light path
Zero adjustment:	reagent blank

	Blank	Sample
Pipette into cuvette		
Bidest. Water	100 µl	---
HDL supernatant	---	100 µl
CHOL reagent	1000 µl	1000 µl

Mix, measure after incubating at +37°C for 5 min or 10 min at +20 to +25°C. Read absorbance against reagent blank within 60min.

Calculation:

Wavelength	Macro	Semi Micro
Hg 546 nm	280 x ΔA	327 x ΔA
500 nm	188 x ΔA	220 x ΔA

Clinical interpretation:

		Lipid disorder
Cholesterol	< 200 mg/dl (5.1 mmol/l)	No
Triglycerides	< 200 mg/dl (2.3 mmol/l)	
Cholesterol	200 – 300 mg/dl 5.1 – 8.6 mmol/l	Yes, if HDL chol. <35mg/dl
Cholesterol	> 200 mg/dl (5.1 mmol/l)	Yes
Triglycerides	> 200 mg/dl (5.1 mmol/l)	

Expected values:

	No risk	Moderate risk	High risk
Men			
mg/dl	> 55	35 - 55	< 35
mmol/l	> 1.45	0.90 – 1.45	< 0.90
Women			
mg/dl	> 65	45 – 65	< 45
mmol/l	> 1.68	1.15 – 1.68	< 1.15

National Cholesterol Education Program (NCEP) guidelines:

< 35 mg/dl Low HDL-Cholesterol (major risk factor for CHD)

>60 mg/dl High HDL-Cholesterol („negative“ risk factor for CHD) HDL-Cholesterol is affected by a number of factors, e.g. smoking, exercise, hormones, sex and age.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes the HDL-cholesterol results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Notes:

The values obtained are reliable, provided that

- no chylomicrons are present in the sample
- the triglyceride concentration does not exceed 400 mg/dl
- the sample does not show signs of type III hyperlipoproteinemia

In measurement at Hg 546 nm, the spectral properties of hemoglobin simulate elevated HDL-cholesterol values which can be ignored up to 200 mg Hb/100 ml.

The supernatant obtained on centrifugation must be clear. If the sample has a high triglyceride content (above 1000mg/dl), lipoprotein precipitation may be incomplete (cloudy supernatant), or part of the precipitate may float on the surface. In these cases, dilute the specimen 1:1 with 0.9 % NaCl solution and repeat the precipitation step.

The result of the cholesterol assay must then be multiplied by 2. High concentrations of ascorbic acid may result in artificially low values.

Imprecision:

Reproducibility was determined using controls in an internal protocol.
The following results were obtained:

Between day			
Sample	Mean mg/dl	SD mg/dl	% CV
Control serum 1	-	-	-
Control serum 2	43.6	0.84	1.92
Control serum 3	69.8	1.74	2.50

Within run			
Sample	Mean mg/dl	SD mg/dl	% CV
Control serum 1	46.57	0.53	1.15
Control serum 2	78.49	1.91	2.44
Control serum 3	-	-	-

Quality control:

Human Control Serum:

Contronorm® L 5 x 2 ml #1302
Controptath® L 5 x 2 ml #1303

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Disposal:

Please note the legal regulations.

Literature:

1. Assmann G. At. what levels of total low-or high-density lipoprotein cholesterol should diet/drug therapy be initiated European guidelines. Amer J Cardiol 1990;65:11F
2. Assmann G Schriewer H Schmitz G et al Qualification of high density lipoprotein cholesterol by precipitation with phosphotungstic acid/Mg/Cl₂ Clin Chem 1983;29:2026-2030
3. AVP Fettstoffwechselstörungen, Therapieempfehlungen 1, 1st 1996:2-16
4. Burstein M. Scholnick HR. Morfix R. Rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. J Lipid Res 1970;11:583-595
5. Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459
6. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
7. Harris N. Galpchian V, Rifai N Three routine methods for measuring high density lipoprotein cholesterol compared with the Reference method. Clin Chem 1996;42:738-743
8. Hatch FT, Lees RS. Practical methods for plasma lipoprotein analysis. Adv Lipid Res 1968;6:1-68
9. Kakuyama T, Kimura S Hashiguchi Y, Fully automated determination of HDL-cholesterol from human serum with Hitachi 911 (Abstract) Clin Chem 1994;40:1104
10. Mustro J. Lawlor JF. HDL-Cholesterol online separation and analysis utilizing an automated chemistry analyser (Abstract) Clin Chem 1993;39:1125
11. Narayan KA, Kummerow FA- Disk electrophoresis of human serum lipoprotein. Nature 1965;205:246-248
12. Nauck M, März W. Jarausch J et al Multicenter evaluation of a homogeneous assay for HDL- cholesterol without sample pretreatment Clin Chem 1997;43:1622-1629
13. Okazaki M, Shiraishi K, Ohno Y et al Heterogeneity of human high density lipoprotein. on high performance liquid chromatography J. Biochem 1982;92:517-524
14. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two different Analytical Methods. J Clin Chem Clin Biochem 1983;21: 709-720
15. Pisani T. Gebbski CP. Leary ET, et al Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127
16. Sugiuchi H. Uji Y, Okabe H. Irie T et al Direct Measurement of High-Density Lipoprotein Cholesterol in Serum with polyethylene Glycol-Modified enzymes and Sulfated α -Cyclo-dextrin. clin chem. 1995;41:717 – 72
17. Matsuzaki Y, Kawaguchi E. Norita Y et al Evaluation of Two Kinds of Reagents for Direct Determination of HDL-Cholesterol. J. Anal Bio Sc 1996;19:419-427
18. Thomas L (ed.). Labor und Diagnose, 4th edition. Marburg: Die Medizinische Verlagsgesellschaft, 1992:208
19. Second report of the Expert Panel on Detection Evaluation and Treatment of High Blood Cholesterol in Adults NIH Publication No 93-3096 September 1993

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Literatur:

1. Assmann G. At. what levels of total low-or high-density lipoprotein cholesterol should diet/drug therapy be initiated? European guidelines. Amer J Cardiol 1990;65:11F
2. Assmann G Schriewer H Schmitz G et al Qualification of high density lipoprotein cholesterol by precipitation with phosphotungstic acid/Mg/Cl₂ Clin Chem 1983;29:2026-2030
3. AVP Fettstoffwechselstörungen, Therapieempfehlungen 1, 1. Auflage 1996:2-16
4. Burstein M. Scholnick HR. Morfix R. Rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. J Lipid Res 1970;11:583-595
5. Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459
6. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
7. Harris N. Galpichian V, Rifai N Three routine methods for measuring high density lipoprotein cholesterol compared with the Reference method. Clin Chem 1996;42:738-743
8. Hatch FT, Lees RS. Practical methods for plasma lipoprotein analysis. Adv Lipid Res1968;6:1-68
9. Kakuyama T, Kimura S Hashiguchi Y, Fully automated determination of HDL-cholesterol from human serum with Hitachi 911 (Abstract) Clin Chem 1994;40:1104
10. Mastro J. Lawlor JF. HDL-Cholesterol online separation and analysis utilizing an automated chemistry analyser (Abstract) Clin Chem 1993;39:1125
11. Narayan KA, Kummerow FA- Disk electrophoresis of human serum lipoprotein. Nature 1965;205:246-248
12. Nauck M, März W. Jarausch J et al Multicenter evaluation of a homogeneous assay for HDL- cholesterol without sample pretreatment Clin Chem 1997;43:1622-1629
13. Okazaki M, Shiraishi K, Ohno Y et al Heterogeneity of human high density lipoprotein. on high performance liquid chromatography J. Biochem 1982;92:517-524
14. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two different Analytical Methods. J Clin Chem Clin Biochem 1983;21: 709-720
15. Pisani T. Gebski CP. Leary ET, et al Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127
16. Suguchi H. Uji Y, Okabe H. Irie T et al Direct Measurement of High-Density Lipoprotein Cholesterol in Serum with polyethylene Glycol-Modified enzymes and Sulfated α -Cyclo-dextrin. clin chem. 1995;41:717 – 72
17. Matsuzaki Y, Kawaguchi E. Norita Y et al Evaluation of Two Kinds of Reagents for Direct Determination of HDL-Cholesterol. J. Anal Bio Sc 1996;19:419-427
18. Thomas L (Hrsg.). Labor und Diagnose, 4. Auflage. Marburg: Die Medizinische Verlagsgesellschaft, 1992:208
19. Second report of the Expert Panel on Detection Evaluation and Treatment of High Blood Cholesterol in Adults NIH Publication No 93-3096 September 1993