Order information:

Catalog No.		Contents					
H11081	Hit I (ILab*)	R0	1 x	60 ml	R1	1 x	30 ml
		R2	1 x	15 ml	R3	1 x	15 ml
H11083	Hit 917 (AU*)	R0	1 x	60 ml	R1	1 x	30 ml
		R2	1 x	15 ml	R3	1 x	15 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 361-400 (user defined method)
 Hitachi 917: ACN 901-905 (user defined method)
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunological assay for the in vitro quantitative determination of hemoglobin A1c in whole blood on automated clinical chemistry analyzers. HbA1c determinations are utilized in the long-term monitoring of glycaemia.

Summary:

Hemoglobin (Hb) consists of four protein chains with four heme portions, and is the red-pigmented protein located in the erythrocytes. Its main function is the transport of oxygen and carbon dioxide in blood. Each Hb molecule is able to bind four oxygen molecules. Hb consists of a variety of subfractions and derivatives. Among this heterogeneous group of hemoglobins HbA1c is one of the glycosylated hemoglobins, a subfraction formed by the attachment of various sugars to the Hb molecule. HbA1c is formed in two steps by the nonenzymatic reaction of glucose with the N-terminal amino group of the β-chain of normal adult Hb (HbA). The first step is reversible and yields labile HbA1c. This slowly rearranges in the second reaction step to yield stable HbA1c.

In the erythrocytes, the relative amount of HbA converted to stable HbA1c increases with the average concentration of glucose in the blood. Due to the erythrocyte's life span of approximately 100 to 120 days, HbA1c reflects the average blood glucose level during the preceding 2 to 3 months. HbA1c is thus suitable to monitor long-term blood glucose control in individuals with diabetes mellitus. More recent glucose levels have a greater influence on the HbA1c level.¹

The approximate relationship between HbA1c and mean blood glucose value during the preceding 2 to 3 months has been analyzed by several studies. A recent study obtained the following correlation:

DCCT / NGSP Standardization

Mean plasma glucose [mmol/l] = 1.98*HbA1c (%) - 4.29 or

Mean plasma glucose [mg/dl] = 35.6*HbA1c (%) - 77.3

IFCC Standardization (not applicable for US customers, recalculated acc. to Ref. 8)

Mean plasma glucose [mmol/l] = 1.73*HbA1c (%) + 0.20 or

Mean plasma glucose [mg/dl] = 31.2*HbA1c (%) + 3.51

The risk of diabetic complications, such as diabetic nephropathy and retinopathy increases with poor metabolic control. In accordance with its function as an indicator for the mean blood glucose level, HbA1c predicts the risk for the development of diabetic complications in diabetes patients.

For routine clinical use, testing every 3 to 4 months is generally sufficient. In certain clinical situations, such as pregnancy diabetes, or after a major change in therapy, it may be useful to measure HbA1c in 2 to 4 week intervals.

Test principle:

Immunoturbidimetric test according DCCT/NGSP protocol.

Both the concentration of HbA1c and the concentration of total hemoglobin are measured. The reported HbA1c result is calculated as a % of the total hemoglobin concentration.

The HbA1c and total hemoglobin values generated in this assay are intended for use in the calculation of the HbA1c/total hemoglobin ratio (%HbA1c) and must not be used individually for diagnostic purposes.

(a) Sample Pre-treatment

The first step of the procedure involves the pre-treatment of the whole blood sample. This lyses red blood cells and causes hydrolysis of the hemoglobin by the action of a protease enzyme in the Hemoglobin Denaturant Reagent.

(b) Determination of Total Hemoglobin

The Total Hemoglobin reagent is used to determine the concentration of total hemoglobin. The method involves the conversion of all the hemoglobin derivatives into haematin in an alkaline solution of a non-ionic detergent as described by Wolf *et al* (1984).

The reaction is initiated by the addition of the pre-treated sample to the total hemoglobin reagent, resulting in a green solution. The conversion of different hemoglobin species into alkaline haematin with one defined absorption spectrum allows the endpoint measurement of total hemoglobin at 600 nm.

(c) Determination of HbA1c


The determination of HbA1c is based on a latex agglutination inhibition assay. The agglutinator, which consists of a synthetic polymer containing multiple copies of the immunoreactive portion of HbA1c, causes agglutination of latex coated with HbA1c specific mouse monoclonal antibodies.

In the absence of HbA1c in the sample, the agglutinator in the HbA1c R3 Reagent and the antibody-coated micro particles in the HbA1c R2 Reagent will agglutinate, resulting in an increase in absorbance.

The presence of HbA1c in the sample will slow the rate of agglutination as it competes with the HbA1c agglutinator for antibody binding sites on the latex.

Hence, the increase in absorbance is inversely proportional to the concentration of HbA1c in the sample.

An increase in absorbance due to agglutination is measured at 700 nm and the extent of agglutination is used to calculate the concentration of HbA1c from a Calibration Curve. The percentage HbA1c is then calculated using the g/dl HbA1c and Total Hemoglobin values.

Working solution concentration

R0: Hemoglobin Denaturant

Porcine Pepsin
 Buffer pH 2.4
 Preservative

R1: Total Hemoglobin Reagent

Sodium Hydroxide, pH 13 0.4% w/v
 Triton 2.5% w/v
 Octylphenoxypolyethoxyethanol 2.5% w/v

R2: HbA1c Antibody Reagent

HbA1c Antibody (mouse) coupled particles <0.1% w/v
 Bovine Serum Albumin
 Buffer pH 8.1
 Non ionic Surfactant 0.6% w/v
 Proclin 150 0.1% w/v

R3: HbA1c Agglutinator Reagent

HbA1c hapten covalently attached to the polymer
 Bovine Serum Albumin
 Buffer pH 2.0
 Proclin 150 0.1% w/v
 Non ionic Surfactant 0.2% w/v

Precautions and warnings:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Do not pipette by mouth.

Total Hb Reagent and Calibrator set Level 1 contain sodium hydroxide which is caustic. Hemoglobin Denaturant contains acidic buffer which is caustic too. In the event of accidental contact, flush affected area with large quantities of water and seek immediate medical attention.

The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.

Storage and stability:

All reagents are ready for use.

All reagents are stable as supplied to expiry date when stored at +2° to +8°C, stored protected from extreme heat, light, or freezing.

Reagents should be mixed thoroughly and equilibrated to system temperature for approximately ½ hour prior to use on the system. The reagents are stable for 30 days on board the analyser at approximately 10°C.

Specimen collection and preparation:

Capillary blood, Potassium-EDTA and Ammonium Heparin whole blood

Stability: at +15 to +25°C 3 days
 at + 2 to + 8°C 5 days
 at - 20°C 6 months (do not refreeze)

Sample preparation:

Capillary blood

Puncture capillary blood and put filled capillary tube into a reagent tube with 400 µl of hemoglobin denaturant reagent R0 (1:41 dilution; 10 µl whole blood sample with 400 µl of hemoglobin denaturant reagent). Close reagent tube and mix gently and thoroughly with vibration swirl

Incubate for a minimum of 5 minutes at room temperature prior to testing.

Stability of sample: at + 2 to + 8°C 48 hours
at +15 to +25°C 8 hours

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Hb Denaturant Reagent R0
- Total Hb Reagent R1
- HbA1c R2 : Antibody Reagent
- HbA1c R3 : Agglutinator Reagent

Additional materials required

- HbA1c/Hb Calibration
- HbA1c Quality Control
- General laboratory equipment

Procedure:

Refer to Hitachi operator's manual for user defined applications.

Hitachi 917: Prior to testing, put the barcodes, according to your chosen application code, on the bottles.

e.g. ACN 901 for Hb testing → barcode 901 for R1 Hb on bottle R1
ACN 902 for A1c testing → barcode 902 for R1 A1c on bottle R2
and barcode 902 for R3 A1c on bottle R3

Place R1 (Hb reagent) on reagent plate 1, R2 (A1c antibody reagent) on reagent plate 1 and R3 (A1c agglutination reagent) on reagent plate 2.

Calculation:

The calculation of the hemoglobin A1c concentration is generated using the following equation:

$$\% \text{HbA1c} = \frac{\text{HbA1c [g/dl]}}{\text{Hb [g/dl]}} \times 100$$

Note: This test is designed for DCCT/NGSP standardisation. In order to convert the results into IFCC units according to the directive of German Medical Association 2008 use the following equation:

$$\text{HbA1c IFCC [mmol/mol]} = 10.93 \times (\% \text{HbA1c NGSP} - 2.15)$$

Limitations – interference:

The assay gives accurate and precise results for a range of total hemoglobin varying between 7g/dl and 23g/dl. Patients with severe anaemias (total hemoglobin < 7g/dl) and those with polythemia (total hemoglobin > 23g/dl) should not be assayed by this method.

Any case of shortened red cell survival such as haemolytic anaemia or other haemolytic disease, pregnancy, recent significant blood loss etc. will result in a decrease in % glycated hemoglobin.

Samples containing hemoglobin variants S and C may result in up to a 40% increase of the expected HbA1c value. Also, those containing variant F (>10%) may yield a lower value than expected. Hence, samples containing variants S,C and F (>10%) should not be compared to published normal or abnormal values.

A sample containing variant E was shown not to interfere.

Criterion: Recovery within ± 10% of initial value.

Icterus: No significant interference up to 30 mg/dl or 513 µmol/l (conjugated and unconjugated bilirubin).

Lipemia (Intralipid): No significant interference up to a triglycerides concentration of 1600 mg/dl or 18 mmol/l. There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factor 2000 IU/ml, Acetylsalicylic acid < 60 mg/dl, Sodium cyanate 50 mg/dl and Urea 500 mg/dl (83 mmol/l) do not interfere.

The labile fraction of glycated hemoglobin (Schiff base) does not interfere as the antibody is specific for the stable ketamine.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Analytical range:

The analytical range for total hemoglobin is 7 g/dl to 23 g/dl.

The analytical range for % HbA1c is the concentration that corresponds to the level 6 HbA1c calibrator (2.06 g/dl HbA1c, 14.7% HbA1c at a total Hemoglobin of 14 g/dl).

Samples with values above 14.7% HbA1c should not be diluted and results should be reported as >14.7% HbA1c.

Sensitivity

The assay is sensitive to a level of 0.3 g/dl HbA1c, 1.38 g/dl total hemoglobin. The % HbA1c sensitivity was found to be 2.571%. This corresponds to 0.3 g/dl HbA1c and 11.667 g/dl total hemoglobin.

Expected values:

Depending upon the assay used, HbA1c is approximately

4 – 6% (DCCT/NGSP) or 20 – 42 mmol/mol (IFCC) in nondiabetics

6 – 8% (DCCT/NGSP) or 42 – 64 mmol/mol (IFCC) in controlled diabetics

As much as 20% (DCCT/NGSP) or 195 mmol/mol (IFCC) in uncontrolled diabetics.

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, hemoglobin A1c results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Imprecision:

Reproducibility was determined using controls and pooled whole blood in an internal protocol. The following results were obtained.

Hb	Within run			Between day		
	Mean g/dl	SD g/dl	CV %	Mean g/dl	SD g/dl	CV %
Level 1	14.440	0.179	1.24	15.445	0.236	1.53
Level 2	15.046	0.185	1.23	14.873	0.238	1.60
Level 3	15.167	0.208	1.37	15.807	0.240	1.52

HbA1c	Within run			Between day		
	Mean g/dl	SD g/dl	CV %	Mean g/dl	SD g/dl	CV %
Level 1	0.818	0.021	2.63	0.824	0.018	2.13
Level 2	1.251	0.021	1.68	1.179	0.021	1.77
Level 3	1.645	0.027	1.64	1.662	0.024	1.43

%HbA1c	Within run			Between day		
	Mean %	SD %	CV %	Mean %	SD %	CV %
Level 1	5.663	0.194	3.43	5.337	0.163	3.06
Level 2	8.316	0.177	2.13	7.927	0.238	3.00
Level 3	10.845	0.249	2.30	10.518	0.251	2.38

Method comparison:

A comparison of the Analyticon HbA1c (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.96x + 0.4532;$$

and a correlation coefficient of 0.9747. 40 patient samples were analysed spanning the range 5.36 % to 12.88 %.

This included 7 normal and 33 abnormal patient samples.

Calibration

Standardization:

HbA1c

S1: 0.9% NaCl

S2-6: Bio Cal[®] HbA1c Calibration set – Level 2 – 6 #14081

We recommend that this assay should be calibrated using Analyticon HbA1c Calibrator Series, levels 2 – 6.

Total Hemoglobin

S1: 0.9% NaCl

S2: Bio Cal[®] HbA1c Calibration set – Level 1 only #14081

We recommend that this assay should be calibrated using 0.9% NaCl solution and Analyticon HbA1c Calibrator Series, level 1 only.

A multipoint calibration is recommended for HbA1c every 30 days, with change of reagent lot/bottle or as indicated by quality control procedures.

Note. Calibrators DO NOT REQUIRE PRETREATMENT.

The calibrators are referenced to an HPLC method for HbA1c and the Drabkin's method for total Hb.

Quality control:

HbA1c Control Set 2 x 0.25 ml #1508

Analyticon HbA1c Controls are recommended for daily quality control.

Two levels of controls should be assayed at least once a day. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error the following steps should be taken:

1. Check instrument settings and light source.
2. Check cleanliness of all equipment in use.
3. Check water, contaminants i.e. bacterial growth may contribute to inaccurate results.
4. Check reaction temperature.

Disposal:

Please note the legal regulations.

References:

1. Cohen P.M. Perspective: measurement of Circulating Glycated Protein to Monitor Intermediate – Term Changes in Glycaemic Control Eur J Clin Chem. Clin. Biochem. 1992;30 (12): 851 – 859
2. The Diabetes Control of complications Trial Research on the Development and Progression of Long – Term Complicity of insulin – dependant Diabetes Mellitus. The New England Journal of Medicine 1993;329 (14):977 – 986.
3. Mayer T.K. and Freedman Z.R.: Protein glycosylation in diabetes mellitus: A review of laboratory measurements and of their clinical utility. Clin. Chem. Acta 127: 147 – 184 (1983)
4. Baynes J.W., Bunn H.F., Goldstein D.E. *et al*: National Diabetes Group: Report of the expert committee on glucosylated hemoglobin. Diabetes Care 7: 602 – 606 (1984).
5. Koenig R.J., Petersen C.M., Kilo C *et al*: Hemoglobin A1c as an indicator of the degree of glucose intolerance in diabetes. Diabetes 25: 230 – 232 (1976).
6. Nathan D.M., Singer D.E., Hurxthal K. and Goodson J.D.: The clinical information value of the glycosylated hemoglobin assay. NE J Med 310: 341 – 346 (1984).
7. McCarren M: DCCT and its works: Intensive therapy reduces the risk of diabetic eye, kidney, and nerve disease. Diabetes Forecast 49 – 51 (September 1993).
8. Larsen M.L., Horder M, and Mogensen E.F.: Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. NE J Med. 323: 1021 – 1025 (1990).
9. Nathan D.M.: Hemoglobin A1c – Infatuation or the real thing? NE J Med. 323: 1062 – 1063 (1990).
10. Wolf HU, Lang W., and Zander R/: Alkaline haematin D – 575, a new tool for the determination of hemoglobin as an alternative to the Cyanhemoglobin method. Clin. Chem. Acta 136: 83 – 104 (1984).
11. National Committee for Clinical Laboratory Standards: Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture – third edition: approved standard. NCCLS Publication H4 – A3 Villanova, PA: NCCLS 1991.
12. Goldstein D.E., Little R.R, Wiedmeyer H.M., *et al*. Glycated Hemoglobin: Methodologies and clinical applications. Clin. Chem. 32: B64 – B70 (1986).
13. Ellis G., Diamonds E.P., Giesbrecht E.E. , *et al*: An automated “high-pressure” liquid chromatographic assay for hemoglobin A1c. Clin. Chem 30: 1746 – 1752 (1984).
14. Burtis C.A. and Ashwood E.R. (eds.) : Tietz textbook of Clinical Chemistry, 2nd edition. W.B. Sandrs Company, Philadelphia, PA, p.2021 (1994)
15. Knowles W.J., Haigh W.B., and Michaud G.C. : A monoclonal antibody-based immunoassay for hemoglobin A1c Diabetes 35 Supplement 94A (1986).

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HbA1c

S1: 0,9% NaCl

S2-6: Bio Cal® HbA1c Calibration set – Level 2 – 6 #14081

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Literatur:

1. Cohen P.M. Perspective: measurement of Circulating Glycated Protein to Monitor Intermediate – Term Changes in Glycaemic Control Eur J Clin Chem. Clin. Biochem. 1992;30 (12): 851 – 859
2. The Diabetes Control of complications Trial Research on the Development and Progression of Long – Term Complicity of insulin – dependant Diabetes Mellitus. The New England Journal of Medicine 1993;329 (14):977 – 986.
3. Mayer T.K. and Freedman Z.R.: Protein glycosylation in diabetes mellitus: A review of laboratory measurements and of their clinical utility. Clin. Chem. Acta 127: 147 – 184 (1983)
4. Baynes J.W., Bunn H.F., Goldstein D.E. *et al*: National Diabetes Group: Report of the expert committee on glycosylated hemoglobin. Diabetes Care 7: 602 – 606 (1984).
5. Koenig R.J., Petersen C.M., Kilo C. *et al*: Hemoglobin A_{1c} as an indicator of the degree of glucose intolerance in diabetes. Diabetes 25: 230 – 232 (1976).
6. Nathan D.M., Singer D.E., Hurxthal K. and Goodson J.D.: The clinical information value of the glycosylated hemoglobin assay. NE J Med 310: 341 – 346 (1984).
7. M^cCarren M: DCCT and its works: Intensive therapy reduces the risk of diabetic eye, kidney, and nerve disease. Diabetes Forecast 49 – 51 (September 1993).
8. Larsen M.L., Horder M, and Mogensen E.F.: Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. NE J Med. 323: 1021 – 1025 (1990).
9. Nathan D.M.: Hemoglobin A_{1c} – Infatuation or the real thing? NE J Med. 323: 1062 – 1063 (1990).
10. Wolf HU, Lang W., and Zander R/: Alkaline haematin D – 575, a new tool for the determination of hemoglobin as an alternative to the Cyanhemoglobin method. Clin. Chem. Acta 136: 83 – 104 (1984).
11. National Committee for Clinical Laboratory Standards: Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture – third edition: approved standard. NCCLS Publication H4 – A3 Villanova, PA: NCCLS 1991.
12. Goldstein D.E., Little R.R, Wiedmeyer H.M., *et al*. Glycated Hemoglobin: Methodologies and clinical applications. Clin. Chem. 32: B64 – B70 (1986).
13. Ellis G., Diamonds E.P., Giesbrecht E.E. , *et al*: An automated “high-pressure” liquid chromatographic assay for hemoglobin A_{1c}. Clin. Chem 30: 1746 – 1752 (1984).
14. Burtis C.A. and Ashwood E.R. (eds.) : Tietz textbook of Clinical Chemistry, 2nd edition. W.B. Sanders Company, Philadelphia, PA, p.2021 (1994)
15. Knowles W.J., Haigh W.B., and Michaud G.C. : A monoclonal antibody-based immunoassay for hemoglobin A_{1c} Diabetes 35 Supplement 94A (1986).