BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B7201	200 / 600	R1	6 x	20 ml
		R2	6 x	8 ml
B7203	300 / 600*	R1	6 x	20 ml
		R2	6 x	8 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Immunturbidimetric assay for the in vitro quantitative determination of IgA in human serum and plasma.

Summary:

IgA accounts for 13% of the plasma immunoglobulins and serves to protect the skin and mucosa against microorganisms. It is capable of binding toxins, and develops in combination with lysozyme antibacterial and antiviral activity. IgA is the predominant immunoglobulin in bodily secretions such as colostum, saliva and sweat. Secretory IgA provides defense against local infections and is important in binding food antigens in the gut. In serum, IgA exists in monomeric, dimeric and trimetric forms, where as in bodily secretions it exists exclusively in dimetric form with an additional chain (secretory component).

Increased polyclonal IgA levels may occur in chronic liver diseases, chronic infections, autoimmune disorders (rheumatoid arthritis systemic lupus erythematosus), sarcoidosis and Wiscott-Aldrich syndrome. Monoclonal IgA increases in IgA myeloma. Decreased synthesis of IgA is observed in congenital and acquired immunodeficiency diseases such as Bruton type agammaglobulinemia. Reduced levels of IgA can be caused by protein-losing gastroenteropathies and loss through skin from burns.

Use of specific antibodies for quantitation of serum proteins has become a valuable diagnostic tool. Light-scattering properties of antigen/antibody aggregates were first observed by Pope and Hearnly in 1938, and later confirmed by Gitlin and Edelhoch. Ritchie employed turbidimetric measurements to quantitate specific proteins. Quantitation of immunoglobulins can also be done using nephelometric techniques. Polymeric enhancement with polyethylene glycol (PEG) to improve sensitivity and increase the rate of antigen/antibody complex formation has been described by Lizana and Hellsing. This IgA assay is based on the principle of immunological agglutination.

Test principle:

Immunturbidimetric assay
Anti-IgA antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity and reduces the risk of samples containing excess antigen producing false negative results.

Reagent concentration:

R1:
TRIS* buffer pH 7.8 100 mmol/l
NaCl 150 mmol/l
Polyethylene glycol 4.5 %
Preservative

R2:
Anti-human IgA antibody (goat): dependent on titer
TRIS* buffer pH 7.8 80 mmol/l
NaCl 150 mmol/l
Preservative
*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability at +2°C to +8°C: R1: 90 days
R2: 90 days

Specimen:

Collect serum using standard sampling tubes
Na-heparin-, Li-heparin- or EDTA-plasma
Stability: 7 days at +20°C to +25°C
3 months at +4°C to +8°C
6 months at - 20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to an index I of 90 (approximate conjugated and unconjugated bilirubin: 90 mg/dl)
Hemolysis: No significant interference up to an index H of 1200 (approximate hemoglobin concentration: 1200 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 700 (approximate triglycerides concentration: 1400 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
There is no cross-reaction between IgA, IgG and IgM under the assay conditions.
Sera with ill-defined characteristics should be subject to protein-electrophoresis to identify a possible antigen excess (e.g. gammopathy)
Rheumatoid factors < 100 IU/mL do not interfere.
A high-dose hook effect may occur at concentrations above 5000 mg/dl (50 g/l).
The effect of drugs or their metabolic products on the test is not known in all cases.
In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided
• Working solutions as described above
Additional materials required
• Calibrators and controls as indicated below
• 0.9% NaCl

Measuring/reportable range:

50 - 600 mg/dl (0.5 - 6 g/l)
Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution.

Expected values:

IFCC 70 – 400 mg/dl
CRM 470* 0.7 – 4.0 g/l

* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.
For diagnostic purposes the IgA results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 2 mg/dl (0.02 g/l)
The lower detection limit represents the lowest measurable IgA concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls (n = 20). The following results were obtained:

Run to run			
Sample	Mean (mg/dl)	SD(mg/dl)	CV %
Sample 1	218.48	4.627	2.1
Sample 2	242.19	4.874	2.0
Sample 3	439.48	11.349	2.6

Within run			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	115.00	1.346	1.2
Sample 2	228.37	3.395	1.5
Sample 3	342.94	13.581	4.0

Method comparison:

A comparison of the Analyticon Turbitex IgA (y) with a commercial obtainable assay (x) gave with 20 samples the following result (mg/dl).

$$y = 0.8649x + 36.281; \quad r = 0.9944$$

Quality control:

Protein Control Level 1 3 x 1 ml #7661
Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The IgA method was calibrated against the CRM 470 standard.

Calibration Type: Logit-Log
S1: 0.9% NaCl
S2-S6: Bio Cal® P Calibration Set 5 x 1 ml #1475

Calibration frequency

Full recalibration is recommended

- after lot change
- as required following quality control procedures

Text passages with grey background were changed in the latest revision of this package insert.

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