### Order information:

Catalog No.	Contents			
7781	R1	3 x	20 ml	R2 3 x 7 ml
H7201 Hit I / 917 (ILab* / AU*)	R1	6 x	20 ml	R2 6 x 8 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911: ACN 777  
Hitachi 917: ACN 737

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of IgA in human serum and plasma.

### Summary:

IgA accounts for 13% of the plasma immunoglobulins and serves to protect the skin and mucosa against microorganisms. It is capable of binding toxins, and develops in combination with lysozyme antibacterial and antiviral activity. IgA is the predominant immunoglobulin in bodily secretions such as colostrum, saliva and sweat. Secretory IgA provides defence against local infections and is important in binding food antigens in the gut. In serum, IgA exists in monomeric, dimeric and trimetric forms, where as in bodily secretions it exists exclusively in dimeric form with an additional chain (secretory component).

Increased polyclonal IgA levels may occur in chronic liver diseases, chronic infections, autoimmune disorders (rheumatoid arthritis systemic lupus erythematosus), sarcoidosis and Wiscott-Aldrich syndrome. Monoclonal IgA increases in IgA myeloma. Decreased synthesis of IgA is observed in congenital and acquired immunodeficiency diseases such as Bruton type agammaglobulinemia. Reduced levels of IgA can be caused by protein-losing gastroenteropathies and loss through skin from burns.

Use of specific antibodies for quantisation of serum proteins has become a valuable diagnostic tool. Light-scattering properties of antigen/antibody aggregates were first observed by Pope and Hearnly in 1938, and later confirmed by Gitlin and Edelhoch. Ritchie employed turbidimetric measurements to quantify specific proteins. Quantitation of immunoglobulins can also be done using nephelometric techniques. Polymeric enhancement with polyethylene glycol (PEG) to improve sensitivity and increase the rate of antigen/antibody complex formation has been described by Lizana and Hellsing. This IgA assay is based on the principle of immunological agglutination.

### Test principle:

Immunoturbidimetric assay

- Sample and addition of R1 (buffer)
- Addition of R2 (anti-IgA antibody/buffer) and start of reaction

Anti-IgA antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity and reduces the risk of samples containing excess antigen producing false negative results.

### Reagent concentration:

#### R1:

TRIS\* buffer pH 7.8 100 mmol/l  
NaCl 150 mmol/l  
Polyethylene glycol 4.5 %  
Preservative

#### R2:

Anti-human IgA antibody (goat) dependent on titer  
TRIS\* buffer pH 7.8 80 mmol/l  
NaCl 150 mmol/l  
Preservative

\*TRIS= Tris (hydroxymethyl)-aminomethane

### Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability at +2°C to +8°C: R1: 90 days  
R2: 90 days

### Specimen:

Collect serum using standard sampling tubes

Na-heparin-, Li-heparin- or EDTA-plasma

Stability: 7 days at +20°C to +25°C  
3 months at +4°C to +8°C  
6 months at - 20°C

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to a bilirubin concentration of 70 mg/dl.

Hemolysis: No significant interference up to a hemoglobin concentration of 1100 mg/dl.

Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.

There is no cross-reaction between IgA, IgG and IgM under the assay conditions.

Sera with ill-defined characteristics should be subject to protein-electrophoresis to identify a possible antigen excess (e.g. gammopathy)

A high-dose hook effect may occur at concentrations above 5000 mg/dl (50 g/l).

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

- Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

<b>Manual Procedure:</b>		
Wavelength:	600 nm	
Temperature:	+37°C	
Cuvette:	1 cm	
Zero adjustment:	against reagent blank	
Sample/Calibrator	Blank	Sample/Calibrator
R1	750 µl	12 µl 750 µl
Mix. Incubate 5 min. at +37°C and read A <sub>1</sub> . Then add:		
R2	300 µl	300 µl
Mix. Incubate 5 min. and read A <sub>2</sub> .		
<b>Calculation:</b>		
$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$		
The concentration of IgA in patient sera has to be calculated from $\Delta A$ using mathematic functions such as logit/log or can be read from a calibration curve obtained using 6 levels of standards in the concentration range of 0 to 800 mg/dl IgA. Saline solution (0.9% NaCl) is recommended for determining the zero point.		

### Measuring/reportable range:

50-800 mg/dl (0.50 – 8.00 g/l) maximum reportable range is dependent on the highest indicated standard concentration.

At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1 + 2). Multiply the result by the appropriate factor (e.g. 3) or determine with rerun function.

### Expected values:

IFCC 70 – 400 mg/dl  
CRM 470\* 0.7 – 4.0 g/l

\* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes the IgA results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 4 mg/dl (0.04 g/l)

The lower detection limit represents the lowest measurable IgA concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls between day (n = 20).

The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	141.9	2.41	1.70
Sample 2	278.7	3.96	1.42
Sample 3	412.5	4.57	1.11

Reproducibility within run was determined using controls (n = 20).

The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	134.0	1.11	0.83
Sample 2	169.5	2.20	1.30
Sample 3	268.2	3.81	1.42

### Method comparison:

A comparison of the Analyticon Turbitex IgA (y) with a commercial obtainable assay (x) gave with 40 samples the following result (mg/dl).

$$y = 1.023 x - 2.45; \quad r = 0.998$$

### Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Standardization: The IgA method was calibrated against the CRM 470 standard.

S1:	0.9% NaCl		
S2-S6:	Bio Cal <sup>®</sup> P	3 x 1 ml	#1470
	Bio Cal <sup>®</sup> P Set	5 X 1 ml	#1475

### Calibration frequency

Multicalibration

Full recalibration is recommended

- after lot change
- as required following quality control procedures

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