Turbitex® IgG-2

IMMUNOGLOBULIN G

Order information:

Catalog No.	Contents
7881	R1 3 x 20 ml R2 3 x 7 ml
H7151 Hit 1 / 917 (ILab* / AU*)	R1 6 x 20 ml R2 6 x 8 ml

(*) Kit contains only reagent barcodes for Hitachi systems.

System information:

Hitachi 911/917: ACN 674 serum/plasma (Standard application)
ACN 673 cerebrospinal fluid (Sensitive appl.)

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Please note: For measurement of CSF on Hitachi 911 and 912 analyzers please replace the reagent bottle labels by the special labels provided in the reagent kit.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of IgG in human serum and plasma and cerebrospinal fluid (CSF).

Summary:

IgG molecules are monomers composed of two light chains, (kappa or lambda), and two gamma heavy chains. Approximately 80% of serum immunoglobulin is IgG; its main task lies in defence against microorganisms, direct neutralisation of toxins and induction of complement fixation. IgG is the only immunoglobulin that can cross the placental barrier and provide passive immune protection for the foetus and newborn. This protection gradually catabolises until the infant's own immunological system begins at about six months of age. Near adult levels are reached by 18 months. Polyclonal IgG increases may be present in systemic lupus erythematosus, chronic liver diseases (infectious hepatitis and Laënnec's cirrhosis), infectious diseases, and cystic fibrosis. Monoclonal IgG increases in myeloma. Decreased synthesis of IgG is found in congenital and acquired immunodeficiency diseases and selective IgG subclass deficiencies, such as Bruton type agammaglobulinemia. Decreased IgG concentrations are seen in protein-losing enteropathies, nephrotic syndrome and through the skin from burns. Increased IgG metabolism is found in Wiskott-Aldrich syndrome, myotonic dystrophy and with anti-immunoglobulin antibodies.

The determination of IgG in cerebrospinal fluid (CSF) is used for evaluation of infections involving the central nervous system (CNS), neoplasms or primary neurologic diseases (in particular, multiple sclerosis). Increased CSF IgG concentrations may occur because of either increased permeability of the blood-brain barrier or local/intrathecal production of IgG, or both. Malfunction of the blood-brain barrier can be reliably quantified by means of the albumin CSF/serum ratio. An elevated albumin ratio is an indication of a disorder of the blood-brain barrier. If IgG and albumin are measured in CSF and serum simultaneously, differentiation between IgG originating from blood and IgG originating from intrathecal production is possible.

Use of specific antibodies for quantitation of serum proteins has become a valuable diagnostic tool. Light-scattering properties of antigen/antibody aggregates were first observed by Pope and Healey in 1938, and later confirmed by Gitlin and Edelhoch. Ritchie employed turbidimetric measurements to quantify specific proteins. Quantitation of immunoglobulins can also be done using nephelometric techniques. Polymeric enhancement with polyethylene glycol (PEG) to improve sensitivity and increase the rate of antigen/antibody complex formation has been described by Lizana and Hellsing.

This IgG assay is based on the principle of immunological agglutination.

Test principle:

Immunoturbidimetric assay

- Sample and addition of R1 (buffer)
- Addition of R2 (anti-IgG antibody/buffer) and start of reaction

Anti IgG antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically.

Reagent concentration:

R1:

Tris buffer* pH 7.8 200 mmol/l
NaCl 85 mmol/l
Polyethylene glycol 3,5 %
Preservative <0,1 %

R2:

Anti-human IgG antibody (goat); dependent on titer
TRIS buffer* pH 7.8 200 mmol/l
NaCl 400 mmol/l
Preservative <0,1 %

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at 2-8°C

Onboard stability at 2-8°C: R1: 90 days
R2: 90 days

Specimen:

Collect serum using standard sampling tubes.

Na-/Li-heparin or EDTA-plasma

Stability: at +15°C - +25°C 4 months
at +2°C - +8°C 8 months
at -20°C 8 months

For analysis on the Roche/Hitachi 902, manually dilute each specimen 1+15 (i.e. 1:16) with 0.9% NaCl.

CSF samples should be as fresh as possible and used undiluted.

Stability: at +15°C - +25°C 1 day
at +2°C - +8°C 7 days
at -20°C not recommended

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

The IgG concentration in serum is much higher than in CSF samples. To achieve high sensitivity the CSF application (IGGC2) is run with 30 µL of undiluted sample. Therefore this application is prone to sample carry-over of IgG. To avoid sample carry-over it is recommended to wash sample probe with 100 µL 1 M NaOH/Multi-clean using the special wash/sample probe evasion option in the software or on the wash-rack. If possible, use batch mode.

Criterion: Recovery within ±10% of initial value.

Standard application (IGG-2):

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl or 1026 µmol/l).

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dl or 621 µmol/l).

Lipemia (Intralipid): No significant interference up to an L index of 1500. (approximate triglycerides concentration: 3000 mg/dl)

There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors < 1200 IU/ml do not interfere.

Falsely low results due to antigen excess may occur at 440 g/l (44000 mg/dl) in polyclonal specimens. There is no cross-reaction between IgG and IgA or IgM under the assay conditions.

As with other turbidimetric or nephelometric procedures, this test may not provide accurate results in patients with monoclonal gammopathy, due to individual sample characteristics which can be assessed by electrophoresis.

Sensitive application (IGGC2):

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl or 1026 µmol/l).

Hemolysis: No significant interference up to an H index of 800 (approximate hemoglobin concentration: 800 mg/dl or 497 µmol/l).

Falsely low results due to antigen excess may occur above 1400 mg/L.

There is no cross-reaction between IgG and IgA or IgM under the assay conditions.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0,9% NaCl

Manual Testing procedure:		
Wavelength:	340 nm	
Temperature:	37°C	
Cuvette:	1 cm light path	
Zero adjustment:	Aqua dest.	
Serum/Plasma		Liquor
R1	750 µl	750 µl
Sample/Calibrator	21 µl	30 µl
Mix. Incubate for 5 min. at 37°C and read A ₁ . Then add:		
R2	255 µl	60 µl
Aqua dest.	---	90 µl
Mix. Incubate for 5 min. at 37°C and read A ₂ .		
Calculation:		
$\Delta A = [(A_2 - A_1) \text{ Sample or Calibrator}] - [(A_2 - A_1) \text{ Aqua dest.}]$		
The concentration of IgG in patient sera has to be calculated from ΔA using mathematic functions such as logit/log or can be read from a calibration curve obtained using 6 levels of standards in the concentration range of 0 to 500 mg/dl IgG. Saline solution (0.9% NaCl) is recommended for determining the zero point.		

Measuring/reportable range:

Standard application: 300-5000 mg/dl
Sensitve application (CSF): 0.2-20 mg/dl (2.0-200mg/L)

At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1+1). Multiply the result by the appropriate factor (e.g. 2) or determine with rerun-function.

Expected values:

Reference values according to CRM 470 Protein Standardization
Serum: Adults 700-1600 mg/dl
CSF: Adults 10-30 mg/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the IgG results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Serum: 16 mg/dl (0,16 g/l)
CSF: 0,1 mg/dl (1 mg/l)

The lower detection limit represents the lowest measurable IgG concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

Imprecision:

Reproducibility was determined using controls between day (n = 20). The following results were obtained:

Standard application:

Sample	Within run			Run/Run		
	Mean mg/dl	SD mg/dl	CV %	Mean mg/dl	SD mg/dl	CV %
Sample 1	746.04	7.18	0.96	205	12.5	6.1
Sample 2	582.43	11.03	1.89	1032	14.86	1.44
Sample 3	1241.7	14.24	1.15	2545	43.77	1.72

Sensitive application:

Sample	Within run			Run/Run		
	Mean mg/L	SD mg/L	CV %	Mean mg/L	SD mg/L	CV %
Sample 1	43.73	0.7	1.61	13.21	0.86	6.52
Sample 2	154.51	2.58	1.67	327.7	6.78	2.07
Sample 3	85.31	0.84	0.99	166.6	2.12	1.27

Method comparison:

Standard application:

A comparison of the Analyticon Turbitex IgG-2 assay (x) with a commercially available assay (y) with 103 human sera samples gave the following correlation (g/L):

$$y = 1.075x - 0.977 \quad r = 0.996$$

Sensitive application (IGGC2):

A comparison of the Analyticon Turbitex IgG-2 assay (x) with a commercially available assay (y) with 54 human liquor samples gave the following correlation (mg/L):

$$y = 1.042x - 3.903 \quad r = 0.988$$

Calibration:

Standardisation: The IgG method has been standardised against the reference preparation CRM 470.

Standard application (IgG-2):

Roche/Hitachi 904, 911, 912, 917:

S1: 0.9% NaCl

S2-S6: BioCal P 3 x 1 ml #1470

Multiply the lot-specific calibrator value by the factors below to determine the standard concentrations for the 6-point calibration curve.

S2: 0.10 S4: 0.50 S6: 2.50
S3: 0.25 S5: 1.00

Sensitive application (IgG-2 CSF):

Roche/Hitachi 904, 911, 912, 917:

S1: 0.9% NaCl

S2-S6: C.f.a.s. PUC from Roche

Multiply the lot-specific C.f.a.s. PUC calibrator value by the factors below to determine the standard concentrations for the 6-point calibration curve.

S2: 0.042 S4: 0.166 S6: 1.000
S3: 0.083 S5: 0.333

Calibration frequency

Full recalibration is recommended:

- after lot change
- as required following quality control procedures

Quality control:

Standard application:

Protein Control Level 1 3 x 1 ml #7661

Protein Control Level 2 3 x 1 ml #7662

Sensitive application:

For quality control, use Precinorm PUC and Precipath PUC from Roche or Liquicheck Spinal Fluid Control from Biorad. Other suitable control material can be used in addition.

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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- Kaplan LA, Pesce AJ, eds. Clinical Chemistry, Theory, Analysis and Correlation, 3. Auflage. Mosby Inc, 1996.