BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B7301	200 / 600	R1	6 x	20 ml
		R2	6 x	8 ml
B7303	300 / 600*	R1	6 x	20 ml
		R2	6 x	8 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Immunturbidimetric assay for the in vitro quantitative determination of IgM in human serum and plasma.

Summary:

IgM normally consists of 10 heavy μ -chains and 10 κ or λ type light chains which are always identical within a molecule. There is also a J-chain linking all the μ -chains together, so that simply speaking, IgM has a pentameric structure when compared to that of IgG. IgM is the largest immunoglobulin molecule (MW=900000), but makes up only 6% of the plasma immunoglobulins.

IgM is the first specific antibody to appear in the serum after infection. It is capable of activating complement, thus helping to kill bacteria. IgM levels sink after the infection has subsided at a relatively rapid rate compared to IgG. This fact is used to advantage in the differential diagnosis of acute and chronic infections by comparing specific IgM and IgG titers. If IgM is prevalent the infection is acute, whereas if IgG predominates the infection is chronic (e.g. rubella, viral hepatitis).

Increased polyclonal IgM levels are found in viral, bacterial, and parasitic infections, liver diseases, rheumatoid arthritis, scleroderma, cystic fibrosis, and heroin addiction. Monoclonal IgM is increased in Waldenström's macroglobulinemia. Increased loss of IgM is found in protein-losing enteropathies and in burns. Decreased synthesis of IgM occurs in congenital and acquired immunodeficiency syndromes.

Use of specific antibodies for quantitation of serum proteins has become a valuable diagnostic tool. Light-scattering properties of antigen/antibody aggregates were first observed by Pope and Healey in 1938, and later confirmed by Gitlin and Edelhoch. Ritchie employed turbidimetric measurements to quantitate specific proteins. Quantitation of immunoglobulins can also be done using nephelometric techniques. Polymeric enhancement with polyethylene glycol (PEG) to improve sensitivity and increase the rate of antigen/antibody complex formation has been described by Lizana and Hellsing.

Test principle:

Immunturbidimetric assay.

Anti-IgM antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity, and reduces the risk of samples containing excess antigen producing false negative results.

Reagent concentration:

R1:

TRIS* buffer pH 7.8 100 mmol/l
NaCl 150 mmol/l
polyethylene glycol 4.5 %
preservative

R2:

Anti-human IgM antibody (goat) dependent on titer
TRIS* buffer pH 7.8 80 mmol/l
NaCl 150 mmol/l
preservative

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability at +2°C to +8°C: R1: 90 days
R2: 90 days

Specimen:

Collect serum using standard sampling tubes

Na-heparin-, Li-heparin- or EDTA-plasma

Stability: 7 days at +20°C to +25°C
3 month at +4°C to +8°C
6 month at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.

Icterus: No significant interference up to a bilirubin concentration of 100 mg/dl.

Hemolysis: No significant interference up to a hemoglobin concentration of 1200 mg/dl.

Lipemia (Intralipid): No significant interference up to an index L of 675 (approximate triglycerides concentration: 1350 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factors < 56 IU/mL do not interfere.

There is no cross-reaction between IgM, IgG and IgA under the assay conditions.

Sera with ill-defined characteristics should be subject to protein-electrophoresis to identify a possible antigen excess (e.g. gammopathy)

A high-dose hook effect may occur at concentrations above 10000 mg/dl IgM.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring/reportable range:

20 – 800 mg/dl (0.2 – 8 g/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution.

Expected values:

IFCC 40 – 230 mg/dl
CRM 470* 0.4 – 2.3 g/l

* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes the IgM results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 2 mg/dl (0.02 g/l)

The lower detection limit represents the lowest measurable IgM concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls (n = 20). The following results were obtained:

<u>Within run</u>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	107.13	2.083	1.9
Sample 2	213.40	2.398	1.1
Sample 3	315.05	7.634	2.4

<u>Between day</u>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	102.65	2.178	2.1
Sample 2	86.38	2.432	2.8
Sample 3	156.64	3.259	2.1

Method comparison:

A comparison of the Analyticon Turbitex IgM (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 0.9808x - 2.2387; \quad r = 0.9953$$

Quality control:

Protein Control Level 1 3 x 1 ml #7661
Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The IgM method was calibrated against the CRM 470 reference material.

Calibration Type: Logit-Log

S1: 0.9% NaCl

S2-S6: Bio Cal® P Calibration Set 5 x 1 ml #1475

Calibration frequency:

Full recalibration is recommended

- after lot change
- as required following quality control procedures

Text passages with grey background were changed in the latest revision of this package insert.

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