Order information:

Catalog No.	Contents
7981	R1 3 x 20 ml R2 3 x 7 ml
H7301 Hit I / 917 (ILab*/AU*)	R1 6 x 20 ml R2 6 x 8 ml

(*) Kit contains only reagent barcodes for Hitachi systems.

System information:

Hitachi 911/917: ACN 789
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of IgM in human serum and plasma.

Summary:

IgM normally consists of 10 heavy μ -chains and 10 κ or λ type light chains which are always identical within a molecule. There is also a J-chain linking all the μ -chains together, so that simply speaking, IgM has a pentameric structure when compared to that of IgG. IgM is the largest immunoglobulin molecule (MW=900000), but makes up only 6% of the plasma immunoglobulins.

IgM is the first specific antibody to appear in the serum after infection. It is capable of activating complement, thus helping to kill bacteria. IgM levels sink after the infection has subsided at a relatively rapid rate compared to IgG. This fact is used to advantage in the differential diagnosis of acute and chronic infections by comparing specific IgM and IgG titers. If IgM is prevalent the infection is acute, whereas if IgG predominates the infection is chronic (e.g. rubella, viral hepatitis).

Increased polyclonal IgM levels are found in viral, bacterial, and parasitic infections, liver diseases, rheumatoid arthritis, scleroderma, cystic fibrosis, and heroin addiction. Monoclonal IgM is increased in Waldenstroem's macroglobulinemia. Increased loss of IgM is found in protein-losing enteropathies and in burns. Decreased synthesis of IgM occurs in congenital and acquired immunodeficiency syndromes.

Use of specific antibodies for quantisation of serum proteins has become a valuable diagnostic tool. Light-scattering properties of antigen/antibody aggregates were first observed by Pope and Healey in 1938, and later confirmed by Gitlin and Edelhoch. Ritchie employed turbidimetric measurements to quantify specific proteins. Quantitation of immunoglobulins can also be done using nephelometric techniques. Polymeric enhancement with polyethylene glycol (PEG) to improve sensitivity and increase the rate of antigen/antibody complex formation has been described by Lizana and Hellsing.

This IgA assay is based on the principle of immunological agglutination.

Test principle:

Immunoturbidimetric assay
• Sample and addition of R1 (buffer)
• Addition of R2 (anti-IgM antibody/buffer) and start of reaction

Anti-IgM antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity, and reduces the risk of samples containing excess antigen producing false negative results.

Reagent concentration:

R1:
TRIS* buffer pH 7.8 100 mmol/l
NaCl 150 mmol/l
Polyethylene glycol 4.5 %
Preservative

R2:
Anti-human IgM antibody (goat) dependent on titer
TRIS* buffer pH 7.8 80 mmol/l
NaCl 150 mmol/l
Preservative

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability at +2°C to +8°C: R1: 90 days
R2: 90 days

Specimen:

Collect serum using standard sampling tubes
Na-heparin-, Li-heparin- or EDTA-plasma
Stability: 7 days at +20°C to +25°C
3 month at +4°C to +8°C
6 month at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.
Icterus: No significant interference up to a bilirubin concentration of 80 mg/dl).
Hemolysis: No significant interference up to a hemoglobin concentration of 11000 mg/dl.
Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.
There is no cross-reaction between IgM, IgG and IgA under the assay conditions.
Sera with ill-defined characteristics should be subject to protein-electrophoresis to identify a possible antigen excess (e.g. gammopathy).
A high-dose hook effect may occur at concentrations above 10000 mg/dl IgM.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:		
Wavelength:	340 nm	
Temperature:	+37°C	
Cuvette:	1 cm	
Zero adjustment:	against reagent blank	
	Blank	Sample/Calibrator
Sample/Calibrator	---	10 μ l
R1	750 μ l	750 μ l
Mix. Incubate 5 min. at 37°C and read A ₁ . Then add:		
R2	250 μ l	250 μ l
Mix. Incubate 5 min. and read A ₂ .		
Calculation:		
$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$		
The concentration of IgM in patient sera has to be calculated from ΔA using mathematic functions such as logit/log or can be read from a calibration curve obtained using 6 levels of standards in the concentration range of 0 to 490 mg/dl IgM. Saline solution (0.9% NaCl) is recommended for determining the zero point.		

Measuring/reportable range:

30 – 850 mg/dl (0.30 – 8.50 g/l)*
At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1+1). Multiply the result by the appropriate factor (e.g. 2) or determine with rerun function.

*The maximum reportable range is dependent on the highest indicated standard concentration.

Expected values:

IFCC 40 – 230 mg/dl
CRM 470* 0.4 – 2.3 g/l

* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the IgM results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 5 mg/dl (0.05 g/l)
The lower detection limit represents the lowest measurable IgM concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls between day (n = 20). The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	86.74	2.61	3.01
Sample 2	167.23	4.51	2.70
Sample 3	246.43	5.49	2.23

Method comparison:

A comparison of the Analyticon Turbitex IgM (y) with a commercial obtainable assay (x) gave the following result (mg/dl):
 $y = 0.960 x + 6,057$; $r = 0.995$

Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The IgM method was calibrated against the CRM 470 reference material.

S1:	0.9% NaCl		
S2-S6:	Bio Cal [®] P	3 x 1 ml	#1470
	Bio Cal [®] P Set	5 X 1 ml	#1475

Calibration frequency

Multicalibration

Full recalibration is recommended

- after lot change
- as required following quality control procedures

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