Order information:

Catalog No.	Contents						
9721	R1	2 x	20 ml	R2	2 x	15 ml	
H8401	Hit I (ILab*)	R1	6 x	47 ml	R2	6 x	30 ml
H8403	Hit 917 (AU*)	R1	6 x	20 ml	R2	6 x	13 ml
AU8404	AU	R1	6 x	60 ml	R2	6 x	40 ml

(*) Kit contains only reagent barcodes for Hitachi systems

System information:

Hitachi 911/917: ACN 731
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:


Enzymatic assay for the in vitro quantitative determination of Lipase in human serum and plasma.

Summary:

Lipases are glycoproteins with a molecular weight of 47000 Daltons. They are defined as triglyceride hydrolases which catalyze the cleavage of triglycerides to diglycerides with subsequent formation of monoglycerides and fatty acids. In addition to a-amylase, pancreatic lipases have for many years been undeniably the most important clinical chemistry parameters for the differential diagnosis of diseases of the pancreas. The lipase activity determination has gained increasing international recognition because of its high specificity and rapid response. After acute pancreatitis the lipase activity increases within 4-8 hours, reaches a peak after 24 hours and decreases after 8 to 14 days. However, there is no correlation between the lipase activity determined in serum and the extent of damage to the pancreas. Numerous methods have been described for the determination of lipase which determine the decrease in substrate turbidimetrically or nephelometrically or determine degradation products. This method is based on the cleavage of a specific chromogenic lipase substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester emulsified with bile acids. The pancreatic enzyme activity is determined specifically by the combination of bile acid and colipase used in this assay. Virtually no lipase activity is detected in the absence of co-lipase. Colipase only activates pancreatic lipase, but not other lipolytic enzymes found in serum. The high amount of cholates ensure that the esterases present in the serum do not react with the chromogenic substrate due to the highly negative surface charge.

Test principle:

Enzymatic colorimetric assay
• Sample and addition of R1 (buffer/colipase/cholate)
• Addition of R2 (emulsion/chromogenic substrate/cholate) and start of reaction:

The chromogenic lipase substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester is cleaved by the catalytic action of alkaline lipase solution to form 1,2-O-dilauryl-rac-glycerol and an unstable intermediate, glutaric acid-(6-methylresorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The colour intensity of the red dye formed is directly proportional to the lipase activity and can be determined photometrically.

Reagent concentration:

R1:	
Buffer/colipase/cholate BICIN buffer*: pH 8.0	50 mmol/l
colipase (porcine pancreas)	1 mg/l
Na-deoxycholate	1.6 mmol/l
calcium chloride	10 mmol/l
detergent	
preservative	
R2:	
Emulsion/chromogenic substrate/cholate	
Tartrate buffer pH 4.0	10 mmol/l
1,2-O-dilauryl-rac-glycero-3 glutaric acid-(6-methylresorufin) ester	0.27 mmol/l
taurodeoxycholate	8.8 mmol/l
detergent	
preservative	
*BICIN = N,N-bis(2-hydroxyethyl)-glycine	

Preparation and Stability:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability:
R1: 28 days opened and refrigerated on the Hitachi-analyser
R2: 28 days opened and refrigerated on the Hitachi-analyser
Store protected from light after opening.

Specimen:

Collect serum using standard sampling tubes
Li-, Na- or NH₄- heparin plasma.
Stability: 7 days at +20°C to +25°C
7 days at +4°C to +8°C
1 year at -20°C

EDTA-, oxalate-, fluoride- or citrated plasma lead to decreased results (inhibition of lipase activity). Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial values.
(Only Hitachi): Use the evasion software with SMS for the individual assays in order to avoid cuvette carry-over of triglycerides, cholesterol, HDL cholesterol plus and LDL-cholesterol plus to lipase and probe carry-over of triglycerides and cholesterol to lipase as well as from lipase to calcium.
Icterus: No significant interference up to an index I of 60 (approximate bilirubin concentration: 60 mg/dl).
Hemolysis: No significant interference up to an index H of 500 (approximate haemoglobin concentration: 500 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl).
There is poor correlation between turbidity and triglycerides concentration. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Procedure:	
Wavelength:	570 nm
Temperature:	+37°C
Cuvette:	1 cm
Zero adjustment:	against air
Sample/Calibrator	12 µl
R1	750 µl
Mix and incubate for 5 min. Then add:	
R2	450 µl
Mix and read initial absorbance A ₁ . After 2 minutes read absorbance A ₂ .	
Calculation:	
$\frac{A_2 \text{ sample} - A_1 \text{ sample}}{A_2 \text{ calibrator} - A_1 \text{ calibrator}} \times \text{calibrator conc.} = \text{Lipase conc.}$	

Measuring /reportable range:

Roche/Hitachi 911/917:
Measuring range: 3 - 300 U/l (0.05-5.00 µkat/l)
Determine samples with lipase activity > 300 U/l via the rerun function. Manually dilute the original samples having activities > 3300 U/l (55.00 µkat/l) with 0.9% NaCl or distilled/deionized water (e.g. 1 + 50). Multiply the result by the appropriate dilution factor (e.g. 51).

Provisional expected values:

Adults: < 60 U/l (<1.00 µkat/l)
In order to relate findings to the lipase reference range obtained using the turbidimetric method, multiply the results (U/l or µkat/l) by the factor 3.2 to obtain an approximate comparison.
Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

Analytical sensitivity (lower detection limit):

Detection limit: 3 U/l (0.05 µkat/l)
The lower detection limit represents the lowest measurable lipase activity that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Sample	Within run		
	Mean U/l	SD U/l	CV %
Control serum 1	34	0.39	1.16
Control serum 2	71	0.51	0.72
Control serum 3	58	0.39	0.67

Sample	day/day		
	Mean U/l	SD U/l	CV %
Control serum 1	35	0.23	0.65
Control serum 2	63	0.55	0.87
Control serum 3	72	0.90	1.26

Method comparison:

A comparison of the Analyticon Fluitest LIP Col. (y) with a turbidimetric assay gives the following result:

$$y = 0.304x + 6.563; \quad r = 0.986$$

Quality control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Standardization: The method was standardized against a nephelometric method.

Calibration frequency:

Full recalibration is recommended

- after lot change
 - as required following quality control procedures
- Calibration verification: Not necessary.

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