BioLyzer® Order information:

Catalog No.	BioLyzer	Contents
B7851	200 / 600	R1 6 x 20 ml
B7853	300 / 600*	R1 6 x 20 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

For the quantitative determination of lactate in plasma or cerebrospinal fluid.

Summary:


Anaerobic glycolysis markedly increases blood lactate and causes some increase in pyruvate levels, especially with prolonged exercise. The common cause for increased blood lactate and pyruvate is anoxia resulting from such conditions as shock, pneumonia and congestive heart failure. Lactic acidosis may also occur in renal failure and leukaemia. Thiamine deficiency and diabetic ketoacidosis are associated with increased levels of lactate and pyruvate.

Lactate levels in cerebrospinal fluid (CSF) are increased in bacterial meningitis. Increased (CSF) levels also occur in hypocapnia, hydrocephalus, brain abscesses, cerebral ischemia and any clinical condition associated with reduced oxygenation of the brain and / or increased intracranial pressure.

The method presented here uses an enzymatic reaction to convert lactate to pyruvate. The hydrogen peroxide produced by this reaction is then used in an enzymatic reaction to generate a colored dye.

Test principle:

Lactate-oxidase (LOD) cleaves lactate into pyruvate and H₂O₂, which reacts in the presence of peroxidase (POD) with 4-aminoantipyrine and TBHB to a red quinoneimin dye. The increase of color is proportional to the lactate concentration.

Reagent concentration:

R1:	
Tris buffer pH 7.5	150 mmol/l
Lactate oxidase (LOD)	> 0.3 kU/l
Peroxidase (POD)	> 1.0 kU/l
4-aminoantipyrine	0.3 mmol/l
TBHB*	2.5 mmol/l

*Tribrom-3-hydroxybenzoic acid

Preparation and stability:

R1: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability: R1: 60 days

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Test results were obtained on a **BioLyzer® 300** analyzer.

Specimen:

- Do not use serum.
- CSF or Plasma with one of the following anticoagulants: Potassium EDTA, Fluoride EDTA, Fluoride Heparin, or Fluoride Oxalate.
- Use plasma from blood collected by standard venipuncture technique. Centrifuge within 15 minutes of drawing specimen.
- Cerebrospinal fluid (CSF) may be used as obtained.
- Blood samples should only be drawn prior to the administration of Metaminizole.

Please note:

1. The lactate level increases during rapidly physical exercise. The time required for return to normal lactate values depends on the physical fitness of the subject. Thirty minutes at rest is usually sufficient for this purpose.
2. Blood samples should be drawn from a stasis-free vein. However, minimal hemostasis (less than 30 seconds) will not affect lactate levels. Avoid the use of a tourniquet, if possible.
3. Glycolysis in blood samples can rapidly increase lactate levels. Cells contribute to the glycolysis and their quick removal is essential for accurate lactate analysis. Heparinized plasma is acceptable, but precautions must be taken to retard glycolysis by keeping the whole blood on ice and then separating the plasma from the cells within 15 minutes of collection.

Stability:

Plasma	2 hours	at +20°C to +25°C
	24 hours*	at +2°C to +8°C
CSF	3 hours	at +20°C to +25°C
	24 hours	at +2°C to +8°C
	1 month	at -20°C

* Determine Potassium EDTA Plasma within 2 hours.

Centrifuge samples containing precipitates before performing the assay.

Limitations - interference:

Criterion: Recovery within ± 10% of initial values.

Plasma

Ascorbic acid: Interference at all concentrations.

Icterus: No significant interference up to an index I of 12.5 (approximate conjugated and unconjugated bilirubin concentration 12.5 mg/dl)

Hemolysis: No significant interference up to an Index H of 125 (approximate hemoglobin concentration: 125 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 67.5 (approximate triglycerides concentration: 125 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metamizole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and Controls as indicated below
- 0.9% NaCl

Conversion into SI-units: Lactate [mg/dl] x 0.11 = [mmol/l]

Measuring / reportable range:

1.0 – 200 mg/dl (0.11 – 22.2 mmol/l)*

Determine samples with higher concentrations via the rerun function using 0.9% NaCl as diluent.

* BioLyzer 200: 1.0 – 150 mg/dl (0.11 – 16.5 mmol/l)

Analytical sensitivity (lower detection limit):

0.2 mg/dl (0.02 mmol/l)

The lower detection limit represents the lowest measurable lactate concentration that can be distinguished from zero.

Expected values:

According to Tietz:

Plasma:	4.5 – 19.8 mg/dl	0.5 – 2.2 mmol/l	venous
	4.5 - 14.4 mg/dl	0.5 – 1.6 mmol/l	arterial
CSF:	10 – 60 mg/dl	1.1 – 6.7 mmol/l	neonate
	10 – 40 mg/dl	1.1 – 4.4 mmol/l	3-10 days old
	10 – 25 mg/dl	1.1 – 2.8 mmol/l	> 10 days old
	10 – 22 mg/dl	1.1 – 2.4 mmol/l	adult

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the lactate results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Imprecision:

Reproducibility was determined using controls between day (n = 20) and within run (n = 20). The following results were obtained:

Plasma	Within run			
	Sample	Mean mg/dl	SD mg/dl	CV %
	Sample 1	14.1	0.103	0.7
	Sample 2	40.8	0.282	0.7
	Sample 3	100.3	0.852	0.8

Plasma	Run to run			
	Sample	Mean mg/dl	SD mg/dl	CV %
	Sample 1	13.9	0.182	1.3
	Sample 2	40.6	0.890	2.2
	Sample 3	101.9	0.979	1.0

Method comparison:

A comparison of the Analyticon Fluitest[®] Lactate determination (y) with another commercial available assay (x) gave the following correlation (mg/dl):
 $y = 1.3344x + 1.891$ $r^2 = 0.9953$

Quality Control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio-Cal [®] E	10 x 3 ml	#1430
Bio Cal [®]	20 x 3 ml	#1420

Calibration frequency:

Two point calibration is recommended:

- After lot change
- as required following quality control procedures

Calibration verification: Not necessary

Disposal:

Please note the legal regulations.

Literature:

1. Barhan D, Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system. *Analyst* 97.1972:142.
2. Bergmeyer HU (ed), New York, NY: Academic Press Inc; 1974;1464.
3. Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 4th ed. St.Louis. Elsevier Saunders; 2006.
4. Gutmann I, Wahlefeld A. *Methods of Enzymatic Analysis*. 2nd ed.
5. Noll F. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU ed.
6. Tietz NW. *clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382 – 383.
7. Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. *Ann. Clin. Biochem.* 1969;6:24.
8. Westgard JO, Lahmeyer BL, Birnbaum ML. *Clin Chem.* 1972; 18:1334-1338.
9. Heil W, Rick W. Lactatbestimmung im Liquor cerebrospinalis zur Differenzierung zwischen bakterieller und viraler Meningitis. *Internist* 1984; 25:320-324.

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