Order information:

Catalog No.	Contents
3011	R1 6 x 20 ml R4 1 x 5 ml
H7851 Hit I (ILab*/AU*)	R1 6 x 20 ml
AU7851 AU	R1 6 x 20 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 361 – 400 (user-defined method)
 Hitachi 917: ACN 901 – 905 (user-defined method)
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

For the quantitative determination of lactate in plasma or cerebrospinal fluid.

Summary:

Anaerobic glycolysis markedly increases blood lactate and causes some increase in pyruvate levels, especially with prolonged exercise. The common cause for increased blood lactate and pyruvate is anoxia resulting from such conditions as shock, pneumonia and congestive heart failure. Lactic acidosis may also occur in renal failure and leukaemia. Thiamine deficiency and diabetic ketoacidosis are associated with increased levels of lactate and pyruvate. Lactate levels in cerebrospinal fluid (CSF) are increased in bacterial meningitis. Increased (CSF) levels also occur in hypocalcaemia, hydrocephalus, brain abscesses, cerebral ischemia and any clinical condition associated with reduced oxygenation of the brain and / or increased intracranial pressure. Lactate measurements that evaluate the acid-base status are used in the diagnosis and treatment of lactic acidosis (abnormally high acidity in the blood). In recent years, enzymatic methods for the determination of lactate have gained favour over colorimetric and titrimetric methods. Enzymatic methods are generally simple and provide greater specificity, accuracy, and reproducibility. The first enzymatic method described for the determination of lactate was based on the transfer of hydrogen from lactate to potassium ferricyanide by lactate dehydrogenase. However, the procedure was cumbersome and did not receive wide acceptance. Subsequent methods involved the UV measurement of the formation of NADH. In 1974, Gutmann and Wahlefeld described a lactate procedure that measures the NADH formed by the oxidation of lactate catalysed by LD, using hydrazine as a trapping agent for pyruvate. A method described by Noll is also based on the catalytic action of LD but includes ALT in the reaction mixture to more rapidly remove the pyruvate formed from the conversion of lactate. The method presented here uses an enzymatic reaction to convert lactate to pyruvate. The hydrogen peroxide produced by this reaction is then used in an enzymatic reaction to generate a colored dye. This method offers longer reagent stability than the previous UV enzymatic methods.

Test principle:

Lactate-oxidase (LOX) cleaves lactate into pyruvate and H₂O₂, which reacts in the presence of peroxidase (POD) with 4-aminoantipyrine and TBHB to a red quinoneimin dye. The increase of color is proportional to the lactate concentration.

Reagent concentration:

R1:	
Tris buffer pH 7.5	150 mmol/l
Lactate oxidase (LOD)	> 0.3 kU/l
Peroxidase (POD)	> 1.0 kU/l
4-aminoantipyrine	0.3 mmol/l
TBHB*	2.5 mmol/l

R4: (#3011)	
Lactate	30 mg/dl

*Tribrom-3-hydroxybenzoic acid

Specimen:

- Do not use serum.
- CSF or Plasma with one of the following anticoagulants: Potassium EDTA, Fluoride EDTA, Fluoride Heparin, or Fluoride Oxalate.
- Use plasma from blood collected by standard venipuncture technique. Centrifuge within 15 minutes of drawing specimen.
- Cerebrospinal fluid (CSF) may be used as obtained.
- Blood samples should only be drawn prior to the administration of Metaminizole.

Please note:

- The lactate level increases rapidly during physical exercise. The time required for return to normal lactate values depends on the physical fitness of the subject. Thirty minutes at rest is usually sufficient for this purpose.
- Blood samples should be drawn from a stasis-free vein. However, minimal hemostasis (less than 30 seconds) will not affect lactate levels. Avoid the use of a tourniquet, if possible.

- Glycolysis in blood samples can rapidly increase lactate levels. Cells contribute to the glycolysis and their quick removal is essential for accurate lactate analysis. Heparinized plasma is acceptable, but precautions must be taken to retard glycolysis by keeping the whole blood on ice and then separating the plasma from the cells within 15 minutes of collection.

Stability:

Plasma	2 hours	at +20°C to +25°C
	24 hours*	at +2°C to +8°C
CSF	3 hours	at +20°C to +25°C
	24 hours	at +2°C to +8°C
	1 month	at -20°C

* Determine Potassium EDTA Plasma within 2 hours.
 Centrifuge samples containing precipitates before performing the assay.

Preparation and stability:

R1: Ready for use. R4: Ready for use.
 Unopened kit components: Up to the expiration date at +2°C to +8°C
 Onboard stability: R1: 60 days at +2°C to +8°C

Notes:

For in vitro diagnostic use.
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
 Exercise the normal precautions required for handling all laboratory reagents.
 Test results were obtained on a Hitachi 911 analyzer.

Limitations - interference:

Criterion: Recovery within ± 10% of initial values.

Plasma

Ascorbic acid: Interference at all concentrations.
 Icterus: No significant interference up to an index I of 12.5 (approximate conjugated and unconjugated bilirubin concentration 12.5 mg/dl)
 Hemolysis: No significant interference up to an Index H of 250 (approximate hemoglobin concentration: 250 mg/dl).
 Lipemia (Intralipid): Interference at all concentrations. There is poor correlation between turbidity and triglycerides concentration.
 False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminizole (Dipyron). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and Controls as indicated below
- 0.9% NaCl

Manual Procedure:		
Wavelength:	Hg 546 nm (540 nm)	
Temperature:	+37°C	
Cuvette:	1 cm	
Zero adjustment:	against reagent blank	
Sample/Calibrator	Blank	Sample/Calibrator
R1	1000 µl	1000 µl
Mix and incubate 5 min. at +37°C. Read absorbance of standard and sample within 30min against reagent blank.		
Calculation:		
ΔA_{sample}	x conc. Calibrator = Lactate [mg/dl]	
$\Delta A_{\text{calibrator}}$		

Conversion into SI-units: Lactate [mg/dl] x 0,111 = [mmol/l]

Measuring /reportable range:

1.0 – 150 mg/dl (0.11 – 16.65 mmol/l)
 Samples with extremely high lactate concentrations may give falsely low results in the measuring range.
 Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1+3). Multiply the result by the appropriate dilution factor (e.g. factor 4).

Analytical sensitivity (lower detection limit):

0.37 mg/dl (0.041 mmol/l)
 The lower detection limit represents the lowest measurable lactate concentration that can be distinguished from zero.

Expected values:

According to Tietz:

Plasma:	4.5 – 19.8 mg/dl	0.5 – 2.2 mmol/l	venous
	4.5 - 14.4 mg/dl	0.5 – 1.6 mmol/l	arterial
CSF:	10 – 60 mg/dl	1.1 – 6.7 mmol/l	neonate
	10 – 40 mg/dl	1.1 – 4.4 mmol/l	3-10 days old
	10 – 25 mg/dl	1.1 – 2.8 mmol/l	> 10 days old
	10 – 22 mg/dl	1.1 – 2.4 mmol/l	adult

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the lactate results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Imprecision:

Reproducibility was determined using controls between day (n = 20) and within run (n = 20). The following results were obtained:

Plasma	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	14.4	0.202	1.4
Sample 2	41.3	0.505	1.2
Sample 3	105	1.178	1.1

Plasma	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	14.6	0.267	1.8
Sample 2	42.5	1.057	2.5
Sample 3	106.3	2.163	2.0

Method comparison

A comparison of the Analyticon Fluitest® Lactate determination (y) with another commercial available assay (x) gave the following correlation (mg/dl):

$$y = 1.004 x + 4.120; \quad r^2 = 0.9923$$

Quality Control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio-Cal®	20 x 3 ml	#1420
Bio-Cal® E	10 x 3 ml	#1430

R4: Calibrator provided in kit

Calibration frequency:

Two point calibration is recommended:

- After lot change
- as required following quality control procedures

Calibration verification: Not necessary

Disposal:

Please note the legal regulations.

Literature:

1. Bablok W et al. A general regression procedure for method transformation. J Clin Chem Biochem. 1988;26:783-790.
2. Barhan D, Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system. Analyst 97.1972:142.
3. Bergmeyer HU (ed), New York, NY: Academic Press Inc; 1974;1464.
4. Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 4th ed. St.Louis. Elsevier Saunders; 2006.
5. Gutmann I, Wahlefeld A. Methods of Enzymatic Analysis. 2nd ed.
6. Noll F. Methods of Enzymatic Analysis. 2nd ed. Bergmeyer HU ed.
7. Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. J Clin Chem Clin Biochem. 1983;21:709-720.
8. Tietz NW. clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382 – 383.
9. Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. Ann. Clin. Biochem. 1969;6:24.
10. Westgard JO, Lahmeyer BL, Birnbaum ML. Clin Chem. 1972; 18:1334-1338.
11. Heil W, Rick W. Lactatbestimmung im Liquor cerebrospinalis zur Differenzierung zwischen bakterieller und viraler Meningitis. Internist 1984; 25:320-324.

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