BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B6251	200 / 600	R1	1 x	20 ml
		R2	1 x	8 ml
B6253	300 / 600*	R1	1 x	20 ml
		R2	1 x	8 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of Myoglobin in human serum and plasma.

Summary:

Myoglobin is a hemoprotein having a molecular weight of approximately 17800 daltons. It transports and binds oxygen reversibly in muscle cells. It is found predominantly in striated muscle tissue (skeletal and cardiac muscle). Myoglobin is liberated from damaged heart muscle cells such as occurs during acute myocardial infarction. An increase in myoglobin concentrations in blood can generally be detected 2 to 4 hours after the onset of pain, which is earlier than other cardiac markers such as CK, CK-MB or troponin. Depending on the therapeutic reperfusion measures taken, the myoglobin concentration reaches its maximum value after 4 to 12 hours and then decreases relatively rapidly to normal levels due to renal elimination (biological half-life: approx. 15minutes). A very rapid increase in the concentration of myoglobin occurs when therapeutic intervention is successful. The gradient of the concentration increase can be taken as an indication of the success of thrombolysis.

The myoglobin determination is of particular value in exclusion diagnosis for myocardial infarction: If there is no increase in the myoglobin concentration 6 hours after the onset of pain and after a repeat determination within a further 4 hours, then acute myocardial damage can essentially be excluded.

Increases in concentration of myoglobin not due to infarction may be a result of muscle trauma, crush syndrome, myopathy, muscle strain/stress, shock, rhabdomyolysis or decreased elimination due to renal failure.

Various nephelometric and turbidimetric methods are available for the determination of myoglobin. This myoglobin assay is based on the principle of immunological agglutination with latex particle reaction enhancement.

Test principle:

Immunoturbidimetric assay

• Sample and addition of R1 (buffer)

• Addition of R2 (anti-myoglobin antibody latex/buffer) and start of reaction

Latex-bound anti-myoglobin antibodies react with antigen in the sample to form an antigen/antibody complex which can be determined turbidimetrically after agglutination.

Reagent concentration:

R1:

Glycine buffer pH 9.0 170 mmol/l

NaCl 100 mmol/l

EDTA 50 mmol/l

Preservative

R2:

Latex particles coated with

Anti-human myoglobin antibodies (rabbit) 0,12 %

Glycine buffer pH 7.3 170 mmol/l

NaCl 100 mmol/l

Preservative

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Mix thoroughly before first use and once per week thereafter.

Unopened kit components: Up to the expiration date at 2-8°C

Onboard stability at 2-8°C: R1: 90 days

R2: 90 days

Specimen:

Collect serum using standard sampling tubes

Li-heparin-, Na citrate- or EDTA-Plasma

Stability: Analyze samples immediately or store frozen at -20°C for a maximum of 60 days.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an I index of 60 (approx. conjugated and unconjugated a bilirubin concentration of 60mg/dl).

Hemolysis: No significant interference up to a H index of 500 (approx. 500mg/dl hemoglobin).

Lipemia (Intralipid): No significant interference up to a L index of 750 (approx. 1500mg/dl triglyceride). There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factors up to 100 IU/ml do not interfere.

A high-dose hook-effect may occur at myoglobin concentrations > 25000 ng/ml.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring /reportable range:

3 – 400 ng/ml (µg/l)

Determine samples having higher activities via the rerun function using 0.9% NaCl.

Expected values:

Men: 23 - 72 ng/ml (µl/l)

Women: 19 - 51 ng/ml (µl/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes the myoglobin results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 3 ng/ml

The lower detection limit represents the lowest measurable myoglobin concentration that can be distinguished from zero.

Imprecision:

Reproducibility between day was determined using controls (n = 20). The following results were obtained:

Sample	Between day		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	36.7	1.05	2.86
Sample 2	70.9	1.31	1.85
Sample 3	274	5.79	2.11

Sample	Within run		
	Mean (ng/ml)	SD (ng/ml)	CV %
Sample 1	62.6	0.87	1.38
Sample 2	314	2.12	0.67
Sample 3	542	4.21	0.78

Method comparison:

A comparison of the Analyticon Turbitex® MYO (y) with a commercial obtainable assay (x) gave with 85 samples the following result:

$$y = 1.046x - 1.123; \quad r = 1.000$$

Quality control:

Human Control Serum

MYO Control 2 x 1 ml #1563

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The method was standardized against a nephelometric method.

Calibration Type: Logit-Log

S1: 0.9% NaCl

S2-5: Bio Cal® MYO Set 4 x 1 ml #14635

Calibration frequency:

Full recalibration is recommended

• after lot change

• as required following quality control procedures

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