Fluitest® MG XB

MAGNESIUM

Order information:

Catalog No.	Contents
3908	R1 2 x 100 ml R4 1 x 5 ml

Intended use:

In vitro test for the quantitative determination of magnesium in human serum, plasma and urine.

Summary:

Magnesium along with potassium is a major intracellular cation. Mg²⁺ is a cofactor of many enzyme systems. Thus, all ATP-dependent enzymatic reactions require Mg²⁺ as a cofactor in the ATP-magnesium complex. Approximately 69% of magnesium ions are stored in bone. The rest are part of the intermediary metabolism about 70% being present in free form while the other 30% is bound to proteins (especially albumin), citrates, phosphate, and other complex formers. The Mg²⁺ serum level is kept constant within very narrow limits (0.65-1.05 mmol/l). Regulation takes place mainly via the kidneys, especially via the ascending loop of Henle.

This assay is used for diagnosing and monitoring hypomagnesemia (magnesium deficiency) and hypermagnesemia (magnesium excess). Numerous studies have shown a correlation between magnesium deficiency and changes in calcium-, potassium- and phosphate homeostasis which are associated with cardiac disorders such as ventricular arrhythmias that cannot be treated by conventional therapy, increased sensitivity to digoxin coronary artery spasms, and sudden death. Additional concurrent symptoms include neuromuscular and neuropsychiatric disorders. Hypermagnesemia is found in acute and chronic renal failure, magnesium excess, and magnesium release from the intracellular space.

In addition to atomic absorption spectrometry (AAS), complexometric methods can also be used to determine magnesium.

The method described here is based on the reaction of magnesium with xylydyl blue in alkaline solution containing EGTA to mask the calcium in the sample

Test principles:

Colorimetric endpoint method

The color intensity of the magnesium complex is proportional to the magnesium concentration and can be used for the photometric determination.

Reagent concentration:

R1:	
Tris buffer pH 11.0	250 mmol/l
Xylydylblue	1 mmol/l
Detergent	15 g/l
R4:	
Magnesium	2.00 mg/dl (0.82 mmol/l)

Preparation and stability:

Reagent is ready for use

Unopened kit components: Up to the expiration date at 2-8°C

Avoid direct light. Keep bottles closed after use.

Specimen:

Collect serum using standard sampling tubes.

Li-heparin plasma

Stability:	7 days	at +20°C to +25°C
	7 days	at + 4°C to + 8°C
	1 year	at -20°C

Hemolyzed specimens are unacceptable.

Urine

Stability:	7 days	at +20°C to +25°C
	7 days	at + 4°C to + 8°C
	1 year	at -20°C

Collect urine sample in a metal-free container with no preservative. Acidify to approximately pH 3-4 prior to assay. Urine samples are manually diluted with distilled/deionized water (e.g. 1+4). Multiply the result by the appropriate factor (e.g. 5).

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitation - interference:

Criterion: Recovery within ± 10% of initial values

Serum/plasma

Icterus: No significant interference up to an index I of 10 corresponding to an approximate bilirubin concentration of 10 mg/dl.

Hemolysis interferes due to magnesium released from erythrocytes. Lipemia (Intralipid): No significant interference up to an index L of 475 (approximate triglycerides concentration: 950 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request

Materials provided

• Working solutions as described above.

Additional materials required

• Calibrators and controls as indicated below

Manual procedure:

Wavelength:	546 nm (520nm)
Temperature:	+25°C
Cuvette:	1 cm
Zero adjustment:	only one reagent blank per series

	Blank	Calibr./Standard	Sample
Reagent/R1	1000 µl	1000 µl	1000 µl
Aqua dest.	10 µl	---	---
Standard/R4	---	10 µl	---
Sample	---	---	10 µl

Mix well and incubate for 5 minutes. Read absorbance of sample and standard against reagent blank within 60min.

Calculation:

$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Calibr./Stand.}}$ x concentration of calibr./Standard = Magnesium conc.

Measuring:

Serum/plasma

Measuring range: 0.04 mmol/l - 2.055 mmol/l (5 mg/dl)

In case of higher results, dilute sample 1:4 with saline solution and repeat test.

Multiply result by 4

Reference values:

Serum	1.9 - 2.5 mg/dl	(0.78 - 1.0 mmol/l)
Cerebrospinal fluid	2.4 - 3.1 mg/dl	(0.98 - 1.3 mmol/l)
Urin	1 - 10 mg/dl	(0.41 - 4.1 mmol/l)
24 h Urin	50 - 200 mg/dl	(20 - 80 mmol/l)

Reference values according to Tietz

Serum/plasma:	1.3-2.1 mEq/l	(1.59 - 2.56 mg/dl)
2-4 days:	1.2-1.8 mEq/l	(1.46 - 2.20 mg/dl)
5 months-6 years:	1.4-1.9 mEq/l	(1.71 - 2.29 mg/dl)
6-12 years:	1.4-1.7 mEq/l	(1.71 - 2.07 mg/dl)
12-20 years:	1.3-1.8 mEq/l	(1.59 - 2.20 mg/dl)
Urine:	1.0-24.0 mEq/24 hours	(12.2 - 292 mg/24 hours)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the magnesium results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.04 mmol/l (0.09 mg/dl)

The lower detection limit represents the lowest measurable magnesium concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Within run			
Sample	Mean mmol/l	SD mmol/l	CV %
sample 1	1.06	0.012	1.13
sample 2	1.28	0.013	1.02
sample 3	1.79	0.010	0.56

Between day			
Sample	Mean mmol/l	SD mmol/l	CV %
sample 1	1.03	0.023	2.28
sample 2	1.27	0.023	1.28
sample 3	1.81	0.062	3.43

Method comparison:

A comparison of the Analyticon Fluitest MG (y) with a commercial obtainable assay (x) gave with samples the following result:

$$y = 0.903x + 0.121; r = 0.986$$

Quality control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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