# NH<sub>3</sub>

## AMMONIA

### Order information:

Catalog No.	Contents
H1701	Hit 1 / 917 (ILab*/ AU*) R1 2 x 20 ml R2a 2 x for 5 ml R2b 2 x 5 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911/917: ACN 478  
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

For the quantitative determination of ammonia (NH<sub>3</sub>) in plasma on automated clinical chemistry analyzers.

### Summary:

In 1859, Berthelot described a reaction between ammonia and an alkaline solution of phenol hypochlorite suitable for the determination of ammonia. The assay, however, proved to be subject to interferences, and several alternatives have been proposed to eliminate the problems inherent in the method.

In 1963, Kirsten, et al introduced an enzymatic method for ammonia determination based on the action of glutamate dehydrogenase. Although the enzymatic method proved to be highly specific and utilized direct evaluation based on the molar absorptivity of NADH, several problems, including difficulties in stabilizing the end reaction, were encountered.

The method presented here is based on Da Fonseca-Wollheim's modification of the Kirsten reaction. The original enzymatic method is improved by the addition of ADP to the reaction mixture, the use of NADPH in place of NADH to eliminate interference from the reaction of endogenous LDH with endogenous pyruvate, and the substitution of plasma for deproteinized supernatant.

Ammonia measurements are used in the diagnosis and treatment of severe liver disorders such as cirrhosis, hepatitis, and Reye's Syndrome.

### Test principle:

Enzymatic kinetic colorimetric assay  
In the reaction catalyzed by glutamate dehydrogenase (GLDH), ammonia reacts with  $\alpha$ -ketoglutarate and NADPH to form glutamate and NADP.

- Sample and addition of R1
- Addition of R2 and start of reaction:

The inclusion of ADP in the reaction mixture causes an acceleration of the rate of conversion and stabilizes the GLDH in the indicated pH range. The amount of NADPH oxidized during the reaction is equivalent to the amount of ammonia in the specimen and can be measured photo-metrically by the resulting decrease in absorbance.

### Reagent concentration:

<b>R1:</b>	
Triethanolamine buffer pH 8.6	151 mmol/l
$\alpha$ -Ketoglutarate	16.6 mmol/l
ADP	$\geq 1.2$ mmol/l
Preservatives	

<b>R2a:</b>	
NADPH	$\geq 458$ $\mu\text{mol/l}$
GLDH (EC 1.4.1.3; bovine liver; +25°C)	$\geq 24.3$ U/ml
Triethanolamine buffer pH 8.6	151 mmol/l
$\alpha$ -Ketoglutarate	16.6 mmol/l
ADP	$\geq 1.2$ mmol/l
Preservatives	

<b>R2b:</b>	
Triethanolamine buffer pH 8.6	151 mmol/l
$\alpha$ -Ketoglutarate	16.6 mmol/l
ADP	$\geq 1.2$ mmol/l
Preservatives	

### Preparation and stability:

R1: Ready for use.  
R2: Reconstitute one Bottle R2a (GLDH) with the content of one bottle R2b and let stand 20 minutes at room temperature, occasionally swirling gently.

Unopened kit components: at +2° to +8°C up to the expiration date.  
Reconstituted GLDH: at +2° to +8°C 42 days when stored tightly closed.  
Onboard: R1: 28 days opened and refrigerated on the analyzer when protected from contamination  
R2: 14 days opened and refrigerated on the analyzer when protected from contamination

### Specimen:

EDTA plasma, nonhemolyzed. Draw the specimen from a stasis-free vein and centrifuge in a tightly closed tube as soon as possible.  
Stability: Place specimen on ice and assay immediately.  
Separated plasma: 3 hours at +4°C in a tightly closed container.

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

Warning! R2 contains hazardous components.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

NOTE: Do not leave sample, calibrator, or quality control cups open for longer than 15 minutes. Results may be affected. Repeat assays must be performed on freshly poured cups, due to evaporation of ammonia.

### Limitations - interference:

Criterion: Recovery within  $\pm 10\%$  of initial value.

Icterus: No significant interference up to an index I of 60

(approximate conjugated and unconjugated bilirubin: 60 mg/dl)

Hemolysis: No significant interference up to an index H of 50 (approximate haemoglobin concentration: 50 mg/d).

Lipemia (Intralipid): No significant interference up to an index L of 250 (approximate triglycerides concentration: 500 mg/dl) There is poor correlation between turbidity and triglycerides concentration. Absorbance flagging may also be observed as a result.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician of an absorbance increase.

### Testing procedure:

#### Materials provided

- Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

### Calculation:

Roche/Hitachi systems automatically calculate the ammonia concentration of each sample.

Conversion factor:  $\mu\text{g/dl} \times 0.587 = \mu\text{mol/l}$

### Measuring /reportable range:

Plasma: 10-1000  $\mu\text{g/dl}$  (5.87-587  $\mu\text{mol/l}$ )

#### Specimen dilution

Manually dilute specimens above the reportable range with freshly distilled or deionized ammonia-free water (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2).

NOTE: Do not report results generated from automatic rerun unless a fresh sample is poured.

### Expected values:

Male: 25-94  $\mu\text{g/dl}$  (14.7 - 55.3  $\mu\text{mol/l}$ )

Female: 19-82  $\mu\text{g/dl}$  (11.2 - 48.2  $\mu\text{mol/l}$ )

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, ammonia results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit):

Detection limit: 10  $\mu\text{g/dl}$  (5.87  $\mu\text{mol/l}$ )

The lower detection limit represents the lowest measurable substance concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

### Imprecision:

Within-Run precision was determined using human-based controls in an internal protocol run (n = 20). Lot-to-lot precision was determined with three different lots of human-based controls. The following results were obtained:

Sample	Within run		
	Mean $\mu\text{g/dl}$	SD $\mu\text{g/dl}$	CV %
Calibrator A/E/C*	519	3.43	0.66
Control A/E/C* Normal	115	3.89	3.37
Control A/E/C* Abnormal	391	3.78	0.97
Sample	Lot to lot		
	Mean Recovery %	SD %	CV %
Control A/E/C* Normal	99.3	5.75	5.79
Control A/E/C* Abnormal	94.3	3.14	3.33

\*Ammonia/Ethanol/CO<sub>2</sub>

**Method comparison:**

A comparison of the ammonia determination on a Hitachi 911 analyzer using the Analyticon ammonia assay (y) with another ammonia assay (x) gave the following correlation ( $\mu\text{g/dl}$ ):

$$y = -0.598 + 0.995 x; \quad r = 0.999$$

**Quality Control:**

Suitable control material should be used.

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

**Calibration:**

Standardization: The method was standardized against a gravimetric standard.

S1: 0.9% NaCl

S2: suitable calibrator for Ammonia

e.g. the Ammonia/Ethanol/CO<sub>2</sub> Calibrator from Roche

**Calibration frequency:**

Two point calibration is recommended:

- every 24 hours
- after reagent bottle change
- after reagent lot change
- as required following quality control procedures

**Disposal:**

Please note the legal regulations.

**Literature:**

1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988; 26:783-790.
2. Berthelot MPE. Repert Chem Appl. 1859:282.
3. Da Fonseca-Wollheim F. Z Klin Chem Klin Biochem. 1973;1 1 :421 ,426. Dokumentation Roche Diagnostics.
4. Hohorst HJ. Biochem Z. 1956;328:507.
5. Kaplan LA, Pesce AJ. Clinical Chemistry Theory, Analysis and Correlation. St. Louis, Mo: CV Mosby Co; 1984:1231.
6. Kirsten E, et al. Biochem Z. 1963;337:312.
7. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983; 21:709-720.
8. Prellwitz W, et al. Med Welt. 1976;27:1277.
9. Siegel JM, Montgomery GA. Arch Biochem/Biophys. 1956;82:288.
10. Tietz NW. Clinical Guide to Laboratory Tests. 3<sup>rd</sup> ed. Philadelphia, PA: W.B. Saunders Co; 1995:44.

Text passages with grey background were changed in the latest revision of this package insert.

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