P-5-P

PYRIDOXAL-5-PHOSPHATE

Order information:

Catalog No.	Contents
1101	R0 6 x 3 ml

Intended use:

Additional reagent for the quantitative determination of AST (GOT) and ALT (GPT) in human serum and plasma with P-5-P activation.

Summary:

Vitamin B₆ occurs in three natural forms: pyridoxine, pyridoxal and pyridoxamine. All three of these moieties are converted to pyridoxal-5-phosphate (P-5-P or PLP) in the body by an enzymatic reaction with ATP. P-5-P serves as a highly versatile coenzyme for many types of reaction: decarboxylation, racemization, transamination and other metabolic functions. P-5-P undergoes amination to pyridoxamine phosphate. The amino group is then transferred to an α -keto acid, returning pyridoxamine phosphate to its pyridoxal-5-phosphate form. Pyridoxal-5-phosphate is essential to these aminotransferase reactions. In normal individuals, part of the aminotransferase enzymes are not saturated with P-5-P. In some hepatic disorders, in vitamin B₆ deficiency, and following renal dialysis, the aminotransferase enzymes may be significantly stripped of P-5-P, thus rendering the enzymes inactive. By adding P-5-P, the transaminases are stabilised and false decreased results of samples are avoided. The addition of pyridoxal-5-phosphate to reaction media for catalytic determinations is widely accepted (IFCC/DGKC/SFBC).

Test principle:

Refer to appropriate insert or application sheet for test-specific information.

Reagent concentration:

R0:	
Good's Buffer, pH 9,6	0,1 mol/L
Pyridoxal-5-Phosphate (P-5-P)	13,8 mmol/L

Preparation and stability:

R0: Ready for use

The reagent is stable up to the imprinted expiry date when stored at +2°C to +8°C. Protect from light. Do not freeze!

Reagent handling:

The reagent is ready to mix with Fluitest GOT/AST or GPT/ALT reagent R1. Take 100 parts R1 and 1 part P-5-P-solution. Mix well.

(e.g. 10ml R1 + 100 μ l P-5-P).

Close bottle tight right after use!

Refer to appropriate package insert for further details.

P-5-P activated working solution stability:

6 days	at +2 - +8°C
24 hours	at +15 - +25°C

Assay:

Materials provided

• Working solutions as described above

Additional materials required

• Fluitest GOT (AST) / GPT (ALT)

Testing procedure:

Refer to Fluitest GOT (AST) / GPT (ALT) package insert.

Reference values for tests with P-5-P activation:

According to DGKC:

+37°C	U/L	μ kat/L
GOT AST/P5P Männer	10 – 50	0,17 – 0,83
GOT AST/P5P Frauen	10 – 35	0,17 – 0,58
GPT ALT/P5P Männer	10 – 50	0,17 – 0,83
GPT ALT/P5P Frauen	10 – 35	0,17 – 0,58

According to Tietz:

+37°C	U/L	μ kat/L
GOT AST/P5P Männer	15 – 40	0,26 – 0,68
GOT AST/P5P Frauen	13 – 35	0,22 – 0,60
GPT ALT/P5P Männer	10 – 40	0,17 – 0,68
GPT ALT/P5P Frauen	7 – 35	0,12 – 0,60

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

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