Fluitest® P-AMYL IFCC

PANCREATIC α -AMYLASE

Order information:

Catalog No.	Contents
H2101 Hit I (ILab*)	R1 6 x 47 ml R2 6 x 11 ml
H2103 Hit 917 (AU*)	R1 6 x 60 ml R2 6 x 14 ml

System information:

Hitachi 911: ACN 029
Hitachi 917: ACN 571
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Enzymatic in vitro test for the quantitative determination of pancreatic α -amylase in human serum, plasma and urine.

Summary:

The α -amylases (1,4- α -D-glucanohydrolases, EC 3.2.1.1) catalyze the hydrolytic degradation of polymeric carbohydrates such as amylose, amylopectin and glycogen by cleaving 1,4- α -glucosidic bonds. In polysaccharides and oligosaccharides, several glycosidic bonds are hydrolyzed simultaneously. Maltotriose, the smallest subunit, is converted into maltose and glucose, but very slowly.

Two types of α -amylases can be distinguished, the pancreatic type (P-type) and the salivary type (S-type). Whereas the P-type can be attributed almost exclusively to the pancreas and is therefore organ-specific, the S-type can originate from a number of sites. As well as appearing in the salivary glands it can also be found in tears, sweat, human milk, amniotic fluid, the lungs, testes and the epithelium of the fallopian tube.

Because of the sparsity of specific clinical symptoms of pancreatic diseases, α -amylase determinations are of considerable importance in pancreatic diagnostics. The determination of pancreas-specific α -amylase instead of total α -amylase is of advantage here.


The determination of pancreatic α -amylase is suitable for the diagnosis and monitoring of acute pancreatitis and acute attacks during chronic pancreatitis. In terms of clinical sensitivity and specificity, the diagnostic value of pancreatic- α -amylase is comparable to that of lipase, the generally recognized pancreas-specific enzyme. The sensitivity of pancreatic α -amylase is 38% higher than that of total α -amylase in the diagnosis of acute pancreatitis when - as commonly used - three times the upper normal limit is taken as the criterion. A variety of methods have been described for determining pancreatic α -amylase: radio- and enzyme-immunoassays as well as the partial inhibition of salivary α -amylase by an inhibitor derived from wheat germ and calculation of the pancreatic α -amylase from the remaining and total amylase activities.

The kinetic method described here is based on inhibition of the activity of human salivary α -amylase by two different monoclonal antibodies and the well-proven cleavage of 4,6-ethylidene-(G₇)-1,4-nitrophenyl-(G₁)- α -D-maltoheptaoside (Ethylidene Protected Substrate = EPS) by pancreatic α -amylase followed by hydrolysis of all the degradation products to p-nitrophenol with the aid of α -glucosidase (100% chromophore liberation). The results of this method correlate with those obtained by HPLC.

Test principle: (simplified)

Enzymatic colorimetric assay, according to the IFCC method for total amylase.

In the first incubation step, the activity of human salivary (α -amylase is inhibited by two 2 different monoclonal antibodies (without affecting the pancreatic- α -amylase). In the second reaction step defined oligosaccharides such as 4,6-ethylidene-(G₇) p-nitrophenyl-(G₁)- α -D-maltoheptaoside (ethylidene-G₇PNP) are cleaved under the catalytic action of pancreatic α -amylases. The G₂PNP G₃PNP and G₄PNP fragments so formed are completely hydrolyzed to p-nitrophenol and glucose by α -glucosidase.

(PNP = p-nitrophenol; G = glucose)

The color intensity of the p-nitrophenol formed is directly proportional to the α -amylase activity and is measured photometrically.

Reagent Concentration:

R1:	
Monoclonal antibodies to s-amylase	40 mg/l
Hepes-buffer, pH 7.15	52.5 mmol/l
NaCl	87 mmol/l
MgCl ₂	12.6 mmol/l
CaCl ₂	0.075 mmol/l
α -Glucosidase, mod.	≥ 4 KU/l
R2:	
Hepes- buffer, pH 7.15	52.5 mmol/l
4,6-ethylidene-G ₇ PNP	22 mmol/l

Preparation and stability:

R1: Ready for use. R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability: R1: 28 days
R2: 28 days

Specimen:

Serum/plasma

Collect serum using standard sampling tubes
Heparinized or EDTA-plasma.
Stability: 7 days at +20°C+25°C
1 month at +2°C+8°C

Urine

Collect without additives.
Stability: 2 days at +20°C+25°C
10 days at +2°C+8°C

Pancreatic- α -Amylase is unstable in acid urine. Assay promptly or adjust pH to alkaline range (about pH 7) before storage.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

The residual activity of salivary α -amylase is approx. 3%. In rare cases, very high activities of salivary α -amylase can hence lead to elevated values being measured for pancreatic α -amylase.

A slight change in the yellow coloration of solution 2 does not interfere with the performance of the test.

Do not pipette by mouth, and ensure that the reagent does not come into contact with the skin. (Saliva and sweat contain α -amylase!)

Criterion: Recovery within $\pm 10\%$ of initial value.

Icterus: No significant interference up to an index I of 65 (approximate conjugated and unconjugated bilirubin: 65 mg/dl)

Hemolysis: No significant interference up to an index H of 150 (approximate haemoglobin concentration: 150 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 850 (approximate triglycerides concentration: 1700 mg/dl). There is poor correlation between turbidity and triglycerides concentration

Highly turbid and grossly lipemic samples may cause absorption flags.

Patients with macroamylase may have elevated p-amylase results. The elevation is not due to an insufficient inhibition of salivary amylase in the serum immune complex. It is caused by a higher than normal level of p-amylase since the immune complex is not subject to glomerular filtration. This elevated p-amylase is not diagnostic for pancreatitis. However, measurement of an elevated p-amylase in urine is confirmatory of pancreatitis, pancreatic trauma, or pancreatic carcinoma as the amylase released is not completely bound by the immune complex and thus subject to glomerular filtration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl - solution

Manual Procedure:		
Wavelength:	405 nm (400-420nm)	
Temperature:	+37°C	
Cuvette:	1 cm light path	
Zero adjustment:	against air or distilled water	
	Serum / plasma	Urine
R1	1000 μ l	1000 μ l
Sample	40 μ l	20 μ l
Mix and incubate 3 minutes. Then add:		
R2	200 μ l	200 μ l
Mix and incubate 2 minutes, read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes. Determine the mean absorbance change per minute $\Delta A/\text{min}$ and use this for the calculation.		
Calculation:		
<i>Serum / plasma:</i>		
$\Delta A/\text{min} \times 3.130 = \text{Activity (U/l)}$		
<i>Urine:</i>		
$\Delta A/\text{min} \times 6.160 = \text{Activity (U/l)}$		

Measuring / reportable range:

Measuring range: 3-1500 U/l (0.05-25.00 µkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl - solution or distilled/deionized water (e.g. 1 + 4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

Expected values:

	U/l	µkat/l
Serum/plasma	13 – 53	0.22 – 0.88
Spontaneously voided urine	≤ 350	≤ 5.83
Pancreatic-α-amylase/ Creatinin Quotient	≤ 205 U/g	≤ 3.42 µkat/g

Pancreatic α-Amylase/ creatinine quotient

To allow for fluctuations in the pancreatic α-amylase activity in urine, it is advisable to determine the α-amylase / creatinine quotient. To do this, determine the pancreatic α-amylase activity and creatinine concentration in spontaneously voided urine.

$$\text{Quotient [U/g] or } \mu\text{kat/mmol} = \frac{\text{pancreatic } \alpha\text{-amylase (U/l or } \mu\text{kat/l)}}{\text{creatinine [g/l or mmol/l]}}$$

Amylase / Creatinine Clearance Ratio (ACCR)

The ACCR is calculated from amylase activity and creatinine concentration. Both the serum and urine samples should be collected at the same time.

$$\text{ACCR [\%]} = \frac{\text{Urine amylase [U/l]} \times \text{serum creatinine [mg/l]}}{\text{Serum amylase [U/l]} \times \text{urine creatinine [mg/l]}} \times 100$$

ACCR approximately equal to 2-5%

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the pancreatic α-amylase results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 3 U/l or 0.05 µkat/l

The lower detection limit represents the lowest measurable α-amylase concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 20). The following results were obtained:

Sample	Between day			Within run		
	Mean U/l	SD U/l	CV %	Mean U/l	SD U/l	CV %
Sample 1	40.2	0.93	2.30	37.3	0.83	2.23
Sample 2	106	2.19	2.06	100	1.03	1.02
Sample 3	162	3.00	1.85	478	3.06	0.64

Method comparison:

A comparison of the Analyticon Fluitest P-AMYL (y) with a commercial obtainable assay (x) gave the following result: (U/l)

$$y = 1.032 x + 1.1911 x; r = 1.000$$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S 1: 0.9% NaCl

S 2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency

Two point calibration is recommended:

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Literature:

- Gerber M, Wulff K. Fortschritte in der spezifischen Bestimmung der Pankreas-α-amylase. Laboratoriumsmedizin 1988;12:110-113.
- Greiling H, Gressner AM ed. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd edition, Stuttgart/New York: Schattauer Verlag, 1995.
- Junge W, Troge B, Klein G et al. Evaluation of a New Assay for Pancreatic Amylase: Performance Characteristics and Estimation of Reference Intervals. Clin Biochem 1989; 22:109-114.
- Keller H ed. Klinisch-chemische Labordiagnostik für die Praxis, 2nd Stuttgart/New York: Georg Thieme Verlag, 1991;354-361.
- Kruse-Jarres JD, Hafkenscheid JCM, Hohenwallner W et al. Evaluation of a New α-Amylase Assay Using 4.6-Ethylidene-(G₇)-1-4-nitrophenyl-(G₁)-α-D-maltoheptaoside as Substrate. J Clin Chem Clin Biochem 1989;27:103-113.
- Kurrie-Weitenhiller A, Hölzel W, Engel D et al. Method for the determination of total and pancreatic α-amylase based on 100% cleavage of the protected substrate ethylidene-4-nitrophenyl-maltoheptaoside. Clin Chem 1996;42:98.
- Lott J. Inflammatory diseases of the pancreas CRC Critical Reviews in Clinical Laboratory Sciences 1982;17:201--228.
- Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983; 21:709-720.
- Rizzotti P, Klein G. Evaluation of a Specific immunoinhibition Method for the Determination of Pancreatic α-Amylase. Eur J Clin Chem Clin Biochem 1994;32:97-106. Salt WB II, Schenker S. Amylase-its clinical significance: a review of the literature [Review]. Medicine 1976; 55:269-281.
- Tietz NW, Huang WY, Rauh DF et al. Laboratory tests in the differential diagnosis of hyperamylasemia. Clin Chem 1986; 32: 301-307.
- Tietz NW, Burlina A, Gerhardt W et al. Multicenter Evaluation of a Specific Pancreatic Isoamylase Assay Based on a Double Monoclonal Antibody Technique. Clin Chem 1988;34:2096-2102.
- Waldenström J, Bouman A et al. Evaluation of the Assays for Total and Pancreatic α-Amylase based on 100% Cleavage of Et-G₇-PNP at 6 European Clinical Centres (Poster Medlab 97). Basel, Switzerland: 12th IFCC European Congress of Clinical Chemistry, August 17-22 1997

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Literatur:

1. Gerber M, Wulff K. Fortschritte in der spezifischen Bestimmung der Pankreas- α -amylase. Laboratoriumsmedizin 1988;12:110-113
2. Greiling H, Gressner AM (Hrsg.). Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer Verlag, 1995.
3. Junge W, Troge B, Klein G et al. Evaluation of a New Assay for Pancreatic Amylase: Performance Characteristics and Estimation of Reference Intervals. Clin Biochem 1989; 22:109-114.
4. Keller H (Hrsg.). Klinisch- chemische Labordiagnostik für die Praxis, 2. Auflage. Stuttgart/New York: Georg Thieme Verlag, 1991;354-361.
5. Kruse-Jarres JD, Hafkenscheid JCM, Hohenwallner W et al. Evaluation of a New α -Amylase Assay Using 4,6-Ethylidene-(G₇)- 1-4-nitrophenyl-(G₁)- α -D-maltoheptaoside as Substrate. J Clin Chem Clin Biochem 1989;27:103-113.
6. Kurrel-Weitenhiller A, Hölzel W, Engel D et al. Method for the determination of total and pancreatic α -amylase based on 100% cleavage of the protected substrate ethylidene-4-nitrophenyl-maltoheptaoside. Clin Chem 1996;42:98.
7. Lott J. Inflammatory diseases of the pancreas CRC Critical Reviews in Clinical laboratory Sciences 1982;17:201--228.
8. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurement from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720
9. Rizzotti P, Klein G. Evaluation of a Specific immunoinhibition Method for the Determination of Pancreatic α - Amylase. Eur J Clin Chem Clin Biochem 1994;32:97-106
10. Tietz NW, Huang WY, Rauh DF et al. Laboratory tests in the differential diagnosis of hyperamylasemia. Clin Chem 1986; 32: 301-307. Junge W.,
11. Tietz NW, Burlina A, Gerhardt W et al. Multicenter Evaluation of a Specific Pancreatic Isoamylase Assay Based on a Double Monoclonal Antibody Technique. Clin Chem 1988;34:2096-2102.
12. Waldenström J, Bouman A et al. Evaluation of the Assay for Total and Pancreatic α -Amylase based on 100% Cleavage of Et-G7-PNP at 6 European Clinical Centres (Poster Medlab 97). Basel, Switzerland : IFCC European Congress of Clinical Chemistry, August 17-22, 1997