### Order information:

Catalog No.		Contents					
H8501	Hit I (ILab*)	R1	12 x	50 ml	R2	6 x	43 ml
H8503	Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	30 ml
AU8504	AU	R1	5 x	20 ml	R2	5 x	12 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems

### System information:

Hitachi 911: ACN 161  
 Hitachi 917: ACN 714  
 Hitachi 917 (Urine): ACN 716

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

In vitro test for the quantitative determination of phosphorus in human serum, plasma and urine

### Summary:

88% of the phosphorus contained in the body is localized in bone in the form of calcium phosphate as the apatite  $Ca^{2+}[Ca_3(PO_4)_2]^{2-}$ . The remainder is involved in intermediary carbohydrate metabolism and in physiologically important substances such as phospholipids, nucleic acids and ATP. Phosphorus occurs in blood in the form of inorganic phosphate and in organically bound phosphoric acid. The small amount of extracellular organic phosphorus is found almost exclusively in the form of phospholipids.

The ratio of phosphate to calcium in the blood is approximately 6:10. An increase in the level of phosphorus causes a decrease in the calcium level. The mechanism is influenced by interactions between parathormone and vitamin D. Hypoparathyroidism, vitamin D intoxication and renal failure with decreased glomerular phosphate filtration give rise to hyperphosphatemia. Hypophosphatemia occurs in rickets, hyperparathyroidism and Fanconi's syndrome. The preferred method for the determination of inorganic phosphorus is based on the formation of ammonium phosphomolybdate with subsequent reduction to molybdenum blue. Reagent stability problems often occur with this method. The method presented here is based on the reaction of phosphate with ammonium molybdate to form ammonium phosphomolybdate without reduction. The addition of an accelerator gives rise to a more rapid rate of reaction and the application of sample blanking yields more precise results.

### Test principle:

Inorganic phosphate forms an ammonium phosphomolybdate complex having the formula  $(NH_4)_2[PO_4(MoO_3)_{12}]$  with ammonium molybdate in the presence of sulfuric acid. The complex is determined photometrically in the ultraviolet region (340 nm).

### Reagent Concentration:

<b>R1:</b>	
H <sub>2</sub> SO <sub>4</sub>	280 mmol/l
NaCl	154 mmol/l
Detergent	2%
<b>R2:</b>	
H <sub>2</sub> SO <sub>4</sub>	280 mmol/l
Ammoniummolybdate	3 mmol/l
NaCl	154 mmol/l

### Preparation and stability:

R1: Ready for use.  
 R2: Ready for use.  
 Unopened Kits are stable up to the expiry date when stored at +2 to +8°C.  
 Onboard stability: R1: 63 days  
 R2: 63 days

### Specimen:

Collect serum using standard sampling tubes

Heparinized or EDTA plasma

Stability: 1 day at +20°C to +25°C  
 4 days at + 2°C to + 8°C

Urine: Stable for 8 hours at +15°C to +25°C

Centrifuge samples containing precipitate before performing the assay.

### Limitations - interference:

Criterion: Recovery within ± 10% of initial value.

Icterus: No significant interference up to an index I 92 (approximate bilirubin concentration: 92 mg/dl).

Hemolysis: Significant positive interference at an index H 1100 (approximate hemoglobin concentration: 1100 mg/dl). Note: This interference results from inorganic phosphates produced by the action of phosphatases on organic phosphates, both of which are released from the red cells upon hemolysis.

Lipemia (Intralipid): No significant interference up to an index L of 825 (approximate triglycerides concentration: 1650 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below  
 • 0.9% NaCl

### Manual Testing Procedure:

Wavelength:	340 nm (334nm)	
Reaction temperature:	+25°C/ +30°C/ +37°C	
Cuvette:	1 cm light path	
Zero adjustment:	One reagent blank for each series	

	Blank	Sample/ Calibrator
R1	750 µl	750 µl
R2	330 µl	330 µl
Sample/ Calibrator	---	15 µl

Mix and incubate 5 minutes. Read the absorbance of specimen against reagent blank within 60min. (ΔA)

### Calculation:

ΔA Sample x Calibrator conc. = Phosphate conc.

SI Units: (mg/dl) x 0.323 = mmol/l

### Measuring/reportable range:

Serum

0.3 - 20 mg/dl (0.10 - 6.46 mmol/l)

Determine samples containing higher phosphate concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled or deionized water (e.g. 1:4). Multiply the result by the appropriate dilution factor (e.g. 4).

### Expected values:

2.7-4.5 mg/dl (0.87-1.45 mmol/l)

Expected values for children are given in „Pediatric reference ranges“, 3<sup>rd</sup> edition, S.J. Soldin, C. Brugnara, I.M. Hicks, AACCC Press.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the phosphorus results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit)

0.3 mg/dl (0.1 mmol/l)

The lower detection limit represents the lowest measurable phosphorus concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Sample	Within run		
	Mean (mmol/l)	SD (mmol/l)	CV %
Controlserum 1	1.25	0.016	1.36
Controlserum 2	1.57	0.012	0.76
Controlserum 3	2.11	0.014	0.66

Sample	Between Day		
	Mean (mmol/l)	SD (mmol/l)	CV %
Controlserum 1	1.31	0.025	1.91
Controlserum 2	1.68	0.034	2.02
Controlserum 3	1.90	0.022	1.16

### Method comparison:

A comparison of the Analyticon Fluitest PHOS (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.030 x + 0.042; r = 0.991$$

# Fluitest® PHOS

INORGANIC PHOSPHORUS

## Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## Calibration for Hitachi systems:

S1: 0.9% NaCl

S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420

## Calibration frequency:

Two-point calibration is recommended:

- after reagent lot change
- as required following quality control procedures

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