Order information:

Catalog No.	Contents	
1914	R1 1 x 100 ml	R2 1 x 100 ml
	R4 1 x 5 ml	

Intended use:

In vitro test for the quantitative determination of phosphorus in human serum, plasma and urine

Summary:

88% of the phosphorus contained in the body is localized in bone in the form of calcium phosphate as the apatite $\text{Ca}^{2+}[\text{Ca}_3(\text{PO})]^{2-}_3$. The remainder is involved in intermediary carbohydrate metabolism and in physiologically important substances such as phospholipids, nucleic acids and ATP. Phosphorus occurs in blood in the form of inorganic phosphate and in organically bound phosphoric acid. The small amount of extracellular organic phosphorus is found almost exclusively in the form of phospholipids.

The ratio of phosphate to calcium in the blood is approximately 6:10. An increase in the level of phosphorus causes a decrease in the calcium level. The mechanism is influenced by interactions between parathormone and vitamin D. Hypoparathyroidism, vitamin D intoxication and renal failure with decreased glomerular phosphate filtration give rise to hyperphosphatemia. Hypophosphatemia occurs in rickets, hyperparathyroidism and Fanconi's syndrome.

The preferred method for the determination of inorganic phosphorus is based on the formation of ammonium phosphomolybdate with subsequent reduction to molybdenum blue. Reagent stability problems often occur with this method. The method presented here is based on the reaction of phosphate with ammonium molybdate to form ammonium phosphomolybdate without reduction. The addition of an accelerator gives rise to a more rapid rate of reaction and the application of sample blanking yields more precise results.

Test principle:

Inorganic phosphate forms an ammonium phosphomolybdate complex with the formula $(\text{NH}_4)_3[\text{PO}_4(\text{MoO}_3)_{12}]$ with ammonium molybdate in the presence of sulfuric acid. The complex is determined photometrically in the ultraviolet region (340 nm).

Reagent Concentration:

R1:	
H ₂ SO ₄	280 mmol/l
NaCl	154 mmol/l
Detergent	2%
R2:	
H ₂ SO ₄	280 mmol/l
Ammoniummolybdate	3 mmol/l
NaCl	154 mmol/l
R4:	
Phosphorus	5 mg/dl (1.62 mmol/l)

Preparation and stability:

For assay series:

A reaction mixture of R1 + R2 in equal volumes may be used. The reaction mixture of R1 + R2 is stable for

8 weeks	at +2 to +8°C and
24 hours	at room temperature.

For single tests:

All solutions are ready to use. Stable up to the expiry date when stored at +2 to +8°C.

Specimen:

Collect serum using standard sampling tubes
Heparinized or EDTA plasma, urine

Stability:	1 day	at +15°C to +25°C
	4 days	at +2°C to +8°C
	1 year	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ± 10% of initial value.

Icterus: No significant interference up to an index I of 92 (approximate bilirubin concentration: 92 mg/dl).

Hemolysis: Significant positive interference at an index H of 1100 (approximate hemoglobin concentration: 1100 mg/dl). Note: This interference results from inorganic phosphates produced by the action of phosphatases on organic phosphates, both of which are released from the red cells upon hemolysis.

Lipemia (Intralipid): No significant interference up to an index L of 825 (approximate triglycerides concentration: 1650 mg/dl). There is poor correlation between turbidity and triglycerides concentration. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Manual procedure:

Wavelength:	340 nm, Hg 334nm
Temperature:	+25°C / +30°C / +37°C
Cuvette:	1 cm
Zero adjustment:	against reagent blank

For assays series:

	Blank	Calib./Stand.	Sample
Reaction mixt.	1000 µl	1000 µl	1000 µl
Sample	---	---	10 µl
Standard/R4	---	10 µl	---

For single tests:

	Blank	Calib./Stand.	Sample
Reagent R1	500 µl	500 µl	500 µl
Reagent R2	500 µl	500 µl	500 µl
Aqua dest.	10 µl	---	---
Standard/R4	---	10 µl	---
Sample	---	---	10 µl

Mix. Incubate for 5 minutes. Read absorbance of specimen against reagent blank within 60 minutes (ΔA).

Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times \text{Conc. Calib./Stand.} = \text{Phosphorus Conc}$$

Measuring/reportable range:

Serum

0.30 - 20 mg/dl (0.10 - 6.46 mmol/l)

Determine samples containing higher phosphate concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled or deionized water (e.g. 1+3). Multiply the result by the appropriate dilution factor (e.g. 4).

Expected values:

2.7-4.5 mg/dl (0.87-1.45 mmol/l)

Expected values for children are given in „Pediatric reference ranges“ 3rd edition, S.J. Soldin, C. Brugnara, I.M. Hicks, AACC Press.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the phosphorus results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

0.3 mg/dl (0.1 mmol/l)

The lower detection limit represents the lowest measurable phosphorus concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Sample	Within run		
	Mean mmol/l	SD mmol/l	CV %
Control serum 1	1.17	0.011	0.94
Control serum 2	1.59	0.012	0.75
Control serum 3	2.08	0.019	0.91

Sample	Between Day		
	Mean mmol/l	SD mmol/l	CV %
Control serum 1	1.31	0.025	1.91
Control serum 2	1.68	0.034	2.02
Control serum 3	1.90	0.022	1.16

Fluitest[®] PHOS

INORGANIC PHOSPHORUS

Method comparison:

A comparison of the Analyticon Fluitest PHOS (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.991 x + 0.007; \quad r = 0.996$$

Quality Control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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