Precision Check

PRÄZISIONS-CHECK

BioLyzer® Order information:

Catalog No.	BioLyzer	Contents
B0020	200 / 600	R1 1 x 45 ml
		R2 1 x 22 ml
		R3 1 x 6 ml
		R4 1 x 5 ml
		R5 1 x 10 ml
B0030	300 / 600*	R1 1 x 45 ml
		R2 1 x 22 ml
		R3 1 x 6 ml
		R4 1 x 5 ml
		R5 1 x 10 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Precision test for validation of the BioLyzer 200, 300 and 600.

Summary:

This test should be carried out after transportation, regular maintenance by technician and exchange of parts.

Test principle:

An amido dye reagent/sample is pipetted by tested needle and precision is checked.

Reagent concentration:

R1 (SPT solution):

Buffer, pH 6
Detergents and stabilizers < 0.1 %

R2 (RPT 1 solution):

Buffer, pH 6
Amido black < 0.01%
Detergents and stabilizers < 0.1 %

R3 (RPT 2 solution):

Buffer, pH 6
Amido black < 0.07%
Detergents and stabilizers < 0.1 %

R4 (Sample SPT):

Buffer, pH 6
Amido black < 0.3%
Detergents and stabilizers < 0.1 %

R5 (Sample RPT):

Buffer, pH 6
Detergents and stabilizers < 0.1 %

Preparation and stability:

All solutions are ready to use.
The solutions are stable at + 4°C to + 25°C up to the date of expiration specified.
Avoid direct sunlight. Avoid contamination.

Notes:

For validation use only.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Materials provided
• Working solutions as described above
Additional materials required
• sample cups

Measuring:

Perform measurement according to Quick Guide.

Calculation:

The following values are calculated:
MV = Mean = sum of values / no. of measurements
SD = standard deviation
CV = variation coefficient = SD * 100 / Mean

Reference value:

The variation coefficient (CV) of results should be

	SPT	RPT1	RPT2
B200	< 2 %	< 2 %	-
B600	< 2 %	< 2 %	< 2 %
B300 w/o degasser	< 2 %	< 1 %	< 3 %
B300 with degasser	< 2 %	< 1 %	< 1 %

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Please note the legal regulations.

Quick Guide for Precision Test on Biolyzer 200:

1. Test settings are pre-installed on analyzer.
2. Menu "System Parameters – System Page 2/2": Choose 20 for number of replicate settings.
3. Place SPT reagent and RPT-1 reagent in reagent tray. Perform "RCU scan" in "Run Monitor – Inventory" to load reagent.
4. Select method for SPT test in Menu „Calibration – Calibration“. Press on "Parameter".
5. Choose "Factor" in pop-up window "Calculation" and insert 1000 in the "K" space. Save entered data by clicking „Calculate“.
6. Proceed for RPT test in the same way.
7. Put approx. 1ml SPT sample in sample cup and 1ml RPT sample in another cup. Place them on sample tray e.g. Pos. 1 and 2.
8. In Menu "Run Monitor – Test Selection" click on position for SPT sample e.g. pos 1. Click on "Type – R" for Replicate. Then select SPT test for SPT sample. "Save" test selection.
9. Select position for the RPT sample e.g. pos 2. Click on "Type – R" for Replicate. Then select RPT test for RPT sample. "Save" test selection.
10. "Start" analyzer

Quick Guide for Precision Test on Biolyzer 300:

1. Test settings are pre-installed on the analyzer (screen "chemistry parameters 1". If necessary, delete another test from "parameter-method no". Select the appropriate test from "method name" and save it onto this method no.
2. Place SPT reagent, RPT-1 and RPT-2 reagent in reagent tray. Perform "RCU scan" in "Run – Inventory" to load reagents.
3. Select method for SPT test in Menu „Calibration Registration“
4. Insert data as follows:
 - Sample Type: Serum
 - Replication: Duplicate
 - Check Interval: 0;
 - Calibration Type: Factor
5. Type "1000" in the "K" field and „Save" data.
6. Proceed for RPT1 and RPT2 test in same way.
7. Put approx. 1ml SPT sample in sample cup and 2 ml RPT sample in another cup. Place them on sample tray (outer ring).
8. In Menu "Run – Test Selection" click on position for SPT sample e.g. pos 1. Select "Normal Sample Cup" and "Serum" for sample type. Choose 20 for repeats. Then select SPT test for SPT sample. "Save" test selection.
9. Click on position for RPT sample. Select "Normal Sample Cup" and "Serum" for sample type. Choose 20 for repeats. Then select RPT-1 and RPT-2 test for RPT sample. "Save" test selection.
10. "Start" analyzer

Quick Guide for Precision Test on Biolyzer 600:

1. Test settings are pre-installed on analyzer (screen "chemistry parameters". If necessary, delete another test from "parameter-method no". Select the appropriate test from "method name" and save it onto this method no.
2. Place SPT reagent, RPT-1 and RPT-2 reagents in reagent tray. Perform "RCU scan" in "Run Monitor – Inventory" to load reagents.
3. Select method for SPT test in Menu „Calibration – Reg Calib“
4. Insert data as follows:
 - Sample Type: Common
 - Sampling: Duplicate
 - Check Interval: 0
 - Type: Factor
5. Type "1000" in the "K" field and „Save" data.
6. Proceed for RPT1 and RPT2 test in same way.
7. Put approx. 1ml SPT sample in a sample cup and 2 ml RPT sample in another cup. Place them on the sample tray (outer ring).
8. In Menu "Run Monitor – Test Selection" click on position for SPT sample. Select "Normal Sample Cup" and "Common" for sample type. Choose 20 for repetition. Then select SPT test for SPT sample. "Save" test selection.
9. Click on position for RPT sample. Select "Normal Sample Cup" and "Common" for sample type. Choose 20 for repetition. Then select RPT-1 and RPT-2 test for RPT sample. "Save" test selection.
10. "Start" analyzer

Text passages with grey background were changed in the latest revision of this package insert.

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