Order information:

Catalog No.	Contents
7540 50 Tests	R1 1 x 2.5 ml R2 1 x 0.5 ml R3 1 x 0.5 ml
7500 100 Tests	R1 1 x 5 ml R2 1 x 1 ml R3 1 x 1 ml
7040 Latex reagent	R1 1 x 2.5 ml
7000 Latex reagent	R1 1 x 5 ml
7011 Control Set	R2 1 x 0.5 ml R3 1 x 0.5 ml
7310 Reaction cards	1 x 20 pcs

Intended use:

A rapid slide test for the qualitative and semi-quantitative determination of rheumatoid factor in human serum.

Summary:

Rheumatoid factors are a heterogeneous group of autoantibodies directed against the antigenic determinants on the Fc-region of IgG molecules. They are important in the diagnosis of rheumatoid arthritis, but can also be found in other inflammatory-rheumatic diseases and in various non-rheumatic diseases. They are also found in clinically healthy persons over 60 years of age. Despite these restrictions, the detection of rheumatoid factors is a diagnostic criterion of the American College of Rheumatology for classifying rheumatoid arthritis. These autoantibodies occur in all the immunoglobulin classes, although the usual analytical methods are limited to the detection of rheumatoid factors of the IgM type.

The classic procedure for the quantitation of rheumatoid factors is by agglutination with IgG-sensitized sheep erythrocytes or latex particles. Particular problems of these semi-quantitative methods are the poor between-laboratory precision and reproducibility, together with standardization difficulties. For these reasons, new assay methods such as nephelometry, turbidimetry, enzyme-immunoassays and radioimmunoassays have been developed. The Analyticon RF assay is based on the immunological agglutination principle with enhancement of the reaction by latex particles.

Principle:

The RF reagent is a suspension of polystyrene latex particles sensitized with specially prepared human IgG in order to avoid unspecific agglutination. The RF latex reagent sensitivity has been adjusted to detect a minimum of 8 IU/ml of rheumatoid factors according with the WHO international standard without previous sample dilution.

Reagent Concentration:

R1:
Latex particles coated with anti-RF preservative

R2:
Protein containing solution with rheumatoid factors > 8 IU/ml preservative

R3:
Protein containing solution with rheumatoid factors < 8 IU/ml preservative

Preparation and stability

The reagents and control sera are stable:
Up to the expiration date at +2°C to +8°C Do not freeze!

Specimen:

Collect serum using standard sampling tubes
Use fresh clear serum specimens. Plasma, lipemic serum or microbial contamination may cause erroneous results.
If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 72 h.
For longer storage, freeze the serum.
Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.
The positive and negative controls were prepared from human sera which have been tested using FDA-approved methods and found to be non-reactive for HbsAg and HIV antibodies. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient sample. In the event of exposure the directives of the responsible health authorities should be followed.

Precautions:

Reagents containing sodium acid may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagent by flushing with large amounts of water to prevent azide build up.
Contaminated sera and a longer reaction time than 3 min. may cause false positive agglutination.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Procedure Qualitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently to homogeneity.
3. Place one drop (40 µl) of each test specimen and if reference is necessary, one drop (40 µl) of the pos./neg. control on **separate fields** of the reaction slide.
4. Add one drop Latex Reagent (40 µl) to each test field. Use a fresh spatula dispenser for each field to spread reaction mixture over entire test field.
5. Rotate the slide for 2 minutes and read immediately under optimal light conditions.

Procedure Semi-Quantitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently to homogeneity.
3. Using saline, dilute the specimens 1:2, 1:4, 1:8, 1:16 or as needed.
4. Place one drop (40 µl) of the test specimen followed by the corresponding prepared dilutions and as reference one drop (40 µl) of the pos./neg. control on **separate fields** of the reaction slide.
5. Add one drop of the Latex Reagent (40 µl) to each sample. Use a fresh spatula dispenser for each field to spread reaction mixture over entire test field.
6. Rotate the slide for 2 minutes and read immediately under optimal light conditions.

Results Qualitative Test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2).
A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

Results Semi-Quantitative Test:

The RF concentration of the serum lies between the highest dilution which just shows a positive reaction (e.g. 1:4) and the following dilution which shows a negative reaction (e.g. 1:8). For calculating the serum concentration, multiply the corresponding dilution factor by 8 IU/ml (cut off of the test).
e.g. Last positive signal 1:4 4x 8 IU/ml = 32 IU/ml
First negative signal 1:8 8x 8 IU/ml = 64 IU/ml
The RF conc. of the sample is between 32 and 64 IU/ml.

Interpretation:

Reference values:

Adults: < 8 IU/ml
The rheumatoid factors are immunoglobulins (most of the IgM) with antibody activity. These factors are present in most of patients suffering Rheumatoid Arthritis. There are different rheumatoid factors and does not exist any test capable to detect all of them, due to the fact that some of them act against human IgG, other against animal IgG, and other against both IgG. We recommend the use of WR test specific for detection of rheumatoid factors against animal IgG.

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the rheumatoid factor results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

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