Order information:

Catalog No.	Contents			
5748	R1	8 x	20 ml	R4 1 x 5 ml
5741	R1	4 x	100 ml	R4 1 x 5 ml
H9001 Hit I (ILab*)	R1	12 x	50 ml	
H9003 Hit 917 (AU*)	R1	12 x	60 ml	
AU9003 AU	R1	12 x	60 ml	

(*) Kit contains only reagent barcodes for Hitachi system

System information:

Hitachi 911: ACN 006
Hitachi 917: ACN 781

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:


Enzymatic *in vitro* test for the quantitative determination of triglycerides in human serum and plasma.

Summary:

Triglycerides are esters of the trihydric alcohol glycerol with 3 long-chain fatty acids. They are partly synthesized in the liver and partly ingested in food. The determination of triglycerides is utilized in the diagnosis and treatment of patients having diabetes mellitus, nephrosis, liver obstruction, lipid metabolism disorders and numerous other endocrine diseases. The enzymatic triglycerides assay as described by Eggstein and Kreutz still required saponification with potassium hydroxide. Numerous attempts were subsequently made to replace alkaline saponification by enzymatic hydrolysis with lipase. Bucolo and David tested a lipase/protease mixture; Wahlefeld used an esterase from the liver in combination with a particularly effective lipase from *Rhizopus arrhizus* for hydrolysis. This method is based on the work by Wahlefeld using a lipoprotein lipase from microorganisms for the rapid and complete hydrolysis of triglycerides to glycerol followed by oxidation to dihydroxyacetone phosphate and hydrogen peroxide. The hydrogen peroxide produced then reacts with 4-aminophenazone and 4-chlorophenol under the catalytic action of peroxidase to form a red dye (Trinder endpoint reaction).

Test principle:

Enzymatic colorimetric test:

Reagent concentration:

R1:	
Pipes buffer pH 7.8	50 mmol/l
p-Chlorophenole	2 mmol/l
Lipoprotein lipase	150000 U/l
Glycerolkinase	800 U/l
Glycerol - 3 - P- oxidase	4000 U/l
Peroxidase	440 U/l
4-Aminoantipyrine	0.7mmol/l
ATP	0.3mmol/l
Mg2+	40 mmol/l
Na-cholat	0.20 mmol/l
Potassium-Hexacyanoferrat(II)	1µmol/l

R4: (#5748, #5741)

Glycerol equivalent to a concentration of 200 mg/dl (2.28 mmol/l) Triglycerides.

Preparation and stability:

R1: Ready for use.

R4: Ready for use.

Unopened kit components: Up to the expiry date at +2°C to +8°C.

Onboard stability: R1 14 days (Hitachi)

Stability: 3 weeks at +20°C to +25°C.

Coloration of the reagent (reagent blank at 546 nm, 1cm >0.2) indicates a contamination or damage by storage at higher temperatures.

Specimen:

Collect serum using standard sampling tubes.

Heparinized or EDTA plasma.

Stability:

5 - 7 days at +2°C to +8°C

3 months at -20°C

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metamizole.

Notes:

For *in vitro* diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 22 (approximately 22 mg/dl bilirubin).

Hemolysis: No significant interference up to an index H of 625 (approximate hemoglobin concentration: 625 mg/dl).

Lipemia (Intralipid): Index L correlates with sample turbidity but not with triglyceride level. Extremely lipemic samples (triglycerides greater than 3000 mg/dl) can produce a normal result. Dilute such samples 1:4 with saline (0.9%) and multiply the result by 5. Roche/Hitachi 911: run the test with decrease a sample volume.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metamizole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength:	Hg 546 nm (500-550 nm)	
Temperature:	+37°C	
Cuvette:	1 cm	
Zero adjustment:	Reagent blank/ one reagent blank per series only	

	Blank	Sample/Calib./Stand.
Sample/Calib./Stand.	---	10 µl
R1	1000 µl	1000 µl

Mix and read initial absorbance. Incubate for 5 min at +37°C or for 10 min at +20°C to +25°C. Then read absorbance of sample against reagent blank within 60 min after start.

Calculation:

$$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Calib./Stand}} \times \text{Calib./Stand. conc.} = \text{Triglycerides conc.}$$

Measuring /reportable range:

3 - 1000 mg/dl (0.05 - 11.4 mmol/l)

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl (e.g. 1:4). Multiply the result by the appropriate dilution factor (e.g. 5)

Expected values:

Clinical interpretation as recommended by the European Arteriosclerosis Society.
lipid metabolism disorder

Cholesterol	<200 mg/dl	no
Triglycerides		
Cholesterol	200 – 300 mg/dl	yes if HDL-CHOL < 35 mg/dl
Cholesterol	> 300 mg/dl	yes
Triglycerides	> 200 mg/dl	

Expected values according to NCEP

Normal value < 200 mg/dl (2.30 mmol/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the triglycerides results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Note:

If the free glycerol is to be taken into account, then 10 mg/dl (0.11 mmol/l) must be subtracted from the triglycerides value obtained. For control sera, note the target value declaration of the manufacturer.

Analytical sensitivity (lower detection limit):

Detection limit: 3mg/dl (0.05 mmol/l)

The lower detection limit represents the lowest measurable triglycerides concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples the following results were obtained:

within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Control serum 1	122	1.09	0.89
Control serum 2	150	1.79	1.19
Control serum 3	206	1.44	0.70

between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Control serum 1	121	1.96	1.62
Control serum 2	161	1.84	1.14
Control serum 3	204	2.36	1.16

Method comparison:

A comparison of the Analyticon Fluitest TG (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.002x + 1.136; r = 0.999$$

Quality control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Lipid Control Serum:

Contronorm® L	5 x 2 ml	#1302
Contropath® L	5 x 2 ml	#1303

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardisation: The triglycerides GPO-PAP method was standardized against the ID-MS-method.

S1: 0.9% NaCl

S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420
Bio Cal® L	1 x 2 ml	#1401
	5 x 2 ml	#1402

Calibration frequency:

Two point calibration is recommended

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

Literature:

1. Bucolo G, David H. Clin Chem 1973;19:476
2. Eggstein M, Kreutz F, Klein Wschr 1966 44:262-267
3. Glick M.R., Ryder K.W., Jackson SA. Graphical comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem. 1986;32:470-474
4. Greiling H., Gressner A.M. (ed.). Lehrbuch der Klinischen Chemie und Pathobiochemie. 3rd. Stuttgart/New York: Schattauer Verlag, 1995
5. Shephard M.D.S., Whiting M.J. Clin Chem 1990; Vol 36, No 2,325-329. False Low Estimation of Triglycerides in Lipemic Plasma by the Enzymatic Triglyceride Method with Modified Trinder's Chromogen.
6. Siedel J et al. AACC Meeting Abstract 34. Clin Chem 1993;39:1127
7. Stein E.A., Myers G.L., National Cholesterol Education Program Recommendations for Triglyceridens Measurements: Executive Summary. Clin Chem 1995; 41:1421-1426
8. Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77
9. Trinder P. Ann Clin Biochem 1969;6:24
10. Wahlefeld A.W.; Bergmeyer H.U. (ed.) Methods of Enzymatic Analysis 2nd English ed. New York NY. academic Press Inc. 1974; 1831

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Literatur:

- Bucolo G. David H. Clin Chem 1973;19:476
- Eggstein M. Kreutz F. Klein Wschr 1966 44:262-267
- Glick M.R., Ryder K.W., Jackson SA. Graphical comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem. 1986;32:470-474
- Greiling H., Gressner A.M. (Hrsg.). Lehrbuch der Klinischen Chemie und Pathobiochemie. 3. Auflage. Stuttgart/New York: Schattauer Verlag, 1995
- Shephard M.D.S, Whiting M.J. Clin Chem 1990; Vol 36, No 2,325-329. False Low Estimation of Triglycerides in Lipemic Plasma by the Enzymatic Triglyceride Method with Modified Trinder's Chromogen.
- Siedel J. et al. AACC Meeting Abstract 34. Clin Chem 1993;39:1127
- Stein E.A., Myers G.L., National Cholesterol Education Program Recommendations for Triglyceridens Measurements: Executive Summaty. Clin Chem 1995; 41:1421-1426
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77
- Trinder P. Ann Clin Biochem 1969;6:24