BioLyzer® Order information:

Catalog No.	BioLyzer	Contents
B9551	200 / 600	R1 6 x 100 ml
B9553	300 / 600*	R1 6 x 60 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

In vitro test for the quantitative determination of total protein in human serum and plasma.

Summary:


Plasma proteins are synthesized predominantly in the liver, plasma cells, lymphnodes, spleen and in bone marrow. In the course of disease the total protein concentration and also the percentage represented by individual fractions can significantly deviate from normal values.

Hypoproteinemia can be caused by diseases and disorders such as loss of blood, sprue, nephrotic syndrome, severe burns, salt retention syndrome and Kwashiorkor (acute protein deficiency). Hyperproteinemia can be observed in cases of severe dehydration and illnesses such as multiple myeloma. Changes in the relative percentage of plasma proteins can be due to a change in the percentage of one plasma protein fraction. Often in such cases the amount of total protein does not change. The A/G-ratio is commonly used as an index of the distribution of albumin and globulin fractions. Marked changes in this ratio can be observed in cirrhosis of the liver, glomerulonephritis, nephrotic syndrome, acute hepatitis, lupus erythematosus as well as in certain acute and chronic inflammations. Total protein measurements are used in the diagnosis and treatment of a variety of diseases involving the liver, kidney, or bone marrow, as well as other metabolic or nutritional disorders.

Test principle:

Colorimetric assay.

Divalent copper reacts in alkaline solution with protein peptide bonds to form the characteristic purple-colored biuret complex. Sodium potassium tartrate prevents the precipitation of copper hydroxide and potassium iodide prevents autoreduction of copper.

The color intensity is directly proportional to the protein concentration which can be determined photometrically.

Reagent Concentration:

R1:	
Potassium iodide	30 mmol/l
Potassium sodium tartrate	32 mmol/l
Copper sulphate	18 mmol/l
Sodium hydroxide	200 mmol/l

Preparation and stability:

R1: Ready for use
 Unopened kit components: Up to the expiration date at +2°C to +8°C
 Onboard stability: R1 60 days

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma

Stability: 3 days at +2°C to +8°C
 6 months at -20°C

Separate the serum or plasma within 4 hours from the clot or cells.

The total protein concentration is by 0.4 to 0.8 mg/dl lower when the sample is collected from a patient situated in the recumbent position rather than upright.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial value.

Icterus: No significant interference up to an index I of 80 corresponding to approximately 80 mg/dl of bilirubin.

Hemolysis: No significant interference up to an index H of 350 (approximate hemoglobin concentration: 350 mg/dl).

This interference results from the fact that hemoglobin is reacting as a protein in the Total Protein assay.

Lipemia (Intralipid): No significant interference up to an index L of 105 (approximate triglycerides concentration: 210 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring range:

0.5 - 13 g/dl (5 - 130 g/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl as diluents.

Reference value:

Expected values according to Josephson

Adults: 6.6-8.7 g/dl or 66-87 g/l

For reference ranges for children, please refer to the book „Pediatric reference ranges“ 3rd ed. S.J. Soldin, C. Brugnara, I.M. Hicks AACCPress

Expected values according to Tietz

Unbilical Cord	4.8-8.0 g/dl
Pemature	3.6-6.0 g/dl
Newborn	4.6-7.0 g/dl
1 week	4.4-7.6 g/dl
7 months – 1 year	5.1-7.3 g/dl
1-2 years	5.6-7.5 g/dl
> 3 years	6.0-8.0 g/dl
adult (ambulatory)	6.4-8.3 g/dl

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, total protein results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.1 g/dl (1 g/l)

The detection limit represents the lowest measurable protein concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Sample	Within run		
	Mean g/dl	SD g/dl	CV %
Control serum 1	6.15	0.044	0.7
Control serum 2	4.87	0.031	0.6
Control serum 3	5.06	0.034	0.7

Sample	Run to run		
	Mean g/dl	SD g/dl	CV %
Control serum 1	6.52	0.177	2.7
Control serum 2	5.17	0.173	3.3
Control serum 3	7.85	0.151	1.9

Method comparison:

A comparison of the Analyticon Fluitest TP (y) with a commercial obtainable assay (x) gave following results:

$$y = 1.1503x + 2.5992; \quad r = 0.9914$$

Quality Control:

Human Control Serum

Contronorm® Plus 5 x 5 ml #1205
 20 x 5 ml #1220

Controptath® Plus 5 x 5 ml #1305
 20 x 5 ml #1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency:

Two point calibration is recommended

• after reagent lot change

• as required following quality control procedures

Calibration verification: Not necessary.

Disposal:

Please note the legal regulations.

Literature:

1. Brobeck J.R. (ed.). Physiological Basis of Medical Practice, 9th Baltimore, MD: Wilkins and Wilkins, 1973:4-7.
2. Glick M.R., Ryder K.W., Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474.
3. Josephson B, Gyllenswärd C. The Development of the Protein Fractions and of Cholesterol Concentration in the Serum of Normal Infants and Children. Scandinav J Clin Lab Investigation 1957; 9:29.
4. Koller A. Total serum protein. in Kaplan L.A., Pesce A.J. (ed.). Clinical Chemistry, Theory, Analysis, and Correlation. St. Louis: Mosby Company, 1984:1316-1319
5. Weichselbaum T.E. Amer J Clin Path 1946;16:40.

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