Turbitex® TRANS

TRANSFERRIN

BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B8901	200 / 600	R1	6 x	20 ml
		R2	6 x	7 ml
B8903	300 / 600*	R1	4 x	60 ml
		R2	4 x	20 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of transferrin in human serum and plasma.

Summary:

Transferrin is a glykoprotein with a molecular weight of 79570 daltons. It consists of a polypeptide stand with two N-glycosidically linked oligosaccharide chains and exists in numerous isoforms. The rate of synthesis in the liver can be altered in accordance with the body's iron requirements and iron reserves. Transferrin is the iron transport protein in serum. In cases of iron deficiency, the degree of transferrin saturation appears to be an extremely sensitive indicator of functional iron depletion. The ferritin levels are depressed when there is a deficiency of storage iron. In sideropenia, an iron deficiency can be excluded if the serum transferrin concentration is low, as in inflammations or - less commonly - in cases of ascorbic acid deficiency. In screening for hereditary hemochromatosis, transferrin saturation provides a better indication of the homozygous genotype than does ferritin. The treatment of anemia with erythropoietin in patients with renal failure is only effective when sufficient depot iron is present. The best monitoring procedure is to determine transferrin saturation during therapy. Transferrin saturation in conjunction with ferritin gives a conclusive prediction of the exclusion of iron overloading in patients with chronic liver disease. A variety of methods are available for determining transferrin including radial immunodiffusion, nephelometry and turbidimetry. This transferrin assay is based on the immunological agglutination principle.

Test principle:

Immunoturbidimetric assay
Anti-transferrin antibodies react with the antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity and reduces the risk of samples containing excess antigen producing false-negative results.

Reagent concentration:

R1:	
TRIS/HCl buffer* pH 7.8	100mmol/l
NaCl	150mmol/l
PEG	4.5 %
Preservative	
R2:	
Anti-human transferrin-antibodies (goat): dependent on titer;	
TRIS/HCl-buffer* pH 7.8	80 mmol/l
NaCl	150 mmol/l
Preservative	
*TRIS= Tris (hydroxymethyl)-aminomethane	

Preparation and Stability:

R1: Ready for use
R2: Ready for use
Unopened kit components: Up to the expiration date at 2-8°C
Onboard Stability at 2-8°C: R1 90 days
R2 90 days

Specimen:

Collect serum using standard sampling tubes

Heparin-Plasma

Stability: 8 days at +20°C to +25°C
8 days at +4°C to +8°C
6 months at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.
TIBC is the non-specific binding of iron ions to various serum proteins. Not only transferrin. For this reason, the conversion factor from transferrin to TIBC is not based on the molecular weight of transferrin, but has been experimentally determined.

The following conversion possibilities exist between the transferrin and the total iron-binding capacity:
1 mg transferrin cor. 1.41 µg Fe
mg transferrin/dl x 1.27 = µg TIBC/dl
mg transferrin/dl x 0.2275 = µmol TIBC/l

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to an index I of 100 (approximate conjugated and unconjugated bilirubin concentration: 100 mg/dl).
Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/d).
Lipemia (Intralipid): No significant interference up to an index L of 625 (approximate triglycerides concentration: 1250 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
Rheumatoid factors < 500 IU/ml do not interfere.
A high-dose-hok-effect may occur at transferrin concentrations of above 2000 mg/dl. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Measuring/reportable range:

15 – 550 mg/dl (0.15 – 5.5 g/l)

Determine samples containing higher transferrin concentrations via the rerun function using 0.9% NaCl.

Reference values:

IFCC/CRM 470*	
mg/dl	g/l
200-360	2.0-3.6

Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.
For diagnostic purposes the transferrin results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 1.0 mg/dl (0.01 g/l)

The detection limit represents the lowest measurable transferrin concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 21). The following results were obtained:

Sample	Within run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	123,62	1,415	1,2
Sample 2	246,21	5,397	2,2
Sample 3	358,35	16,021	4,5

Sample	Run to Run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	253,21	15,424	6,1
Sample 2	317,03	20,398	6,4
Sample 3	401,45	38,551	9,6

Method comparison:

A comparison of the Turbitex® TRANS (y) with a commercial obtainable assay (x) assay (x) gave with 20 samples the following result.
 $y = 0.7588 + 37.946x$; $r = 0.9702$

Quality control:

Protein Control Level 1 3 x 1 ml #7661
Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardisation: The transferrin method has been standardized against the reference preparation CRM 470.

Calibration Type: Logit-Log

S1: 0.9% NaCl

S2-S6: Bio Cal® P Calibration Set 5 x 1 ml #1475

Calibration frequency:

Full calibration is recommended:

• after lot change
• as required following quality control procedures

Turbitex[®] TRANS

TRANSFERRIN

Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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*TRIS= Tris (hydroxymethyl)-aminomethan-hydrochlorid

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