Order information:

Catalog No.		Contents			
H8901	Hit I/ 917 (ILab* / AU*)	R1	6 x	20 ml	R2 6 x 7 ml
H8903	Hit 917 (AU*)	R1	4 x	60 ml	R2 4 x 20 ml
AU8901	AU	R1	6 x	20 ml	R2 6 x 7 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 216
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of transferrin in human serum and plasma.

Summary:

Transferrin is a glykoprotein with a molecular weight of 79570 daltons. It consists of a polypeptide strand with two N-glycosidically linked oligosaccharide chains and exists in numerous isoforms. The rate of synthesis in the liver can be altered in accordance with the body's iron requirements and iron reserves. Transferrin is the iron transport protein in serum. In cases of iron deficiency, the degree of transferrin saturation appears to be an extremely sensitive indicator of functional iron depletion. The ferritin levels are depressed when there is a deficiency of storage iron. In sideropenia, an iron deficiency can be excluded if the serum transferrin concentration is low, as in inflammations or - less commonly - in cases of ascorbic acid deficiency. In screening for hereditary hemochromatosis, transferrin saturation provides a better indication of the homozygous genotype than does ferritin. The treatment of anemia with erythropoietin in patients with renal failure is only effective when sufficient depot iron is present. The best monitoring procedure is to determine transferrin saturation during therapy. Transferrin saturation in conjunction with ferritin gives a conclusive prediction of the exclusion of iron overloading in patients with chronic liver disease. A variety of methods are available for determining transferrin including radial immunodiffusion, nephelometry and turbidimetry. This transferrin assay is based on the immunological agglutination principle.

Test principle:

Immunoturbidimetric assay
• sample and addition of R1 (buffer)
• addition of R2 (anti-transferrin antibody/buffer) and start of reaction
Anti-transferrin antibodies react with the antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity and reduces the risk of samples containing excess antigen producing false-negative results.

Reagent concentration:

R1:	
TRIS/HCl buffer* pH 7.8	100mmol/l
NaCl	150mmol/l
PEG	4.5 %
Preservative	
R2:	
Anti-human transferrin-antibodies (goat): dependent on titer;	
TRIS/HCl-buffer* pH 7.8	80 mmol/l
NaCl	150 mmol/l
Preservative	
*TRIS= Tris (hydroxymethyl)-aminomethane	

Preparation and Stability

R1: Ready for use
R2: Ready for use
Unopened kit components: Up to the expiration date at 2-8°C
Onboard Stability at 2-8°C: R1 90 days
R2 90 days

Specimen

Collect serum using standard sampling tubes
Heparin-Plasma
Stability: 8 days at +20°C to +25°C
8 days at +4°C to +8°C
6 months at -20°C
Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to an index I of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl).
Hemolysis: No significant interference up to an index H of 1000 (approximate haemoglobin concentration: 1000 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 600 (approximate triglyceride concentration: 1200 mg/dl). There is poor correlation between turbidity and triglyceride concentration.
Rheumatoid factors < 500 IU/ml do not interfere.
A high-dose-hook-effect may occur at transferrin concentrations of above 2000 mg/dl. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:		
Wavelength:	600 nm	
Temperature:	37°C	
Cuvette:	1 cm	
Zero adjustment:	against blank	
Sample/Calibrator	Blank	Sample/Calibrator
R1	750 µl	10 µl 750 µl
Mix and incubate 5 min. at 37°C and read A ₁ . Then add:		
R2	200 µl	200 µl
Mix and incubate 5 min. at 37°C and read A ₂		
Calculation:		
$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$		
The concentration of transferrin in patient sera has to be calculated from ΔA using mathematic function as logit/log or can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 160 mg/dl transferrin. For zero value is recommended to use saline solution (0.9%)		
TIBC is the non-specific binding of iron ions to various serum proteins, not only transferrin. For this reason, the conversion factor from transferrin to TIBC is not based on the molecular weight of transferrin, but has been experimentally determined.		
The following conversion possibilities exist between the transferrin and the total iron-binding capacity: 1 mg transferrin cor. 1.41 µg Fe mg transferrin/dl x 1.27 = µg TIBC/dl mg transferrin/dl x 0.2275 = µmol TIBC/l		

Measuring/reportable range:

10-500 mg/dl (0.1-5.0 g/l)

At higher concentrations, determine via rerun function or manually dilute the sample with 0,9% NaCl solution, (e.g. 1:1). Multiply the result by the appropriate factor (e.g. Factor 2).

Reference values:

IFCC/CRM 470*
200-360 mg/dl (2.0-3.6 g/l)

*Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the transferrin results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 15 mg/dl (0.15 g/l)
The detection limit represents the lowest measurable transferrin concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 21). The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	124	2.15	1.72
Sample 2	253	4.42	1.75
Sample 3	372	5.75	1.54

Method comparison:

A comparison of the Turbitex[®] TRANS (y) with a commercial obtainable assay (x) assay (x) gave with 40 samples the following result.

$$y = 1.017 + 0.969 x; \quad r = 0.988$$

Quality control:

Protein Control Level 1 3 x 1 ml #7661
Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardisation: The transferrin method has been standardized against the reference preparation CRM 470.

• Manual procedure:

S1: 0.9% NaCl
S2-S6: Bio Cal P 3 x 1 ml #1407
 Bio Cal P Set 5 x 1 ml #1475

• Roche/Hitachi 911/917:

S1: 0.9% NaCl
S2-S6: Bio Cal P 3 x 1 ml #1407

Factors for calculating the concentrations of standards for the 6-point calibration curve using the lot-assigned value of the Bio Cal P:

S2: 0.125 S5: 1.000
S3: 0.250 S6: 2.000
S4: 0.500

Calibration frequency

Full calibration is recommended:

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Literatur:

1. Bablok W et al. General Regression Procedure for Method Transformation. J Clin Chem Biochem 1988; 26:783-790
2. Buffone GJ, Lewis SA, Josefsohn M et al. Chemical and immuno-chemical measurement of total iron-binding capacity compared. Clin Chem 1978; 24:1788-1791
3. Consensus values of the Deutsche Gesellschaft für Laboratoriumsmedizin, the Deutschen Gesellschaft für Klinische Chemie and the Verband der Diagnostica-Industrie e.V. (VDGH). Clin Lab 1995;41:743-748
4. Gottschalk R. Evaluation of total iron binding capacity and transferrin determination under several clinical conditions. Publication in preparation.
5. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag, 1996
6. Haupt H, Baudner S. Behring Inst Mitt 1990; 86: 16-19
7. Heidelberger M, Kendall FE. J Exp Med 1935;62:697
8. Kreuzer HJH. An immunological turbidimetric method for serum transferrin determination. J Clin Chem Biochem 1976;14:401-406
9. Lievens M, Borque de Larrea L, Rus A et al. Determination of transferrin without sample predilution by a turbidimetric assay on Boehringer Mannheim/Hitachi analysis systems. Lab med 1994;18:83-90
10. Lizana J Hellsing K Clin Chem 1974;20:1181
11. Mancini G, Carbonara AO, Heremans JF. Immunochemical Quantitation of antigens by single radial Immunodiffusion. Immunochemistry 1965;2:235-243
12. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two different Analytical Methods. J Clin Chem Clin Biochem 1983;21: 709-720
13. Tietz NW. Fundamentals of Clinical Chemistry, 2nd ed. Philadelphia, Pa: WB Saunders, 1976:278-280
14. Wick M, Pinggera W, Lehmann P eds.. Iron Metabolism, Diagnosis and Therapy of Anemias, 3th Wien/New York: Springer Verlag 1996

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