Order information:

Catalog No.	Contents						
4852	R1	8 x	15 ml	R2	1 x	25 ml	
				R4	1 x	5 ml	
H9801	Hit I (ILab*)	R1	6 x	47 ml	R2	6 x	11 ml
H9803	Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	14 ml
AU9803	AU	R1	6 x	60 ml	R2	6 x	14 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 044
Hitachi 917: ACN 700

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Enzymatic *in vitro* test for the quantitative determination of uric acid in human serum and plasma.

Summary:

Uric acid is the final product of purine metabolism in the human organism. Uric acid measurements are used in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukaemia, psoriasis, starvation or other wasting conditions, and of patients receiving cytotoxic drugs.


The oxidation of uric acid provides for two approaches to the quantitative determination of this purine metabolite. One approach is the reduction of phosphotungstic acid in an alkaline solution to tungsten blue, which is measured photometrically. The method is, however, subject to interferences from drugs and reducing substances other than uric acid. A second approach, described by Praetorius and Poulson, utilizes the enzyme uricase to oxidise uric acid; this method eliminates the interferences intrinsic to chemical oxidation. Uricase can be employed in methods that involve the UV measurement of the consumption of uric acid or in combination with other enzyme to provide a colorimetric method. Another method is the colorimetric method developed by Town, et al. The sample is initially incubated with a reagent mixture containing ascorbate oxidase and a clearing system. In this test system it is important that any ascorbic acid is eliminated in the preliminary reaction; this precludes any ascorbic acid interferences with the subsequent POD indicator reaction. Upon addition of the starter reagent, oxidation of uric acid begins.

The assay described here is a slight modification of the colorimetric method described above. The modifications were described by Siedel.

Test principle:

Kinetic colorimetric assay

- Sample and addition of R1 (buffer/enzyme/DHBSA)
- Addition of R2 (buffer/enzymes/aminoantipyrene) and start of reaction

Uricase cleaves uric acid to form allantoin and hydrogen peroxide.

The peroxide reacts in the presence of peroxidase, DHBSA and aminoantipyrene to form a Quinonediimine dye. The intensity of the red color is proportional to the uric acid concentration and is determined photometrically.

Reagent concentration:

R1:
Phosphate buffer pH 7.4 50 mmol/l
DHBSA* 7 mmol/l
Preservative

R2:
Uricase 6 kU/l
POD 5 kU/l
Aminoantipyrene 1 mmol/l
Preservative

R4: (#4852)
Uric acid 6 mg/dl (356.9 µmol/l)

* 3,5-Dichloro-2-hydroxy-benzenesulfonic acid

Preparation and stability:

R1: ready for use

R2: ready for use

The reagents are stable up to the expiry date when stored at +2°C to +8°C.

Onboard stability: R1: 28 days Protect from light
R2: 28 days Protect from light

Notes:

For *in vitro* diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Specimen:

Serum/plasma

Collect serum using standard sampling tubes.

Heparin, or EDTA-plasma

Stability: 5 days at +2 to +8°C
6 months at -20°C

Urine

Collect urine without using preservatives.

Stability: Assay urinary uric acid as soon as possible. Do not refrigerate.

• Roche/Hitachi analyzers with automatic sample dilution:

Urine samples are diluted 1 : 10 with system water or 0.9% NaCl. This dilution is taken into account in the calculation of the results:

• Roche/Hitachi analyzers without sample dilution:

Dilute urine samples manually with distilled water or 0.9% NaCl (e.g. 1 : 10). Multiply the result by the appropriate dilution factor (e.g. 11).

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of **Metamizole**.

Limitations - interference:

Criterion: Recovery within ±10% of initial values

Icterus: No significant interference up to an index I of 40 (approximate conjugated and unconjugated bilirubin concentration: 40 mg/dl)

Hemolysis: No significant interference up to an H Index of 1000 (approximate hemoglobin concentration: 1000 mg/dl)

Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Elevated levels of ascorbic acid produce false low values.

Uricase reacts specifically with uric acid. Other purine derivatives can inhibit the uric acid reaction.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metamizole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure for substrate start:			
Wavelength:	546 nm		
Temperature:	+25 / +30 / +37°C		
Cuvette:	1 cm light path		
Zero adjustment:	against reagent blank		
	Blank	Standard	Sample
Calibrator	---	15 µl	---
Sample	---	---	15 µl
R1	500 µl	500 µl	500 µl
Mix, incubate for 1 min. at assay temperature, read initial absorbance, then add:			
R2	100 µl	100 µl	100 µl
Mix, start stopwatch simultaneously. Repeat reading after exactly 4 minutes.			
Calculation:			
$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{Uric acid conc.}$			

Measuring /reportable range:

0.2 – 25.0 mg/dl (11.9 – 1487 µmol/l)

Unit conversion: mg/dl x 59.5 = µmol/l
mg/dl x 0.059 = mmol/l

Determine samples with uric acid concentrations > 25.0 mg/dl via the rerun function. On instruments without rerun function, dilute the samples manually with 0.9% NaCl or distilled/ deionized water (e.g. 1 : 1). Multiply the result by the appropriate dilution factor (e.g. factor 2).

Expected values:

Serum/plasma:

Male: 3.4 – 7.0 mg/dl (202.3 – 416.5 µmol/l)

Female: 2.4 – 5.7 mg/dl (142.8 – 339.2 µmol/l)

Urine

(Reference range according to Krieg and Colombo)

Morning urine 37 – 92 mg/dl (2200 – 5475 µmol/l)

24 hour urine 200 – 1000 mg/24h (1200 – 5900 µmol/l/24h)

corresponding to 13 – 67 mg/dl* (773 – 3986 µmol/l*)

*Calculated from a urine volume of 1.5 l/24h

Urine (Reference range according to Tietz)

Average diet: 250 – 750 mg/24h

Low purine diet: Male: < 480 mg/24h

Female: < 400 mg/24h

High purine diet: < 1000 mg/24h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the uric acid results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.2 mg/dl (11.9 µmol/l)

The lower detection limit represents the lowest measurable uric acid concentration that can be distinguished from zero.

Imprecision:

Serum

Reproducibility was determined using human samples and controls (n = 20).

The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Human serum	6.0	0.04	0.8
Contronorm	4.23	0.03	0.7
Controptath	11.49	0.06	0.5

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Human serum	7.21	0.12	1.7
Contronorm	4.43	0.04	0.9
Controptath	11.03	0.13	1.2

Method comparison:

A comparison of the Analyticon Fluitest® UA 5+1 (y) with a commercial obtainable assay (x) gave with 44 samples the following result:

$$y = 1.008 x + 0.0645; \quad r = 0.995$$

Quality Control:

Human Control Serum:

Contronorm® Plus 5 x 5 ml #1205

20 x 5 ml #1220

Controptath® Plus 5 x 5 ml #1305

20 x 5 ml #1320

Human Urine Control:

Urine control Set 8 x 5 ml #1507

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Bio Cal® 20 x 3 ml #1420

Calibration frequency:

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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Grey text passages were changed in the latest revision of this package insert.

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