Fluitest[®] UA

URIC ACID

BioLyzer[®] Order information:

Catalog No.	BioLyzer	Contents
B9801	200 / 600	R1 6 x 47 ml
		R2 6 x 11 ml
B9833	300 / 600*	R1 6 x 55 ml
		R2 6 x 14 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Enzymatic in vitro test for the quantitative determination of uric acid in human serum and plasma.

Summary:

Uric acid is the final product of purine metabolism in the human organism. Uric acid measurements are used in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukaemia, psoriasis, starvation or other wasting conditions, and of patients receiving cytotoxic drugs.


The oxidation of uric acid provides for two approaches to the quantitative determination of this purine metabolite. One approach is the reduction of phosphotungstic acid in an alkaline solution to tungsten blue, which is measured photometrically. The method is, however, subject to interferences from drugs and reducing substances other than uric acid.

A second approach, described by Praetorius and Poulson, utilizes the enzyme uricase to oxidise uric acid; this method eliminates the interferences intrinsic to chemical oxidation. Uricase can be employed in methods that involve the UV measurement of the consumption of uric acid or in combination with other enzyme to provide a colorimetric method.

The assay described here is a slight modification of the colorimetric method. The modifications were described by Siedel. This reaction, the peroxide reacts in the presence of peroxidase, DHBSA and aminoantipyrine to form a quinoneimine dye. The intensity of the red color is proportional to the uric acid concentration and is determined photometrically.

Test principle:

Colorimetric assay

Uricase cleaves uric acid to form allantoin and hydrogen peroxide.

The peroxide reacts in the presence of peroxidase, DHBSA and aminoantipyrine to form a quinoneimine dye. The intensity of the red colour is proportional to the uric acid concentration and is determined photometrically.

Reagent concentration:

R1:	
Phosphate buffer pH 7.4	50 mmol/l
DHBSA*	7 mmol/l
Preservative	
R2:	
Uricase	6 kU/l
POD	5 kU/l
Aminoantipyrine	1 mmol/l
Preservative	

* 3,5-Dichloro-2-hydroxy-benzenesulfonic acid

Preparation and stability:

R1: ready for use
R2: ready for use

The reagents are stable up to the expiry date in the label when stored at +2°C to +8°C.

Onboard stability: R1: 28 days Protect from light
R2: 28 days Protect from light

Specimen:

Serum/plasma

Collect serum using standard sampling tubes.

Heparin, or EDTA-plasma

Stability: 5 days at +2°C to +8°C
6 months at -20°C

Urine

Collect urine without using preservatives.

Stability: Assay urinary uric acid as soon as possible. Do not refrigerate.

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metaminazole.

Note:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial values

Icterus: No significant interference up to an index I of 22 (approximate conjugated and unconjugated bilirubin concentration 22 mg/dl)

Hemolysis: No significant interference up to an index H of 185 (approximate hemoglobin concentration: 185 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 80 (approximate triglycerides concentration: 160 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Elevated levels of ascorbic acid produce false low values.

Uricase reacts specifically with uric acid. Other purine derivatives can inhibit the uric acid reaction.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminazole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Unit conversion: mg/dl x 59.5 = µmol/l
mg/dl x 0.059 = mmol/l

Measuring /reportable range:

0.5 – 30 mg/dl (30 - 1785 µmol/l)

Determine samples with higher concentrations via the rerun function using 0.9% NaCl as diluents.

Expected values:

Serum/plasma:

Male: 3.4 – 7.0 mg/dl (202.3 – 416.5 µmol/l)

Female: 2.4 – 5.7 mg/dl (142.8 – 339.2 µmol/l)

Urine

(Reference range according to Krieg and Colombo)

Morning urine 37 – 92 mg/dl (2200 – 5475 µmol/l)

24 hour urine 200 – 1000 mg/24h (1200 – 5900 µmol/l/24h)

corresponding to 13 – 67 mg/dl* (773 – 3986 µmol/l*)

* Calculated from a urine volume of 1.5 l/24h

Urine (Reference range according to Tietz)

Average diet: 250 – 750 mg/24h

Low purine diet: Male: < 480 mg/24h

Female: < 400 mg/24h

High purine diet: < 1000 mg/24h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the uric acid results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.1 mg/dl (6 µmol/l)

The lower detection limit represents the lowest measurable uric acid concentration that can be distinguished from zero.

Imprecision:

Serum

Reproducibility was determined using human samples and controls (n = 20). The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	4.69	0.035	0.7
Sample 2	11.51	0.08	0.7
Sample 3	5.89	0.049	0.8

Reproducibility was determined using human samples and controls (n = 20). The following results were obtained:

Sample	Run to run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	10.87	0.115	1.1
Sample 2	5.31	0.065	1.2
Sample 3	10.60	0.145	1.4

Fluitest[®] UA

URIC ACID

Method comparison:

A comparison of the Analyticon Fluitest[®] UA (y) with a commercial obtainable assay (x) gave with 148 samples the following result:
 $y = 0.9295x + 0.8719$; $r = 0.9858$

Quality Control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine Control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Calibration Type: Linear

S1:	0.9% NaCl		
S2:	Bio Cal [®] E	10 x 3 ml	# 1430

Calibration frequency:

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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