### Order information:

Catalog No.	Contents	
4652	R1 1 x 110 ml	R2 5 x for 20 ml
	R4 1 x 5 ml	

### Intended use:

Enzymatic in vitro test for the quantitative determination of uric acid in human serum and plasma.

### Summary:

Uric acid is the final product of purine metabolism in the human organism. Uric acid measurements are used in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukaemia, psoriasis, starvation or other wasting conditions, and of patients receiving cytotoxic drugs.


The oxidation of uric acid provides the basis for two approaches to the quantitative determination of this purine metabolite. One approach is the reduction of phosphotungstic acid in an alkaline solution to tungsten blue, which is measured photometrically. The method is, however, subject to interferences from drugs and reducing substances other than uric acid.

A second approach, described by Praetorius and Poulson, utilizes the enzyme uricase to oxidise uric acid; this method eliminates the interferences intrinsic to chemical oxidation. Uricase can be employed in methods that involve the UV measurement of the consumption of uric acid or in combination with other enzymes to provide a colorimetric method.

The assay described here is a slight modification of the colorimetric method. The modifications were described by Siedel. The peroxide reacts in the presence of peroxidase, phenole and amino-antipyrine to form a quinoneimine dye. The intensity of the red color is proportional to the uric acid concentration and is determined photometrically.

### Test principle:

Colorimetric endpoint assay

Uricase cleaves uric acid to form allantoin and hydrogen peroxide.

The increase in absorbance is measured.

### Reagent concentration:

<b>R1:</b>	
Phosphate buffer pH 8.0	50 mmol/l
Chlorophenol	7.5 mmol/l

<b>R2:</b>	
Uricase	300 U/l
POD	> 1 KU/l
4-aminoantipyrine	3 mmol/l

<b>R4:</b>	
Uric acid	6 mg/dl (356.9 µmol/l)

### Preparation and stability:

Dissolve contents of enzyme reagent/R2 with the corresponding volume of buffer /R1.

Gently swirl until completely dissolved. DO NOT SHAKE.

This working reagent is stable (protected from light!):

3 days	at +20 to +25°C
14 days	at + 2 to + 8°C

### Specimen:

#### Serum/plasma

Collect serum using standard sampling tubes.

Heparin, or EDTA-plasma

Stability:	5 days	at +2°C to +8°C
	6 months	at - 20°C

#### Urine

Stability: Assay urinary uric acid as soon as possible. Do not refrigerate.

Stability from urine with using NaOH (pH > 8.0):

4 days	at + 20 to 25°C
--------	-----------------

Dilute urine samples manually with distilled water or 0.9% NaCl - solution (e.g. 1 + 10). Multiply the result by the appropriate dilution factor (e.g. 11)

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metamizole.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial values

Icterus: No significant interference up to an index I of 55 (approximate conjugated and unconjugated bilirubin concentration 55 mg/dl)

Hemolysis: No significant interference up to an index H of 25 (approximate hemoglobin concentration: 25 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 300 (approximate triglycerides concentration: 600 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Elevated levels of ascorbic acid produces false low values.

Uricase reacts specifically with uric acid. Other purine derivatives can inhibit the uric acid reaction.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metamizole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request

#### Materials provided

• Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl - solution

### Manual procedure:

Wavelength:	Hg 546 nm, (490 - 550)
Temperature:	+25 / +30 / +37°C
Cuvette:	1 cm light path
Zero adjustment:	against reagent blank one reagent blank per series only

	Blank	Calibrator	Sample
Sample	---	40 µl	40 µl
Working reagent	1000 µl	1000 µl	1000 µl

Mix, incubate for 5 min. at +37°C or 10 min at +20°C or +25°C.  
Read absorbance of the sample against reagent blank within 30 min.

### Calculation:

By standard:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{Uric acid conc.}$$

### Measuring /reportable range:

0.2 – 40.0 mg/dl (11.9 – 2380 µmol/l)

Determine samples having higher concentrations via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. factor 2).

Unit conversion:

$$\begin{aligned} \text{mg/dl} \times 59.5 &= \mu\text{mol/l} \\ \text{mg/dl} \times 0.059 &= \text{mmol/l} \end{aligned}$$

### Expected values:

#### Serum/plasma:

Male: 3.4 - 7.0 mg/dl (202.3 - 416.5 µmol/l)

Female: 2.4 - 5.7 mg/dl (142.8 - 339.2 µmol/l)

#### Urine (Reference range according to Krieg and Colombo)

Morning urine

37 - 92 mg/dl (2200 - 5475 µmol/l)

24 hour urine

200 - 1000 mg/24h (1200 - 5900 µmol/24h)

corresponding to 13 - 67 mg/dl\* (773 - 3986 µmol/l\*)

\*Calculated from a urine volume of 1.5 l/24h

#### Urine (Reference range according to Tietz)

Average diet: 250 - 750 mg/24h

Low purine diet: Male: < 480 mg/24h

Female: < 400 mg/24h

High purine diet: < 1000 mg/24h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the uric acid results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 0.2 mg/dl (11.9 µmol/l)

The lower detection limit represents the lowest measurable uric acid concentration that can be distinguished from zero.

### Imprecision:

#### Serum

Reproducibility was determined using human samples and controls within run (n = 20). The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample1	4.34	0.029	0.67
Sample2	5.92	0.056	0.95
Sample3	11.63	0.062	0.53

Reproducibility was determined using human samples and controls between day (n = 20). The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample1	4.73	0.097	2.05
Sample2	6.78	0.116	1.71
Sample3	11.03	0.315	2.86

### Method comparison:

A comparison of the Analyticon UA (y) with a commercially available assay (x) gave with 38 samples the following result (mg/dl):

$$y = 0.999x - 0.034; \quad r = 0.997$$

### Quality Control:

#### Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
ControPath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

#### Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

S1:	0.9% NaCl		
S2:	Bio Cal®	20 x 3 ml	#1420
	Bio Cal® E	10 x 3 ml	#1430
	R4: calibrator provided in kit		

### Calibration frequency

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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